

**NEW!**

**smiths medical**  
bringing technology to life

## SafeSheath® Ultra Lite Valved Peel-Away Introducers




**For lower risk of blood loss and air embolisms**

**Deltec**

# SafeSheath® Ultra Lite Valved Peel-Away Introducers

Smiths Medical now offers a new valved peel-away introducer for lower risk of blood loss and air embolisms. Available in a wide variety of sizes and for use with PORT-A-CATH® implantable access systems.



- Low insertion force valve membrane allows smooth insertion and removal of port catheters
- Extruded scoreline peel-away sheath for clean, even peeling during removal
- Smooth guidewire/dilator transition

## Order Information

SafeSheath® introducers are available in Smiths Medical PORT-A-CATH® system trays and also as separate sets.

### PORT-A-CATH® II Systems with SafeSheath® Valved Peel-Away Introducers

	Introducer Size	Catheter Material	OD(mm)	Length(cm)	Connector	Tray Reorder #
Single-lumen	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4066
Single-lumen	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK®	21-4067
Single-lumen Low Profile™	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4068
Single-lumen Low Profile™	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK®	21-4069
Dual-lumen Low Profile™	7 Fr	Polyurethane	2.2	76	WING-LOCK®	21-8069
Dual-lumen Low Profile™	10 Fr	Polyurethane	3.2	76	WING-LOCK®	21-8070

### PORT-A-CATH® POWER P.A.C. Systems with SafeSheath® Valved Peel-Away Introducers

	Introducer Size	Catheter Material	OD(mm)	Length(cm)	Connector	Tray Reorder #
Single-lumen	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4474
Single-lumen Low Profile™	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4475
Single-lumen	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK®	21-4476
Single-lumen Low Profile™	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK®	21-4477
Dual-lumen Low Profile™	7 Fr	Polyurethane	2.2	76	WING-LOCK®	21-8469
Dual-lumen Low Profile™	10 Fr	Polyurethane	3.2	76	WING-LOCK®	21-8470

### Introducers only

#### SafeSheath® Valved Peel-Away Introducer Sets

French Size	Reorder # (5/box)
6 Fr	21-2476
7 Fr	21-2477
8.5 Fr	21-2478
9 Fr	21-2479
10 Fr	21-2480
11 Fr	21-2481

For additional information, visit [www.smiths-medical.com](http://www.smiths-medical.com) > products > vascular access > port accessories

Introducer sheath length – 13 cm Dilator sheath length – 18 cm

**Smiths Medical ASD, Inc.**  
St. Paul, MN 55112, USA  
Phone: 1-214-618-0218  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at:  
[www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

PORT-A-CATH and the Deltec and Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks, or service marks of their respective owners. ©2011 Smiths Medical. All rights reserved. 10/11 VA192200

**smiths medical**