





















CADD™ Ambulatory Tubing Sets






Rev 6/06 19830





Medication Cassette Reservoirs <i>Must be used with CADD® Extension Set</i>							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7001	12/box	Medication Cassette Reservoir with clamp and female luer – 50ml	125 ml/hr	0.2 ml	8"	DEHP
	21-7002	12/box	Medication Cassette Reservoir with clamp and female luer – 100ml	125 ml/hr	0.2 ml	8"	DEHP
	21-7100	12/box	Yellow Medication Cassette Reservoir with clamp and female luer – 100ml	125 ml/hr	0.2 ml	8"	DEHP



Medication Cassette Reservoirs with Flow Stop Free-Flow Protection <i>Must be used with CADD® Extension Set</i>							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7301*	12/box	Medication Cassette Reservoir with flow stop, clamp, and female luer – 50ml	125 ml/hr	0.2 ml	8"	DEHP
	21-7302*	12/box	Medication Cassette Reservoir with flow stop, clamp, and female luer – 100ml	125 ml/hr	0.2 ml	8"	DEHP
	21-7300*	12/box	Yellow Medication Cassette Reservoir with flow stop, clamp, and female luer – 100ml	125 ml/hr	0.2 ml	8"	DEHP

CADD™ Extension Sets designed to be used with a Medication Cassette Reservoir							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7045 21-7060	50/box 12/box	CADD™ Extension Set with male luer, clamp, and integral anti-siphon valve with male luer	125 ml/hr	1.1 ml	30" 30"	TOTM (non-DEHP)
	21-7046 21-7061	50/box 12/box	CADD™ Extension Set with male luer, clamp, and integral anti-siphon valve with male luer	125 ml/hr	1.4 ml	45" 45"	TOTM (non-DEHP)
	21-7047 21-7062	50/box 12/box	CADD™ Extension Set with male luer, clamp, and integral anti-siphon valve with male luer	125 ml/hr	1.7 ml	60" 60"	TOTM (non-DEHP)
	21-7092	12/box	CADD™ Extension Set with male luer, integral anti-siphon valve with male luer, y-extension with one-way female checkvalve, and clamps	125 ml/hr	2.9 ml	68"	TOTM (non-DEHP)
	21-7105	12/box	CADD™ Extension Set with male luer, yellow-striped tubing, clamp, and integral anti-siphon valve with male luer	125 ml/hr	2.3 ml	90"	TOTM (non-DEHP)
	21-7040	25/box	CADD™ Extension Set with male luer, clamp, 0.2-micron air-eliminating filter, and integral anti-siphon valve with male luer	125 ml/hr	4.6 ml	55"	TOTM (non-DEHP)
	21-7106	50/box	CADD™ Extension Set with male luer, 0.2-micron air-eliminating filter, clamp, and integral anti-siphon valve with male luer	125 ml/hr	1.8 ml	60"	TOTM (non-DEHP)

CADD™ Administration Sets with Flow Stop Free-Flow Protection <i>For use with CADD-Legacy® and CADD-Prizm® pumps</i>							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7321	12/box	CADD™ Administration Set with female luer, flow stop, clamp, and one-way checkvalve with male luer (for use with a syringe)	125 ml/hr	2.7 ml	90"	TOTM (non-DEHP)
	21-7322	12/box	CADD™ Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	125 ml/hr	2.3 ml	60"	TOTM (non-DEHP)
	21-7324	12/box	CADD™ Yellow Striped Administration Set with bag spike, flow stop, yellow striped tubing, clamp, and one-way checkvalve with male luer	125 ml/hr	3.2 ml	105"	TOTM (non-DEHP)
	21-7359	12/box	CADD™ Administration Set with male luer, flow stop, clamp, and one-way checkvalve with male luer (for use with 21-6167)	125 ml/hr	2.0 ml	60"	TOTM (non-DEHP)
	21-7390	12/box	CADD™ Administration Set with female luer, flow stop, 7.6cm y-extension, one-way checkvalve with female luer, one-way checkvalve with male luer, and clamps (for use with a syringe)	125 ml/hr	4.0 ml	90"	TOTM (non-DEHP)
	21-7391	12/box	CADD™ Administration Set with bag spike, flow stop, 7.6cm y-extension, one-way checkvalve with female luer, one-way checkvalve with male luer, and clamps	125 ml/hr	4.2 ml	90"	TOTM (non-DEHP)
	21-6167	48/box	250 ml TOTM-plasticized medication reservoir with female luer-activated valve	-----	-----	-----	TOTM (non-DEHP)

CADD™ Administration Sets with Air-eliminating Filters and Flow Stop Free-Flow Protection <i>For use with CADD-Legacy® and CADD-Prizm® pumps</i>							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7336	15/box	CADD™ Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male luer	125 ml/hr	10.0 ml	75"	TOTM (non-DEHP)
	21-7383	15/box	CADD™ Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, spiral tubing, clamp and one-way checkvalve with male luer	125 ml/hr	8.0 ml	73"	TOTM (non-DEHP)
	21-7339	12/box	CADD™ Administration Set with bag spike, flow stop, yellow-striped tubing, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	125 ml/hr	6.2 ml	112"	TOTM (non-DEHP)
	21-7394	12/box	CADD™ Administration Set with bag spike, flow stop, clamp, 0.2-micron air-eliminating filter, and one-way checkvalve with male luer	125 ml/hr	3.3 ml	90"	TOTM (non-DEHP)
	21-7395	12/box	CADD™ Administration Set with female luer, flow stop, clamp, 0.2-micron air-eliminating filter, and one-way checkvalve with male luer	125 ml/hr	3.1 ml	90"	TOTM (non-DEHP)

CADD™ High Volume Administration Sets with Flow Stop Free-Flow Protection <i>For use with the CADD-Prizm® pump only (except in PCA delivery mode)</i>							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-7355	15/box	CADD™ High Volume Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male luer	350 ml/hr	20 ml	103"	TOTM (non-DEHP)
	21-7381	15/box	CADD™ High Volume Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male luer	350 ml/hr	18 ml	101"	TOTM (non-DEHP)
	21-7375	15/box	CADD™ High Volume Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, luer-activated needleless injection site, clamp, and one-way checkvalve with male luer	350 ml/hr	18 ml	101"	TOTM (non-DEHP)
	21-7357	15/box	CADD™ High Volume Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	350 ml/hr	11 ml	75"	TOTM (non-DEHP)

Miscellaneous							
	Reorder No.	Units	Description	Use For Rates Up To:	Priming Volume	Output Tubing Length	Plasticizer
	21-1200	12/box	Bag spike set with male luer lock and clamp	-----	0.5 ml	12"	-----
	21-7450	10/box	3 ml syringe-type micro medication reservoir with latex-free plunger (for use with CADD-MS™ 3 pump only)	-----	-----	-----	-----

*subject to availability

smiths

Smiths Medical MD, Inc.
Medication Delivery and
Patient Monitoring Division
St. Paul, MN 55112, USA
www.smiths-medical.com

CADD, CADD-Legacy, CADD-Prizm, and the Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
©2006 Smiths Medical family of companies. All rights reserved. 7/06 19830

Customer and Clinical Services:

1-800-426-2448

www.smiths-medical.com