

smiths

Portex® RapID™ Spinal Needles

- Wide range of needle sizes
- Clear needle hub designed to allow for easy identification of spinal fluid flashback
- Stainless steel needles and stylets
- Color-coded stylet caps
- Offered in lancet point and pencil point tip configurations
- Available as individually packaged sterile needles or in standard and custom spinal anesthesia trays



The Point of Spinal Anesthesia™

Portex® RapID™ Spinal Anesthesia Needles

Ordering Information

Portex® Lancet Point Spinal Needle



Portex® Lancet Point Spinal Needle		Introducer Needle	
49261-20	20g x 1 1/2" (yellow hub)	None	10 Units Per Case
49262-22	22g x 2 1/2" (black hub)	None	10 Units Per Case
4926-22	22g x 3 1/2" (black hub)	None	10 Units Per Case
49263-25	25g x 3" (blue hub)	20g x 1 1/2"	10 Units Per Case
4926-25	25g x 3 1/2" (blue hub)	20g x 1 1/2"	10 Units Per Case
4926-26	26g x 3 1/2" (gold hub)	20g x 1 1/2"	10 Units Per Case
4926-27	27g x 3 1/2" (smoke hub)	20g x 1 1/2"	10 Units Per Case

Portex® Pencil Point Spinal Needle



Portex® Pencil Point Spinal Needle		Introducer Needle	
4924-22	22g x 3 1/2" (black hub)	18g x 1 1/2"	10 Units Per Case
4924-25	25g x 3 1/2" (blue hub)	20g x 1 1/2"	10 Units Per Case
4924-27	27g x 3 1/2" (smoke hub)	20g x 1 1/2"	10 Units Per Case

Please refer to Instructions for Use supplied with the product for detailed instructions, warnings, and cautions. Federal law restricts this device to sale by, or on the order of a physician.

For a complete listing of standard spinal trays containing Portex® RapID™ spinal anesthesia needles, see the Smiths Medical Pain Management Products catalog.

smiths

Smiths Medical MD, Inc.
Medication Delivery Division
 1265 Grey Fox Road
 St. Paul, MN 55112

800-258-5361 603-352-3812 Fax: 603-352-3703
 www.smiths-medical.com

©2004 Smiths Medical family of companies. All rights reserved. Portex® and RapID™ are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

LT-2165 8/04

