

smiths

**LockIt Plus<sup>™</sup>**  
Catheter Securement Device




Making life less painful

“Seventeen percent of catheters fail because they are prematurely dislodged....”

ASRA, 1999 Gaston Labat Lecture, *Acute Pain: Lessons Learned From 25,000 Patients*, L Brian Ready, MD, FRCPC

## LockIt Plus™ Catheter Securement Device

The LockIt Plus™ device is a catheter fixation product from Smiths Medical. It is designed to minimize catheter migration. This helps ensure that pain control delivered via a catheter is not compromised, providing increased confidence and peace of mind. The LockIt Plus™ device has a low profile to help ensure increased patient comfort, and a large adhesive pad to optimize adherence to the skin.



Open LockIt Plus™  
securement device



Insert catheter  
through device



Close device  
to secure catheter

- > Low profile transparent locking mechanism maintains catheter patency and permits monitoring of catheter exit site
- > Large foam adhesive pad makes application easy and minimizes patient discomfort
- > Medical grade adhesive helps decrease skin irritation
- > Easy to peel backing makes application straightforward, even when wearing gloves
- > Available for catheters used with 16/17 and 18 gauge needles

### Ordering Information

Description	Product Code	
LockIt Plus™ Catheter Securement Device	100/399/216	for catheters used with a 16/17G needle
LockIt Plus™ Catheter Securement Device	100/399/218	for catheters used with an 18G needle

**Customer and Clinical Services**  
**1-800-258-5361**

[www.smiths-medical.com](http://www.smiths-medical.com)

**smiths**

Smiths Medical MD, Inc.  
Medication Delivery and  
Patient Monitoring Division  
St. Paul, MN 55112, USA  
[www.smiths-medical.com](http://www.smiths-medical.com)

LockIt Plus and the Smiths design mark is a trademark of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2006 Smiths Medical family of companies. All rights reserved. 8/06 19354