

SmartSite®

Needle-Free Systems

Smiths Medical offers a full line of Needle-Free I.V. disposables to support your continuum of care.



One of the most widely used needle-free systems in the world. Its ergonomic design makes the valve easy to handle, intuitive and simple to use.

SmartSite® Valves

To Add To Your Systems

Extension Sets

With SmartSite® Valves

Bifurcated Sets

And T-Connectors
With SmartSite®

IV Sets

With SmartSite®

Customized

Specialty Sets
With SmartSite®

1









2

3

4

5

SmartSite®

catalog number	description	product status	approx. priming volume	approx. overall length	units/case
SM5000	SmartSite® needle-free valve with male Luer lock (Non-PVC). 	Non-DEHP Formulation Needle-Free Access Site	0.1 ml	1" (3 cm)	100
SM5002	Small bore extension set with SmartSite® injection site, non-removable slide clamp, male Luer lock. 	Non-DEHP Formulation Needle-Free Access Site	0.3 ml	8" (20 cm)	100
SM5008	Small bore T-connector extension set with SmartSite® injection site, non-removable slide clamp, rotating Luer lock. 	Non-DEHP Formulation Needle-Free Access Site	0.2 ml	5" (13 cm)	100
SM5020	Standard bore extension set with 0.2 micron filter, 1 SmartSite® injection site and pinch clamp, swivel luer lock. 	Non-DEHP Formulation Needle-Free Access Site	3.1 ml	20" (51 cm)	25
SM873	20 drop vented/non-vented spike, gravity administration set with back check valve and 3 SmartSite® injection sites, swivel Luer lock. 	Needle-Free Access Site	19 ml	106" (270 cm)	25
2200E	SmartSite® vial access pin. 	Non-DEHP Formulation Needle-Free Access Site			50
2205E	SmartSite® vial adaptor. 	Non-DEHP Formulation Needle-Free Access Site			100
2309E	SmartSite® bag spike with back check valve. 	Non-DEHP Formulation Needle-Free Access Site			100

SmartSite® Valve Specifications

Features

Smooth, swabable surface

Approximate Priming Volume

0.1 ml

Product Status

Non-DEHP Formulation
Straight Through Fluid Path
Latex Free
Lipid Resistant

SmartSite® Needle-Free Systems

Smiths Medical offers a line of access pins, vial adapters and spikes with SmartSite® valves.

If you are looking for an additional set configuration, please contact your Smiths Medical sales representative or call customer service department at 800-848-1757 (Prompt 1). The SmartSite® Valve is also available in custom designed set configurations. Thank you for your interest in Smiths Medical products.

Smiths Medical
2231 Rutherford Road
Carlsbad, CA 92008
USA

T: 800 848 1757
F: 800 621 2590

©2005 Smiths Medical family of companies. All rights reserved. The ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the names, trademarks or service marks of their respective owners. SmartSite is a trademark of Alaris Medical, Inc.

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.

smiths

Smiths Medical