

Portex® Silikon Larynx Maske Because we listened to you



AIRWAY MANAGEMENT


PORTEX™

Die neue Portex® Silikon Larynx Maske

- Material & Design gewährleisten ein Höchstmaß an Komfort, Funktionalität und Sicherheit
- Zeitgemäß – einzeln und steril verpackte Einwegprodukte werden hygienischen und ökonomischen Ansprüchen gerecht
- Erhältlich in 5 Größen für Kinder und Erwachsene

Because we listened to you...

- Einmalprodukt
- Verbesserte Materialeigenschaften (Silikon)
- Leicht einzuführen
- Atraumatisches Design
- Transparentes Material



BESTELLINFORMATION

Silikon Larynx Maske

Bestellnr.	Beschreibung	Größe	I.D. Tubus	A.D. Tubus	Max. Cuff Volumen in ml
100/222/200	Kleinkinder/Kinder	10-20 kg	2,0	7,0	10
100/222/250	Kinder	20-30 kg	2,5	9,0	14
100/222/300	Kinder	30-50 kg	3,0	10,2	20
100/222/400	Erwachsene	50-70 kg	4,0	10,2	30
100/222/500	Erwachsene	70-100 kg	5,0	11,5	40



THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical Deutschland GmbH Bretonischer Ring 3, 85630 Grasbrunn, Telefon +49 89 242959-0, Telefax +49 89 242959-100

Smiths Medical Österreich GmbH

Campus 21, Europaring A 03 5 02, 2345 Brunn am Gebirge, Telefon +43 1 86670 21560, Telefax +43 1 86670 21566

Smiths Medical Schweiz AG Ringwiesenstrasse 2, 8600 Dübendorf, Telefon +41 43 388 62 00, Telefax +41 43 388 62 11

www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Portex and the Smiths Medical Design Mark are trademarks of the Smiths Medical family of companies.
© 2009 Smiths Medical family of companies. All rights reserved. 06/2009

smiths medical