

Precise Effective Wallace Amniocentesis Needles

Sterile, echo marked, disposable needles for sampling amniotic fluid.



Designed for visible sampling in prenatal diagnosis:

- 1cm of spiral engraving on the distal tip provides a reflective surface for enhanced visualisation under ultrasound guidance
- Self-locating stylets to prevent blockage
- Colour coded hubs for easy identification of needle length and gauge
- Range of sizes: 90mm, 120mm, 150mm

Ordering information

Description	Quantity	Order Code
Amniocentesis Needles		
Amniocentesis needle 120mm 21g	10	12021
Amniocentesis needle 150mm 21g	10	15021
Amniocentesis needle 90mm 22g	10	9022

