



Pneupac 880

'The' MRI anaesthetic machine

This compact, mobile unit has been specifically designed to be compatible within the strong magnetic field of an MRI suite. The 880, with its small footprint, is ideal for use where the anaesthetist has limited access to the patient but patient care must not stop.



smiths



Other features

Pipelines for O₂, N₂O and Air (where fitted) Supply and Vacuum: 4 metres O₂ and N₂O high pressure anti-static colour coded hoses to BS EN739 fitted with NIST Connectors at machine end. (Pipeline inlets and cylinders are dependent on specification of machine).

Dimensions: 1500 H x 560 W x 530 D mm **Weight:** 60kg

O₂ and Air Auxiliary Outlet: (To BS4272 pt3) Dependent on specification

Please note: This machine is designed to work in the MRI magnet room within one metre of the ISO centre of the magnet. It is advised that a survey of the MRI room is carried out to address any position problem that could arise regarding the installation.

Accessories

For a full range of MRI compatible accessories please refer to our **Anaesthesia Accessories** leaflet

THE DETAILS IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY, HOWEVER, RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN

For further information, please call your local PneuPac distributor or PneuPac direct on +44 (0)1582 430000

PneuPac

PneuPac Limited

Bramingham Business Park, Enterprise Way, Luton, Beds, LU3 4BU, UK
Tel: +44 (0)1582 430000, Facsimile: +44 (0)1582 430001
<http://www.pneuPac.co.uk> E mail: info@pneuPac.co.uk

Smiths and PneuPac are registered trademarks of Smiths Group plc. Selectatec is a registered trademark of BOC/Ohmeda
Literature No: 504-3055 Issue 3



smiths

Smiths Medical - A part of Smiths Group plc