

EzPAP[®]

Positive Airway Pressure System

Catalog Number	23-0747	System w/mouthpiece	(10/cs)
	23-0757	EzPAP System w/Disposable Manometer	(10/cs)
	23-1747	System w/pediatric mask	(1/cs)
	23-2747	System w/medium mask	(1/cs)
	23-3747	System w/large mask	(1/cs)
	23-6000	EzPAP Kit	(1/cs)

Intended Use: EzPAP is intended for use in conjunction with a medical need for lung expansion therapy and to prevent or reverse atelectasis.

FEATURES	BENEFITS
<ul style="list-style-type: none"> • Use with Mouthpiece or Mask attachment. • 22 OD patient end 	<ul style="list-style-type: none"> • May be used with mask or mouthpiece • Nebulizer
<ul style="list-style-type: none"> • Scalloped Ambient Air Inlet 	<ul style="list-style-type: none"> • Prevents inadvertent occlusion.
<ul style="list-style-type: none"> • Pressure Port with Cap 	<ul style="list-style-type: none"> • Permits connection to gauge tubing.
<ul style="list-style-type: none"> • Complete procedural kit 	<ul style="list-style-type: none"> • Eliminates collecting components.
<ul style="list-style-type: none"> • Disposable manometer 	<ul style="list-style-type: none"> • Eliminates the need to carry a manometer.