

Flo-MistTM
CONTINUOUS MEDICATION NEBULIZER

DHD 11-4100 Flo-mist Continuous Medication Nebulizer (CMN)

11-4110 Flo-Mist CMN, 72" popple tubing and adult aerosol mask

11-4120 Flo-Mist CMN, 72" popple tubing and pediatric aerosol mask

11-4130 Flo-Mist CMN, 72" popple tubing and adult non-rebreathing mask, valved Y piece with plain reservoir bag

11-4140 Flo-Mist CMN, 72" popple tubing and pediatric non-rebreathing mask, valved Y piece with plain reservoir bag

Packaging: 11-4100 = 25/case

11-4110, 11-4120, 11-4130, 11-4140 = 10/case

Intended Use: The Flo-Mist Continuous Medication Nebulizer is designed to deliver up to 8 hours of high volumes of aerosolized medication and humidification therapy.

Feature	Benefit
Auxiliary port	Allows the clinician to add an auxiliary gas flow without recalculating the medication.
Rescue Valve	Located in the non-rebreather circuit, a fail-safe one-way valve in the valved Y piece that allows the patient to entrain ambient air should their inspiratory demand exceed the system.
Connects directly to flowmeter	Once connected to the flowmeter, Flo-Mist is fully secured and easy to access for refilling.
Up to 8 hours of continuous therapy	Flo-Mist delivers medications as well as humidification for any prescribed length of time up to 8 hours.
Dosing chart included	Convenient dosing chart included in each package for quick reference.
Aerosol and non-rebreather kits	A fast and convenient way to administer therapy in emergency care settings.
Standard 22mm aerosol port	Use your existing tubing or purchase Flo-Mist as a kit.
Colorful graphics and gradients	Easy to read gradient lines calibrated at 5ml, ensuring precise medication delivery.