

Minimize hypothermia and maintain comfort when delivering high volume irrigation.



The Level 1[®] NormoFlo[®] Irrigation System. Deliver body temperature irrigation at optimum flow rates.

The Level 1[®] NormoFlo[®] Irrigation System is perfect for a variety of procedures including:

- Laparoscopy
- Arthroscopy
- Hysteroscopy
- Wound Debridement
- Urology
- Liposuction
- Lavage



Model H-1129 shown.



Pressure system accommodates standard 1-liter bags in rigid chambers and also includes 3-liter pressure cuffs (System H-1129 only).

Reduce bleeding, speed recovery times and lower infection rates.

With the Level 1® NormoFlo® Irrigation System you can deliver large volumes of body temperature irrigation at optimum flow rates. This innovative system helps you effectively maintain your patients' core temperature during the entire length of the procedure. Minimizing hypothermia has been shown to reduce bleeding, speed recovery and lower infection rates.

No change in surgical technique.

The Level 1® NormoFlo® Irrigation System delivers a clear, bubble-free field of vision, which allows surgeons to operate with confidence and with **no** change in surgical technique. It easily interfaces with standard surgical tools, hand pieces, scopes and specialized pumps to provide consistently safe performance that improves patient comfort and delivers better outcomes.

Easy to set up and use.

There are two versions of the Level 1® NormoFlo® Irrigation System:

Model H-1100 is intended for gravity flow procedures and features an electrically operated I.V. pole that can automatically lift up to four 3-liter bags at the touch of a button.

Model H-1129 has an electrically operated I.V. pole and also includes an automated pressure system that comes with 1-liter rigid chambers (maximum pressure; 500mmHg) and 3-liter cuffs (maximum pressure; 300mmHg).

Both systems are easy to use and require minimal training time.



Lift heavy bags at the touch of a button.

Ordering Information

HARDWARE

- H-1100** Level 1® NormoFlo® Irrigation Fluid Warmer includes electrically operated I.V. pole
- H-1129** Level 1® NormoFlo® Irrigation Fluid Warmer includes electrically operated I.V. pole and pressure system

DISPOSABLES

- IR-600** Irrigating set, for all surgical irrigating, gravity flow, or can be pressurized up to 300mm/Hg. 10/case
- IR-500** Irrigating set, for all surgical irrigating, gravity flow, or can be pressurized up to 300mm/Hg. 10/case
- IR-700** Irrigating set, for all surgical irrigating, gravity flow, or can be pressurized up to 300mm/Hg. 9/case
- IR-40** Irrigating set, can be pressurized up to 500mm/Hg. 10/case

DISPOSABLE ACCESSORIES

- IR-6** High Flow "Y" Adapter, converts IR-40, IR-500, and IR-600 into 4 spike configuration. 10/case
- IR-7** High Flow "Y" Adapter, same as IR-6 except incorporates vented spikes. For use with rigid plastic bottles. 10/case

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS.

Smiths Medical ASD, Inc.

Smiths Medical ASD
160 Weymouth Street
Rockland, MA 02370
Phone: 800 848 1757 Fax: 800 621 2590
www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Smiths Medical Design Mark, the Level 1 Design Mark, and NormoFlo are registered trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. Copyright 2007 Smiths Medical family of companies. All rights reserved.

TM407 2007.9