

## Giving you the EDGE on needle safety

Hypodermic Needle-Pro® EDGE™ Safety Device



SAFETY ACCESS



**A new standard in safety**

A leader in the design of innovative needle safety devices that are easy to use and provide cost-effective solutions to needle safety requirements.

In a fast-paced, clinical environment like yours, patient and clinician safety can be easily overlooked. That's why today's strict hospital and federal directives demand accuracy, efficiency and cost-effective safety measures.


Smiths Medical has responded to these requirements by incorporating innovative safety features into our Hypodermic Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device. This unique device gives clinicians an important edge when it comes to safety, convenience and cost-efficiency.

The future of hypodermic needle safety is here. Get the Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device today!



**Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device**

- The Hypodermic Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device has a color-coded safety sheath for rapid identification of needle gauge size
- Available with 18g, 19g, 20g, 21g, 22g, 23g, 25g, 26g, 27g and 30g needles
- Adds no dead space as compared to a traditional needle hub. The needle is fully integrated into the safety device
- Bevel is fixed in the "up" position with the safety sheath to the right. The arrow on the sheath points in the direction of the bevel position
- Simple, one-handed operation of the safety device



Syringe markings are compliant with Joint Commission's Official "Do Not Use" list regarding use of decimal points.

Easy-to-read graduation markings

One-piece design completely integrates safety device

Bevel-oriented needle

Safety sheath is color coded for rapid identification of needle gauge size

**Safety and compliance made easy**



Perform a one-handed technique



Close safety device on any flat surface



Discard in sharps container

**With the Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device, you can make injections more convenient with a pre-attached syringe.**

Features of the Needle-Pro<sup>®</sup> EDGE<sup>™</sup> Safety Device include:

- Easy-to-read graduation markings on the syringe
- Simple, one-handed activation of the safety feature helps you stay safe and compliant
- Helps reduce potential for needlestick injuries by minimizing critical exposure time



Molded arrow points in direction of bevel.

## Hypodermic Needle-Pro® EDGE™ Safety Device

Reference Code	Description	Color-Coded	Units Per Case
401810	18g x 1" Hypodermic Needle-Pro® EDGE™	Pink	1000
401815	18g x 1 ½" Hypodermic Needle-Pro® EDGE™	Pink	1000
401910	19g x 1" Hypodermic Needle-Pro® EDGE™	Cream	1000
401915	19g x 1 ½" Hypodermic Needle-Pro® EDGE™	Cream	1000
402010	20g x 1" Hypodermic Needle-Pro® EDGE™	Yellow	1000
402015	20g x 1 ½" Hypodermic Needle-Pro® EDGE™	Yellow	1000
402110	21g x 1" Hypodermic Needle-Pro® EDGE™	Green	1000
402115	21g x 1 ½" Hypodermic Needle-Pro® EDGE™	Green	1000
402210	22g x 1" Hypodermic Needle-Pro® EDGE™	Black	1000
402215	22g x 1 ½" Hypodermic Needle-Pro® EDGE™	Black	1000
402310	23g x 1" Hypodermic Needle-Pro® EDGE™	Blue	1000
402315	23g x 1 ½" Hypodermic Needle-Pro® EDGE™	Blue	1000
402510	25g x 1" Hypodermic Needle-Pro® EDGE™	Orange	1000
402515	25g x 1 ½" Hypodermic Needle-Pro® EDGE™	Orange	1000
402558	25g x ⅝" Hypodermic Needle-Pro® EDGE™	Orange	1000
402605	26g x ½" Hypodermic Needle-Pro® EDGE™	Brown	1000
402705	27g x ½" Hypodermic Needle-Pro® EDGE™	Gray	1000
4027125	27g x ¼" Hypodermic Needle-Pro® EDGE™	Gray	1000
403005	30g x ½" Hypodermic Needle-Pro® EDGE™	Yellow	1000

## Hypodermic Needle-Pro® EDGE™ Safety Device with 3mL Luer lock syringe

431810	18g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Pink	400
431815	18g x 1 ½" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Pink	400
431910	19g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Cream	400
431915	19g x 1 ½" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Cream	400
432010	20g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Yellow	400
432015	20g x 1 ½" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Yellow	400
432110	21g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Green	400
432115	21g x 1 ½" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Green	400
432210	22g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Black	400
432215	22g x 1 ½" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Black	400
432310	23g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Blue	400
432510	25g x 1" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Orange	400
432558	25g x ⅝" Hypodermic Needle-Pro® EDGE™ with 3 mL Luer lock syringe	Orange	400

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

### Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370 USA

T: 800 848 1757 F: 800 621 2590 www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Smiths Medical, JELCO design marks, Needle-Pro and EDGE are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © 2008 Smiths Medical family of companies. All rights reserved. SS400 2.2008