

TheraPEP[®]

Positive Expiratory Pressure System

Catalog Number	20-1112 TheraPEP system w/mouthpiece	(10/cs)
	20-0010 pressure port	(10/cs)
	20-0022 pressure port, tube & indicator	(10/cs)
	20-0120 valved resistor, pressure port	(10/cs)
	20-1110 valved resistor, pressure port, mouthpiece	(10/cs)
	20-3112 system w/pediatric mask	(1/bx)
	20-5112 system w/small mask	(1/bx)
	20-7112 system w/large mask	(1/bx)
	20-3115 pediatric mask	(10/bx)
	20-5115 small mask	(10/bx)
	20-7115 large mask	(10/bx)

Intended Use: Using TheraPEP, patients can master Positive Expiratory Pressure therapy quickly and self-administer effective secretion removal away from the hospital. TheraPEP is a great alternative for patients who cannot tolerate CPT and those suffering from lung conditions with secretion problems such as, COPD including: (chronic bronchitis & CF) pneumonia, post-operative atelectasis and asthma.

FEATURES	BENEFITS
<ul style="list-style-type: none"> • Six Fixed Orifice Options 	<ul style="list-style-type: none"> • Can accommodate virtually any patient's lung capacity.
<ul style="list-style-type: none"> • Built-in durable pressure indicator 	<ul style="list-style-type: none"> • Resists breakage, unlike fragile, costly manometers. • Provides immediate, visual 360° feedback of prescribed pressure.
<ul style="list-style-type: none"> • 22mm OD ambient end • 22mm OD patient end 	<ul style="list-style-type: none"> • May be used with a mask or mouthpiece, or • Nebulizer
<ul style="list-style-type: none"> • Inspiratory valve 	<ul style="list-style-type: none"> • Allows inhalation and exhalation without removing from mouth.