

# Security Freedom

## Wallace Ring Pessaries



Designed to ease uterine prolapse and other uterine displacements.

### Range of sizes to suit all patients

- Flexible PVC with large cross section
- Sixteen sizes from 50mm to 110mm (outside diameter)
- Individually packed in easy to open peel pouches

### Simple to Fit

- Ring softens when immersed in warm water to aid fitting in a figure of eight configuration

### Easy to Maintain

- Remove and check ring, cervix and vaginal cavity after three months
- Wash ring in soapy water before re-insertion according to instructions for use
- Replace with a new ring after six months

### Ordering information

Description	Quantity	Order Code
<b>Ring Pessaries</b>		
PVC Ring	10	700/300/xxx*
Sizes available (mm) 50, 53, 56, 59, 62, 65, 68, 71, 74, 77, 80, 85, 90, 95, 100, 110		
Packed in individual peel pouches in cartons of 10		
*xxx is the size of the individual ring preceded by 0 if only two digits		

