

Soft Visible

Wallace Sure View® Ultrasound Visible Embryo Replacement Catheters



Sure View® technology - allowing for the whole of the catheter to be viewed under ultrasound ¹¹

Sure View® Catheter under ultrasound in phantom uterus

Competitor ultrasound visible tip

Cross-section of the Sure View® catheter illustrating smooth inner lumen

Sure View® Embryo Replacement Catheter CE118 [18 cm]

Sure View® Embryo Replacement Catheter CE123 [23 cm]

Sure View® Trial Transfer Catheter CE418 [18 cm]

Sure View® Trial Transfer Catheter CE423 [23 cm]

More and more clinicians and embryologists worldwide have witnessed demonstrable improvement in pregnancy rates when they use ultrasound during the embryo replacement procedure.⁷⁻¹³

In cooperation with clinical experts, Smiths Medical International, the designer and manufacturer of the gold standard Wallace embryo replacement catheter, has developed the Sure View® range to improve ultrasound guided embryo transfer. Sure View® combines all of the benefits of the classic Wallace catheter with the unique opportunity to view the entire length of the catheter under ultrasound. The catheter is soft, visible and requires no compromise.

KEY PRODUCT BENEFITS

Visible

- Visible under ultrasound from hub to tip
- Highly visible in simulated uterine environment
- Increased confidence of direction and position of the uterus
- Optimal placement of catheter in the uterus
- More accurate positioning for embryo expulsion

No Compromise

- Combines the benefits of the classic range with state of the art technology
- Exploits the known benefits of ultrasound with the Wallace gold standard catheter
- Compatible with the Wallace malleable stylet for difficult transfers
- Unique ability to view the entire catheter under ultrasound

ORDERING INFORMATION AND COMPATIBLE PRODUCTS

Description	Rigidity	Quantity	Order Code
Sure View® Embryo Replacement Catheter (18cm)	soft	10	CE118
Sure View® Embryo Replacement Catheter (23cm)	soft	10	CE123
Sure View® Trial Transfer Catheter blind end (18cm)	soft	10	CE418
Sure View® Trial Transfer Catheter blind end (23cm)	soft	10	CE423
Stylet for introducing 18cm catheter	rigid	10	1816ST
Stylet for introducing 23cm catheter	rigid	10	1816NST

The Sure View® range is fully compatible with the Classic range.