



Dual Sharp


Wallace Dual Lumen Oocyte Recovery Systems




Available in 16g and 17g with colour coded sleeves



Echo markings right to the tip




Flushing lumen extends to needle tip



Ergonomic handle with orientation marker

Wallace Dual Lumen Oocyte Recovery System

Cross Section of Dual Lumen Needle



Needles

Smiths Medical, designer and manufacturer of the worldwide gold standard Wallace Embryo Replacement Catheter, have used their knowledge and expertise to introduce a range of premium quality, Dual Lumen Oocyte Recovery products.

The Wallace Dual Lumen Needles have been carefully designed in order to maximise the number of oocytes retrieved, whilst minimising operating time and patient discomfort.

KEY PRODUCT BENEFITS

- Echo marked to tip for safer positioning of needle within follicle
- Flushing lumen extends to needle tip to maximise flow rate through dual lumen needle
- Ergonomic handle with orientation marker, allowing easy and controlled handling
- Silicone bung for secure and safe use in test tubes
- 16g and 17g needle size with colour-coded tubing sleeves

ORDER INFORMATION

Description	Quantity	Order Code
Dual lumen 16g needle (33cm) with 950mm aspiration	10	DNS1633-950
Dual lumen 16g needle (33cm) with 750mm aspiration	10	DNS1633-750
Dual lumen 16g needle (33cm) with 500mm aspiration	10	DNS1633-500
Dual lumen 17g needle (33cm) with 950mm aspiration	10	DNS1733-950
Dual lumen 17g needle (33cm) with 750mm aspiration	10	DNS1733-750
Dual lumen 17g needle (33cm) with 500mm aspiration	10	DNS1733-500