

Clear-Cuff® Pressure Infusor



Economical | Durable | Convenient



Clear-Cuff® Pressure Infusor




Durable

Durable, easy-to-clean, polyurethane construction allows the Clear-Cuff Pressure Infusor to be re-used if desired.*



Circumferential bladder construction

The Clear-Cuff Pressure Infusor provides a consistent flow rate for effective fluid management. Its circumferential bladder construction allows for an optimal squeeze of the entire fluid bag, enabling more consistent flow rates.



Beneficial features

The Clear-Cuff Pressure Infusor, made of thick, durable polyurethane, has an inherent transparent quality, providing easy viewing of the fluid bag and the fluid level. Additionally, a special built-in hook allows the fluid bag to be replaced without removing the Clear-Cuff Pressure Infusor from the IV pole. The Clear-Cuff Pressure Infusor also features an extra-large squeeze bulb and pressure gauge for fast and easy inflation.



Product Code	Description	Qty/Case
MX4705	Clear-Cuff 500ml Pressure Infusor	10
MX4710	Clear-Cuff 1,000ml Pressure Infusor	10

*Cleaning the Clear-Cuff Pressure Infusor should be done according to the Instructions for Use.

Medex, Inc.

2231 Rutherford Road
Carlsbad, CA 92008
USA

Medex Holding GmbH

Prinzenallee 15
40549 Düsseldorf
Germany



info@medex.com
www.medex.com

(800) 848 1757 phone
(800) 621 2590 fax

49 (0) 211/52773-0 phone
49 (0) 211/52773-26 fax

© 2005 Medex, Inc. The trademarks and names herein are the property of Medex, Inc.