



**For infusion therapy, Smiths Medical
knows you have a lot of shoes to fill**



**A single large volume pump platform
may not be the best fit for all your patient needs**

smiths medical

bringing technology to life

Smiths Medical has the systems you need to deliver specialized infusions



Two highly-targeted infusion systems with a lot in common

- Small, lightweight, portable pumps
- Medication safety software for dose error reduction
- Configurable for multiple care areas
- The right reservoir options for flexible, efficient drug compounding and filling
- Accurate, reliable medication delivery from specialized pumps – not add-on modules

Smart. Really Smart.

Our medication safety software does much more than just dose error reduction. It enables you to optimize your infusion services from prescription to delivery.

CADD[®]-Solis

Ambulatory Infusion System with
CADD[®]-Solis Medication Safety Software



*Right in-step
with your pain
management
infusion needs*

Smart Medication Safety

- Meets industry-recognized standards for advanced dose error reduction in PCA pumps
- Treatment-based medication safety software reflects best practices
- Facility-defined therapy protocols promote safety and standardization
- On-screen graphs, trend reports and user audit trail

Functional and Flexible

- Easy to load, program, and start for busy pain management teams
- Wide variety of CADD[®] administration sets, featuring 50, 100 and 250 mL medication cassette reservoirs

Medfusion[™] 3500

Syringe Infusion System with
PharmGuard[®] Safety Software



*Keeping
pace with
your syringe
infusion
requirements*

PharmGuard[®] Safety Software

- Designed to input more than 16 profile or care areas
- Dose/kg/time for more accurate drug delivery
- QuickLibrary[™] alternative for fast-response anesthesia programming
- Over 1200 possible profiles storable in pump memory

Therapy Ready without “Add-ons”

- Highly accurate flow rates from 0.01 mL/hr to 1130 mL/hr
- Wide syringe compatibility from 1 mL to 60 mL syringe size
- Easy, simple-to-load syringe system
- Integral rapid occlusion detection and post-occlusion bolus reduction

safety...portability...scalability...connectivity

Large volume pumps serve their purpose for many infusion needs, but your infusion teams serve critical niches that only a highly-targeted, focused infusion solution can effectively deliver. That's where CADD®-Solis and Medfusion™ 3500 infusion systems from Smiths Medical fit in your infusion practice.

When it's time to evaluate infusion systems, consider the full spectrum of infusion needs

Smiths Medical provides a perfect fit for nursing, pharmacy, bio-med, materials management – all the shoes your team needs to fill to meet your patients' specialty infusion needs.

Find out how we can fit into your hospital infusion plan

contact Smiths Medical today

call..... 1-800-426-2448

email bestfit@smiths-medical.com

visit www.smiths-medical.com/bestfit

Smiths Medical infusion systems

the best fit

for all the shoes your team needs to fill
to meet your patients' specialty infusion needs

www.smiths-medical.com/bestfit

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN

Smiths Medical MD, Inc.

St. Paul, MN 55112, USA

Phone: 1-800-426-2448

www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Customer and Clinical Services 1-800-426-2448
www.smiths-medical.com

medfusion™

CADD™

CADD-Solis, CADD, Medfusion, PharmGuard and the Smiths Medical design mark are trademarks of the Smiths Medical family of companies.
The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
©2008 Smiths Medical family of companies. All rights reserved. 11/08 193700

smiths medical