

NEW!

smiths medical
bringing technology to life

Medfusion® 3500 Syringe Pump with PharmGuard® v5 Software

Medfusion® 3500
syringe pump
now recognizes

**Enteral
Syringes**



**For safe, accurate, enteral
nutrition delivery**

medfusion

™

Reduce the risk of IV / enteral misconnections

Medfusion® 3500 syringe pumps with PharmGuard® v5 software now recognize NeoMed®, Baxa® NeoThrive™, and Vygon® enteral syringes

Adding enteral syringe data to PharmGuard® v5 software allows you to:

- Meet Joint Commission requirements for enteral feedings¹
- Differentiate your patients' IV and enteral lines²
- Meet your facility's needs and provide smarter, safer enteral delivery

Supported Enteral Syringes

NeoMed® - 2 sizes - 20 mL, 60 mL - orange
Baxa® NeoThrive™ - 3 sizes - 20 mL, 35 mL, 60 mL - orange
Vygon® C-Gon - 2 sizes - 20 mL, 60 mL - purple

Smart and safe just got smarter and safer

By partnering with the leading enteral syringe companies in the market, Smiths Medical ensures you'll have the best of both worlds - a smarter infusion system and safer enteral delivery

- Orange or purple coloring for easy enteral identification
- Unique tips will not connect to standard I.V. Luers



neothrive
Brought to you by Baxa Corporation.

VYCON

NEOMED
incorporated

Call Smiths Medical today to find out how to incorporate enteral syringe data into your Medfusion® 3500 pumps with PharmGuard® v5 software

1-800-258-5361

Order Information

PharmGuard® Supported Syringes Series 1 67-2449-51

1. Guenter, Peggi, PhD., R.N., C.N.S.N., et al. *Enteral Feeding Misconnections: A Consortium Position Statement*. The Joint Commission Journal on Quality and Patient Safety. May 2008; Volume 34 Number 5
2. Abstract: *Enteral feeding set connectors and adapters*. Association for the Advancement of Medical Instrumentation. ANSI/AAMI ID54:1996/©2005.

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

St. Paul, MN 55112, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361

www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Product(s) described may not be licensed or available for sale in Canada or other countries outside of the United States.

Medfusion, PharmGuard and the Smiths Medical and Medfusion design marks are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2009 Smiths Medical family of companies. All rights reserved. 10/09 193816

Smiths Medical Canada Ltd.

Markham, Ontario, Canada, L3R 4Y8
Phone: 905-477-2000
Toll-Free: 1-800-387-4346

Customer and Clinical Services

U.S. 1-800-258-5361
Canada 1-800-387-4346

www.smiths-medical.com

smiths medical