

Bivona® Customized Adjustable Tracheostomy Tube Template

smiths medical

bringing technology to life

PORTEX™ NEW easy-to-read template
Just follow steps 1-9. Please specify:

PHONE: 1 800 424 8662 FAX: 219 989 7435

WEB: www.smiths-medical.com



1. Shaft Length:

Shaft Length _____mm (30mm min, max. 200mm)



2. Shaft:

Shaft: I.D. _____mm (2.5-9.5mm) or Shaft: O.D. _____mm



3. Connector:

<input type="checkbox"/> 3a. Swivel 2.5-9.5mm	<input type="checkbox"/> 3b. Fixed 2.5-6.0mm
	

4. Neck Flange:

<input type="checkbox"/> 4a. 2.5-5.5mm	<input type="checkbox"/> 4c. 5.0-9.5mm
	
Neck Flange Length: 60mm	Neck Flange Length: 83mm


5. Cuff:

<input type="checkbox"/> 5a. TTS™ Tight to Shaft	<input type="checkbox"/> 5b. Aire-Cuf® Pedi	<input type="checkbox"/> 5b. Aire-Cuf® Adult	<input type="checkbox"/> 5c. Fome-Cuf®	<input type="checkbox"/> 5d. Cuffless
				
2.5-9.5mm	2.5-5.5mm	5.0-9.5mm	2.5-9.5mm	2.5-9.5mm

6. Specialty Shaft:

6b. Hyperflex

Reinforced- Straight



2.5-9.5mm

Not intended for long term use.
This tube provides temporary airway access for a tracheostomized patient when determining the optimal tube length for the patient. This tube must be replaced with a fixed neck flange tracheostomy tube when the optimal length is determined.

7. Optional:

7a. Distal (Bottom) Cuff

Cuff distance from tip to shaft _____mm

Note: Cuff distance from tip cannot be less than the minimum standard distance listed below:
a. Pediatric TTS, Aire & Fome – 3mm
b. Adult Aire & Fome – 5mm
c. Adult TTS – 9mm

8. Approval: NOTE: Please print clearly. Smiths Medical will review the completed template for compliance. Any required modifications will be made with the approval of a clinician.

HOSPITAL/CLINIC

DATE

P.O. NUMBER

CLINICIAN SIGNATURE/CREDENTIALS

PATIENT/CLIENT ID

(Smiths Medical Use Only) Approval: _____ Reorder: _____

9. Select Type of Service:

3-4 weeks (Sterile) Standard Service (NON-Sterile) Express Service (NON-Sterile) Same Day Service (NON-Sterile) – USA Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

The Smiths Medical and Portex design marks; Bivona, TTS, Aire-Cuf, and Fome-Cuf are trademarks of the Smiths Medical family of companies. The Symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © 2012 Smiths Medical family of companies. All rights reserved.