




Pneupac
Airway-Ventilation Systems

Simply Anaesthesia




smiths

Simply Anaesthesia




This compact, lightweight machine has been designed as the perfect solution for emergency anaesthesia. Ideal for mobile medical units, its balanced portability enables it to be carried anywhere, used free-standing or rail mounted.

- On/Off switch
- Selectatec backbar
- O₂ flush
- Ratio flow device
- Guarded flowmeter valves
- Visual and audible O₂ supply failure alarm
- Pipeline supply gauges
- Handles
- O₂ auxiliary outlet
- O₂ and N₂O inlets
- Rail clamps for mounting




This small, rugged machine is perfectly designed as a mobile anaesthetic unit. Its small footprint allows it to be easily transported between departments and it is ideally suited for operating theatres, casualty departments, field hospitals and anaesthetic rooms to provide you with the exact versatility you need.

- On/Off switch
- O₂ flush
- Selectatec backbar
- Ratio flow device
- Guarded flowmeter valves
- Visual and audible O₂ supply failure alarm
- Pipeline supply gauges
- Cylinder supply gauges
- O₂ and Air auxiliary outlets
- O₂, Air and N₂O inlets
- Drip stand (option)
- Monitor shelf
- Drawer options (single as standard)
- Pull handle (not shown on picture above)




The versatile 550 machine is designed to suit applications within operating theatres, induction rooms and casualty departments. It can be used as a simple or sophisticated machine and therefore be adapted for use in any environment.

- On/Off switch
- O₂ flush
- Twin station selectatec backbar
- Ratio flow device
- Guarded flowmeter valves
- Visual and audible O₂ supply failure alarm
- Pipeline supply gauges
- Cylinder supply gauges
- O₂ and Air auxiliary outlet
- O₂, Air and N₂O inlets
- Machine protection valve
- Drip stand (option)
- Eye level instrument shelf
- Stowage hooks
- Drawer options (single as standard)
- Pull handle (not shown on picture above)



This compact, lightweight wall mounted unit is ideal for use in small anaesthetic rooms, casualty departments, maternity suites and plaster rooms. It has been designed to offer unhindered visibility and perform in tough environments where space may be at a premium.

- On/Off switch
- Selectatec backbar
- O₂ flush
- Ratio flow device
- Guarded flowmeter valves
- Visual and audible O₂ supply failure alarm
- Pipeline supply gauges
- O₂ and Air auxiliary outlet
- O₂, Air and N₂O inlets
- Rail clamps for mounting
- Work surface
- Instrument shelf



This compact, mobile unit has been specifically designed to be compatible within the strong magnetic field of an MRI suite. The 880, with its small footprint, is ideal for use where the anaesthetist has limited access to the patient but patient care must not stop.

- On/Off Switch
- O₂ flush
- Single station Selectatec* backbar
- Ratio flow device
- Guarded flowmeter valves
- Visual and audible O₂ supply failure alarm
- Pipeline supply gauges
- O₂, N₂O & Air cylinder supply gauges
- O₂ and Air auxiliary outlets
- O₂, Air and N₂O inlets
- Machine protection valve
- Eye level instrument shelf
- Lower box unit
- Drip stand
- Vaporiser bay

* Selectatec is a registered trademark of BOC/Ohmeda

THE DETAILS IN THIS LEAFLET GIVEN ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY, HOWEVER, RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN

**For further information, please call your local Pneupac distributor or
Pneupac direct on +44 (0)1582 430000**

Pneupac Limited

Bramingham Business Park, Enterprise Way, Luton, Beds, LU3 4BU, UK
Tel: +44 (0)1582 430000, Facsimile: +44 (0)1582 430001
<http://www.pneupac.co.uk> E mail: info@pneupac.co.uk

Smiths and Pneupac are registered trademarks of Smiths Group plc. Selectatec is a registered trademark of BOC/Ohmeda
Literature No: 504-3048 Issue 2



smiths
Smiths Medical - A part of Smiths Group plc