

## CADD-Sentry Pro™ Medication Safety Software – Administrator and Point of Care Installation Requirement Checklist

The CADD-Sentry Pro™ Medication Safety Software consists of the Administrator and Point of Care software applications. This software allows you to create a set of standard pump Protocols to be used with the CADD-Prizm® PCS II Ambulatory Infusion Pump.

Read through the *Instruction For Use* for the Administrator and Point of Care applications, as well as the help files located on the software programs for detailed instructions, warnings and cautions. If you have questions while installing the software, please contact Clinical Services at 1-800-426-2448.

### Administrator Requirements

#### Computer Requirements:

- Microsoft Windows 2000 Server or Windows 2003 Server (when hosting the Database on a Server)
- A PC running Windows 2000 (with SP4 or newer installed) or XP (with SP2 or newer installed)
- CD-ROM drive
- Minimum 70 MB available hard disk space (application)
- Minimum 100 MB available hard disk space (database)
- Minimum 256 MB RAM
- Mouse or other pointing device
- 16-bit color capable graphics card
- Display device with minimum screen resolution of 1024 X 768

#### CADD™ Pump:

- CADD-Prizm® PCS II Ambulatory Infusion pump (with software revision H or higher)

#### Cable Requirements:

- CADD-Prizm® Interface Cable/Null Modem Cable Set
- Available RS232 serial port for connection to the Interface Cable/Null Modem Cable Set

#### Barcode Reader Requirements, if applicable:

- Barcode reader - configurable as a Human Interface Device (HID) code 128 compatible

#### Printer Requirements, if applicable:

- Zebra brand barcode printer (compatible with ZPL or ZPL II programming language – three inch wide by one inch high label stock on continuous feed roll – USB)

Installation Notes:

---



---



---

## CADD-Sentry Pro™ Medication Safety Software – Administrator and Point of Care Installation Requirement Checklist

The CADD-Sentry Pro™ Medication Safety Software consists of the Administrator and Point of Care software applications. This software allows you to create a set of standard pump Protocols to be used with the CADD-Prizm® PCS II Ambulatory Infusion Pump.

Read through the *Instruction For Use* for the Administrator and Point of Care applications, as well as the help files located on the software programs for detailed instructions, warnings and cautions. If you have questions while installing the software, please contact Clinical Services at 1-800-426-2448.

### Point of Care Requirements

#### Computer Requirements:

- A PC running Windows 2000 (with SP 4 or newer installed) or XP (with SP2 or newer installed) Or
- A Tablet PC running Windows XP Tablet PC Edition 2005 (SP1)
- CD-ROM drive
- Minimum 70 MB available hard disk space
- Minimum 256 MB RAM
- Mouse or other pointing device
- 16-bit color capable graphics card
- Display device with minimum screen resolution of 800 X 600

#### CADD™ Pump:

- CADD-Prizm® PCS II Ambulatory Infusion pump (with software revision H or higher)

#### Cable Requirements:

- CADD-Prizm® Interface Cable/Null Modem Cable Set
- Available RS232 serial port for connection to the Interface Cable/Null Modem Cable Set

#### Barcode Reader Requirement, if applicable:

- Barcode reader - configurable as a Human Interface Device (HID) code 128 compatible

Installation Notes:

---



---



---