


# Level 1® Fast Flow Fluid Warmer Disposable Administration Sets and Accessories

## Level 1® Fast Flow Fluid Warmer Disposable Administration Sets


Smiths Medical offers a full line of Disposable Administration Sets for use with all Level 1® Fast Flow Fluid Warmers. These easy to use Disposable Administration Sets:

- Have a Replaceable Gas Vent/Filter Assembly
- Are Latex Free
- Provide a Sterile Fluid Path
- Are Designed for Single Patient Use




### DI-50

- Priming Volume: 56 mL
- Maximum Normothermic Flow Rate: 300 mL/minute
- Maximum Flow Rate: 580 mL/minute (Crystalloid, 300 mmHg, 14 Gauge catheter)
- Quantity: 20/case



### D-60HL / DI-60HL

- Priming Volume: 74 mL
- Maximum Normothermic Flow Rate: 500 mL/minute
- Maximum Flow Rate: 500 mL/minute (Crystalloid, 300 mmHg, 14 Gauge catheter)
- Triple lumen tube prevents patient line cool down
- For use in unpredictable low- to high-flow surgical procedures
- Quantity: 10/case




### D-70 / DI-70

- Priming Volume: 64 mL
- Maximum Normothermic Flow Rate: 600 mL/minute
- Maximum Flow Rate: 600 mL/minute (Crystalloid, 300 mmHg, 14 Gauge catheter)
- Quantity: 20/case

D-Series Disposable Administration Sets: For use in the USA  
DI-Series Disposable Administration Sets: For use in markets outside of the USA  
Flow rate  $\pm$  10%; 20° C  $\pm$  2° C infusate; and 300 mmHg + 0 / -12  
Approximate Priming Volume  $\pm$  10 mL


### Compatible Disposable Administration Sets

Level 1® Fast Flow Fluid Warmers	D-50 DI-50	D-60HL DI-60HL	D-70	D-100 DI-100	D-300 DI-300
Level 1® H-1200 with Integrated Air Detection	•	•	•	•	•
Level 1® H-1000	•	•	•	•	•



### D-100 / DI-100

- Priming Volume: 65 mL
- Maximum Normothermic Flow Rate: 800 mL/minute
- Maximum Flow Rate: 800 mL/minute (Crystalloid, 300 mmHg, 8.5 Fr catheter)
- For use in trauma cases
- Quantity: 10/case



### D-300 / DI-300

- Priming Volume: 90 mL
- Maximum Normothermic Flow Rate: 1000 mL/minute
- Maximum Flow Rate: 1255 mL/minute (Crystalloid, 300 mmHg, 8.5 Fr catheter)
- For use in massive trauma or liver transplant procedures
- Quantity: 10/case

® Fast Flow  
HOTLINE®

® NORMOFLO®  
Equator Warm®  
Probes®


LEVEL 1

Warmth

## Level 1® Fast Flow Fluid Warmer Accessories


Level 1® Fast Flow Fluid Warmer Accessories:

- Are Latex Free
- Provide a Sterile Fluid Path
- Are Designed for Single Patient Use




### F-10 Replacement Gas Vent/Filter Assembly

- Filter Size: 170 micron
- Replacement Gas Vent/Filter Assembly for DI-50, D-60HL/DI-60HL, D-70/DI-70, DI-100, D-100/DI-100
- Change every three hours, when filter becomes clogged, or if air is venting slowly
- Quantity: 20/case




### F-30 Replacement Gas Vent/Filter Assembly

- Filter Size: 170 micron
- Large Bore Tubing
- Replacement Gas Vent/Filter Assembly for D-300/DI-300
- Change every three hours, when filter becomes clogged, or if air is venting slowly
- Quantity: 20/case




### PF-1 Fast Flow Pre-Filter

- Filter Size: 340 micron
- Spikes directly into blood bag
- Extends life of F-10 and F-30 Filters with Gas Vent/Filter Assembly Replacement
- Traps large particles and clots
- Use one per unit of blood
- Quantity: 60/case




### Y-INJ Fast Flow Extension with Injection Site

- (8 in) 20.3 cm fast flow I.V. extension with a Y-injection site
- Large bore patient lead
- Priming volume: 2.2 mL
- Use as a patient lead or for adding a second I.V. line
- Contains needle or needleless injection site
- Quantity: 40/case




### X-36 Fast Flow Extension Line

- (36 in) 91.4 cm fast flow extension line with a roller clamp
- Priming volume: 8 mL
- Quantity: 40/case



### Y-30 Fast Flow Y-Type Extension

- (30 in) 76.2 cm fast flow I.V. extension line and roller clamp with a 15.2 cm (6 in) Y-Type extension and clamp
- Intended to split fluid output from a Level 1® Disposable Administration Set to two I.V. sites
- Priming volume: 9.5 mL
- Quantity: 40/case



### SC-3 Fast Flow 3-Way Stopcock

- (9 Fr) 3.0 mm inner diameter
- Connects to Level 1® Disposable Administration Set to provide syringe access or additional I.V. administration set access to the patient
- Quantity: 40/case

smiths medical

#### Smiths Medical Canada Ltd.

Markham, Ontario, Canada, L3R 4Y8

Phone: 905-477-2000

Toll-Free: 1-800-387-4346

www.smiths-medical.com

Level 1, HOTLINE, NormoFlo, Snuggle Warm and the Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates it is registered in the U.S. Patent and Trademark Office and certain other countries. Please read the instructions for use supplied with the product for detailed instructions, warnings and precautions.

©2010 Smiths Medical family of companies. All rights reserved. TM194804 2010.7