

Safety and accuracy when you need it most

Level 1[®] Neo-Therm[®] Skin Temperature Probes



TEMPERATURE MANAGEMENT

Neo-Therm® Skin Temperature Probes

Neo-Therm® skin temperature probes offer exceptional safety, quality, and reliability. The sensor consists of a 1"-diameter, hypoallergenic, adhesive foam disk with a Mylar® cover.

Neo-Therm® Connecting Cables, available in 3 foot and 10 foot lengths, allow for that important connection; the parents' ability to hold their child while continuously monitoring temperature.

Color-Coded Connectors

- Reduces risk of wrong connections
- Allows for simple and quick setup

Mylar® Cover

- Reflects heat generated by warmer and heat lamps

Thermistor Positioning

- Thermistor is embedded into foam so it lies flat against baby's skin allowing for better contact and more reliable monitoring



Neo-Therm® Neonatal Skin Temperature Probes

Reference Code	Description	Color	Units Per Case
STS-201	For use with: Airshields C100, C200, C550; PM78s	Pink	20
STS-202	For use with: Airshields C86; Warmers ICS; Mobile Units (most older Airshields units)	Yellow	20
STS-203	For use with: Ohmeda NIC 190 ASC Incubators; Ohio Isolette	Beige	20
STS-204	For use with: Ohmeda Incubators IC, GC, and CP; Warmers 2001, 3000, 3300, 3500, and 5000	Light Blue	20
STS-205	For use with: Hill Rom Warmers 200 and 2000 Series; Dräger	White	20
STS-206	For use with: Hill Rom Warmers 1500A, 1250A, and 2000C	Lavendar	20
STS-207	For use with: Dräger/Airshields C2000	Orange	20
STS-209	For use with: Ohmeda Medical/GE Giraffe®	Turquoise	20

Neo-Therm® Connecting Cables

Reference Code	Description	Color	Units Per Case
C201-3	3 ft. (0.91 m) Cable for use with STS-201*	Pink	1
C202-3	3 ft. (0.91 m) Cable for use with STS-202*	Yellow	1
C203-3	3 ft. (0.91 m) Cable for use with STS-203*	Beige	1
C204-3	3 ft. (0.91 m) Cable for use with STS-204*	Light Blue	1
C205-3	3 ft. (0.91 m) Cable for use with STS-205 on Hill Rom devices*	White	1
C205-3HD	3 ft. (0.91 m) Cable for use with STS-205 on HealthDyne devices*	White	1
C205-3D	3 ft. (0.91 m) Cable for use with STS-205 on Dräger devices	White	1
C206-3	3 ft. (0.91 m) Cable for use with STS-206	Lavendar	1
C207-3	3 ft. (0.91 m) Cable for use with STS-207	Orange	1
C209-3	3 ft. (0.91 m) Cable for use with STS-209*	Turquoise	1

* Available in 10 ft. (3.05 m) lengths

Neo-Therm® Replacement Probe Covers

Reference Code	Description	Units Per Case
GP-200	Hydro-Gel Probe Covers	200

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370 USA

T: 800 848 1757 F: 800 621 2590 www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Smiths Medical and Level 1 design marks, Neo-Therm and Temperature Management Solutions are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Patents pending. All products described are covered by one or more of the following U.S. Patent Nos.: 5,066,140; 5,073,034; 5,088,835; 5,161,892; 5,161,893. © 2007 Smiths Medical family of companies. All rights reserved. Literature No. TM409 2007.9