

Are You Prepared To Save Lives?



In the chaos of disaster, people look to you for help. Before disaster strikes, look to Smiths Medical for your emergency preparedness supplies. Smiths Medical products are designed for ease of use, with minimal training, and are suitable for contaminated environments. All products are available for individual purchase, as part of pre-designed kits or your own custom-designed kit. Each kit contains the critical care items needed by first responders for scenarios requiring emergency and transport ventilation and resuscitation for adults, children and infants. Smiths Medical – we're ready.

For more information, visit smiths-medical.com/bci or smiths-medical.com/pneupac.

DISASTER PREPAREDNESS



Disaster Preparedness Kits

Disaster Preparedness Kit #1

Developed for scenarios requiring a single ventilator, this kit has an extended shelf life, and requires no maintenance during that time period.

Kit contains:

- 1 – BCI® Disaster Preparedness Carrying Bag
- 1 – Pneupac® VR1 Emergency and Transport Ventilator
- 1 – Disposable Circuit
- 1 – BCI® Capnocheck® Capnometer
- 25 – Disposable Capnocheck® Adapters

Catalog #WW170000



Disaster Preparedness Kit #2

Developed for scenarios requiring multiple ventilators that will be powered by individual tanks of gas, this kit is shipped in 1 carton for easy storage until needed. All of the products require no annual maintenance.

Kit contains:

- 1 – BCI® Disaster Preparedness Carrying Bag
- 1 – PneuPac® paraPAC Emergency Transport Ventilator
- 1 – Disposable Circuit
- 1 – BCI® Capnocheck® Capnometer
- 25 – Disposable Capnocheck® Adapters

Catalog #WW170001



Disaster Preparedness Kit #3

Developed for scenarios requiring multiple ventilators that will be powered by one large tank of gas, this kit is shipped in 1 carton for easy storage until needed. All of the products included have an extended shelf life, and require no maintenance during that time period.

Kit contains:

- 1 – BCI® Disaster Preparedness Carrying Bag
- 1 – PneuPac® VR1 Emergency and Transport Ventilator
- 1 – Disposable Circuit
- 1 – BCI® Capnocheck® Capnometer
- 25 – Disposable Capnocheck® Adapters
- 1 – Multi-Ventilator Manifold

Catalog #WW170002



Design Your Own Disaster Preparedness Kit

Design your own custom Disaster Preparedness Kit to meet your specific requirements. Include as many or as few items as you need for each kit choosing from the items below.

VENTILATORS

	Catalog #
Pneupac® VR1 with Air Mix	520A1126
Pneupac® paraPAC Emergency Transport Ventilator with Alarms	120003
Pneupac® ventiPAC Emergency Transport Ventilator with Alarms	130003
Pneupac® babyPAC™ Emergency Transport Ventilator with Alarms	150003

CAPNOMETERS

Capnocheck® II Capnograph/Oximeter	8400
Capnocheck® Capnometer with Alarms	WW8110-001SYS

OXIMETERS

BCI® 3301 Hand-Held Pulse Oximeter	3301
Digit® Finger Pulse Oximeter	3420Y



DISPOSABLE CIRCUITS

Disposable Patient Circuit Kit, 10/pkg (paraPAC, ventiPAC, VR1)	122003
Disposable Patient Circuit Kit w/diaphragm, 15/case (babyPAC™)	122004

ACCESSORIES

CBRN Circuit	WW122007
Multi-Ventilator Manifold	WW331003
Disaster Preparedness Carrying Bag	WW341002VR

Capnometers/Capnographs



BCI® Capnocheck® Capnometer

The BCI® Capnocheck® Capnometer is a small portable quantitative capnometer that monitors carbon dioxide concentrations and respiratory rate in one easy to use device.
Catalog #WW8110-001SYS (with alarms)



Capnocheck® II Capnograph/Oximeter

The Capnocheck® II Handheld Capnograph/Oximeter is a full featured, highly portable unit. Its back-lit display assures easy viewing of waveforms, numeric values and on-screen trending data.
Catalog #8400 (with battery and charger)

Oximeters



BCI® 3301 Hand-Held Pulse Oximeter

The easy-to-read display and ergonomic design quickly and reliably delivers pulse oximetry, and pulse rate for patients from neonate to adult.
Catalog #3301



Digit® Finger Pulse Oximeter

A rugged pocket-size oximeter with a large LED display for use in a wide range of environments that provides SpO₂, pulse rate and pulse strength measurements for patients.
Catalog #3420Y

Transport Ventilators



Pneupac® VR1 Emergency and Transport Ventilator

The simplicity, portability and safety of the Pneupac® VR1 Emergency and Transport Ventilator provides the full-featured solution for emergency and transport ventilation resuscitation.

**Catalog #520A1126
(with Air Mix)**



Pneupac® paraPAC Emergency Transport Ventilator

Designed specifically for use by trained emergency personnel, the paraPAC enables greater control of breathing parameters. The dual controls allow easy selection of tidal volume and frequency to match your patient's ventilatory requirements. Suitable for ventilation during CPR and emergency transportation of adults, children and infants.

**Catalog #120003
(with alarms)**



Pneupac® ventiPAC Emergency Transport Ventilator

The ventiPAC is portable, compact, lightweight and easy to use in difficult situations. It operates reliably and provides alarms and monitoring similar to those found throughout the hospital. The clinician is able to alter the inspiratory and expiratory phases of ventilation to allow for patient comfort, and can also provide essential ventilatory support options such as positive end expiratory pressure (PEEP).

**Catalog #130003
(with alarms)**



Pneupac® babyPAC™ Emergency Transport Ventilator

Designed specifically for delicate neonate and infant lungs, the babyPAC™ ventilator has a unique variable gas mixing system which extends cylinder life and allows precise selection of oxygen concentration. The unit incorporates a comprehensive range of features to help control and monitor ventilation parameters at all times.

**Catalog #150003
(with alarms)**

Accessories



CBRN Circuit

Provides patient protection by filtering out all environmental contaminants such as nuclear and biological agents.

Catalog #WW122007



Multi-Ventilator Manifold

Offers Emergency and First Responder personnel the option of simultaneously running up to 8 emergency and transport ventilators off of a single tank of gas.

Catalog #WW331003



BCI® Disaster Preparedness Carrying Bag

A lightweight vinyl bag designed specifically for use by Emergency and First Responder personnel. This easily visible carrying case stores the Pneupac® VR1 Emergency and Transport Ventilator, Circuits, BCI® Capnocheck® Capnometer and Airway Adapters and has room to spare for additional devices and supplies.

Catalog #WW341002VR

Information subject to change without notice. Please read the instructions for use supplied with the product for detailed instructions, warnings and cautions.

For more information, please call Smiths Medical PM, Inc., at 262-542-3100 or 800-558-2345, or your Smiths Medical distributor.

Smiths Medical PM, Inc.

Disaster Preparedness

N7 W22025 Johnson Drive, Waukesha, Wisconsin 53186 USA

Phone: 262-542-3100 Fax: 262-542-0718 Toll-Free: 800-558-2345

www.smiths-medical.com

Smiths Medical – a part of Smiths Group plc

babyPAC, Capnocheck, Digit, Pneupac, BCI and the Smiths Medical design mark are trademarks of the Smiths Medical family of companies.

The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

©2008 Smiths Medical family of companies. All rights reserved. BC2032 Rev 00 09/08



smiths medical