

smiths medical
portex®

BLUselect® Tracheostomy Tubes

Every Breath Counts



BLUSELECT® TRACHEOSTOMY TUBES

EVERY BREATH COUNTS

BLUselect® Tracheostomy Tubes

You've asked, we've listened. After interviewing hundreds of clinicians, we have made improvements to our tracheostomy tubes to meet your patient's needs.

When caring for a patient with a tracheostomy, you should choose a tracheostomy tube that best suits their specific clinical need. The BLUselect® tracheostomy tubes have options for a variety of patients:



Color coding tubes and packages

- Tubes and packaging are color-coded by size
- Color coding intended to ensure the right size tube is used for the patient

Laser printing on flanges

- Product information cannot be altered by secretions or cleaning
- Displays the tube size prescribed to the patient

Solid or checkered markings around tube

- Indicates if the tube is fenestrated or unfenestrated
- Markings indicate whether a fenestrated or unfenestrated inner cannula should be used as prescribed for the patient

Comprehensive tube sizes

- With or without inner cannula
- Cuffed or uncuffed
- Fenestrated or unfenestrated

BLUselect® Suctionaid® tubes

- Suction created in the tube removes secretions from above the cuff to help minimize aspiration
- This suction feature can be used for ventilated patient vocalization

SOFT-SEAL® tubes

- Low pressure, high volume cuff designed to minimize trauma
- Upon deflation, the cuff rests smoothly to the tube to aid in ease of insertion

BLUselect® tracheostomy tubes are not manufactured with DEHP plasticizers.

VISIT OUR WEBSITE FOR MORE INFORMATION
WWW.SMITHS-MEDICAL.COM

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN ALL GEOGRAPHIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. BLUselect, Portex and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries ©2019 Smiths Medical. All rights reserved. TR194560EN-102019

**Rx
ONLY**

MMSPCA-0566

smiths medical