

## General Anaesthesia Airway

Breathing Systems Catalogue



Breathing  
Circuits



Face  
Masks



Circuit  
Accessories





## Contents

|  |    |
|--|----|
| Expandable Length Dual-Limb Circuits .....               | 4  |
| Fixed Length Dual-Limb Circuits.....                     | 7  |
| Fixed Length Single-Limb Circuits .....                  | 10 |
| Custom Circuits .....                                    | 11 |
| VALUE PACK Custom Kits.....                              | 12 |
| Premium Soft Plus™ Face Masks with Inflation Valve ..... | 13 |
| Breathing Bags .....                                     | 14 |
| Bacterial / Viral Filter .....                           | 15 |
| Heat and Moisture Exchangers (HMEs).....                 | 16 |
| Filtered Heat and Moisture Exchangers (HMEFs).....       | 17 |
| Accessories .....  | 18 |
| Index.....   | 19 |

## Expandable Length Dual-Limb Circuits

### ORDERING INFORMATION

| <b>2x Adult 42"-120" (107 cm-305 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |                    |                    |               |                |                |
|--|-----------|--------------------|--------------------|---------------|----------------|----------------|
| Item Number  | Temp Port | Filter(s)          | Gas Sampling Line  | Breathing Bag | Breathing Mask | Units per Case |
| 450950-NL  |           | Breathing          | GSL without filter | 3.0L (NL)     | 5245           | 12             |
| 450951-NL  | Yes       | Breathing (2), HME | GSL without filter | 3.0L (NL)     | 5245           | 12             |
| 450952-NL  |           | Breathing (2)      |                    | 3.0L (NL)     |                | 12             |
| 451040-NL  |           | HEPA               |                    | 3.0L (NL)     |                | 12             |

| <b>2x Adult 42"-120" (107 cm-305 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |                 |               |               |                |
|--|-----------------|---------------|---------------|----------------|
| Item Number  | Tubing          | Filter(s)     | Breathing Bag | Units per Case |
| 370900-NL  | 2x144" (366 cm) | Breathing (2) | 3.0L (NL)     | 10             |

| <b>2x Adult 42"-108" (107 cm-274 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |                |                   |               |                |                |
|--|-----------|----------------|-------------------|---------------|----------------|----------------|
| Item Number  | Temp Port | Filter(s)      | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 450953-NL  |           | Breathing      | GSL no filter     | 3.0L (NL)     | 5245           | 20             |
| 450954-NL  |           | Breathing      |                   | 3.0L (NL)     |                | 20             |
| 451041-NL  |           | HEPA           | GSL no filter     | 3.0L (NL)     |                | 12             |
| 450955-NL  |           | Breathing (2)  |                   | 3.0L (NL)     |                | 20             |
| 450956-NL  |           | Breathing, HME | GSL no filter     | 3.0L (NL)     | 5245           | 12             |

| <b>2x Pediatric 60"-150" (153-381 cm) MRI Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |           |                   |               |                |                |
|--|-----------|-----------|-------------------|---------------|----------------|----------------|
| Item Number  | Temp Port | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 371001-NL  |           |           | GSL no filter     | 2.0L (NL)     |                | 10             |
| 371002   |           |           | GSL no filter     |               |                | 10             |

| <b>2x Adult 120" (304 cm) Adult MRI Length Dual-Limb Breathing Circuit</b> |                                |                |
|--|--------------------------------|----------------|
| Item Number  | Tubing                         | Units per Case |
| 487000   | 2x 20"-40"-60" (51-102-152 cm) | 20             |



## ORDERING INFORMATION

| <b>2x Adult 33"-87" (84-221 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |                |                    |                        |                |                |
|---|-----------|----------------|--------------------|------------------------|----------------|----------------|
| Item Number   | Temp Port | Filter(s)      | Gas Sampling Line  | Breathing Bag          | Breathing Mask | Units per Case |
| 450960-NL   |           | Breathing      | GSL without filter | 3.0L (NL)              | 5245           | 20             |
| 450961-NL   |           | Breathing (2)  |                    | 3.0L (NL)              |                | 20             |
| 450962-NL   |           | Breathing      | GSL without filter | 3.0L (NL)              | 005445         | 20             |
| 450963-NL   |           | Breathing      | GSL without filter | 3.0L (NL)              |                | 20             |
| 451050-NL   |           | HEPA           |                    | 3.0L (NL)              |                | 20             |
| 450964-NL   |           | Breathing (2)  | GSL without filter | 3.0L (NL)              | 5245           | 20             |
| 451051-NL   |           | HME            | GSL without filter | 3.0L (NL)              | 5245           | 20             |
| 450965-NL   | Yes       | HME, Breathing |                    | 3.0L (NL)              | 5245           | 20             |
| 450904  |           | BV             |                    | 3.0L (L)               |                | 20             |
| 650904  |           | BV             |                    | 3.0L (NL)              |                | 20             |
| 450904-SG   |           | BV             |                    | 3.0L (L) Suregrip™ Bag |                | 20             |
| 450907  |           | BV             |                    | 2.0L (L)               |                | 20             |
| 450907-NL   |           | BV             |                    | 2.0L (NL)              |                | 20             |
| 450906  | Yes       | BV             |                    | 3.0L (L)               |                | 20             |
| 450906-NL   | Yes       | BV             |                    | 3.0L (NL)              |                | 20             |
| 451004  |           |                |                    | 3.0L (L)               |                | 20             |
| 451004-NL   |           |                |                    | 3.0L (NL)              |                | 20             |
| <b>Elbow: GSP WYE: Bifurcated with port</b>   |           |                |                    |                        |                |                |
| 454204  |           | BV             |                    | 3.0L (L)               |                | 20             |
| <b>Elbow: GSP WYE: Parallel with Port</b>   |           |                |                    |                        |                |                |
| 455477D   | Yes       | BV             | GSL without filter | 3.0L (L)               |                | 20             |
| 455417D   | Yes       | BV             | GSL with filter    | 3.0L (L)               |                | 20             |
| <b>Elbow: Std. WYE: Bifurcated with port</b>  |           |                |                    |                        |                |                |
| 455300  |           |                |                    | 3.0L (L)               |                | 20             |
| <b>Elbow: - WYE: Swivel</b>   |           |                |                    |                        |                |                |
| 454000  |           |                |                    | 3.0L (L)               |                | 20             |
| 453600  |           | BV             |                    | 3.0L (L)               |                | 20             |
| 4509901D  |           | Breathing (2)  | GSL without filter | 3.0L (L)               |                | 20             |
| <b>Elbow: GSP WYE: Straight with port</b>   |           |                |                    |                        |                |                |
| 453600-NL   |           | BV             |                    | 3.0L (NL)              |                | 20             |



## ORDERING INFORMATION

| <b>2x Pediatric 40"-100" (102-254 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |                    |                    |               |                |                |
|---|--------------------|--------------------|---------------|----------------|----------------|
| Item Number   | Filter(s)          | Gas Sampling Line  | Breathing Bag | Breathing Mask | Units per Case |
| 490900-NL   | 600 HME, Breathing | GSL without filter | 1.0L (NL)     |                | 15             |
| 490901-NL   | 600HME, Breathing  |                    | 2.0L (NL)     |                | 15             |
| 490902-NL   | Breathing (2)      |                    | 1.0L (NL)     |                | 20             |
| 490903-NL   | Breathing (2)      |                    | 2.0L (NL)     |                | 20             |
| 491000-NL   | HEPA               | GSL without filter | 2.0L (NL)     | 5248           | 15             |

| <b>2x Pediatric 40"-100" (102-254 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |  |               |                |
|---|-----------|--|---------------|----------------|
| Item Number   | Filter(s) |  | Breathing Bag | Units per Case |
| 490804  |           |  | 1.0L (L)      | 20             |
| 490804-NL   |           |  | 1.0L (NL)     | 20             |
| 492804  | BV        |  | 1.0L (L)      | 20             |
| 492804-NL   | BV        |  | 1.0L (NL)     | 20             |

| <b>2x Pediatric 20"-50" (51-127 cm) Expandable Length Dual-Limb Breathing Circuit</b><br>Elbow: GSP Wye: Straight |           |  |               |                |
|---|-----------|--|---------------|----------------|
| Item Number   | Filter(s) |  | Breathing Bag | Units per Case |
| 480804  |           |  | 1.0L (L)      | 40             |
| 482804  | BV        |  | 1.0L (L)      | 40             |
| 482804-NL   | BV        |  | 1.0L (NL)     | 40             |
| 482802  | BV        |  | 2.0L (L)      | 40             |
| 482802-NL   | BV        |  | 2.0L (NL)     | 40             |



## Fixed Length Dual-Limb Circuits

### ORDERING INFORMATION

| <b>2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit</b> |           |                        |                |                |  |
|---|-----------|------------------------|----------------|----------------|--|
| Elbow: GSP Wye: Straight  |           |                        |                |                |  |
| Item Number   | Filter(s) | Breathing Bag          | Breathing Mask | Units per Case |  |
| 225-3709-804  | BV        | 3.0L (L)               |                | 20             |  |
| 370936  | (2) BV    | 3.0L (L)               |                | 20             |  |
| 370936-NL   | (2) BV    | 3.0L (NL)              |                | 20             |  |
| 670904  | BV        | 3.0L (NL)              |                | 20             |  |
| 370924  | BV        | 3.0L (L)               | 5045           | 15             |  |
| Elbow: GSP WYE: Swivel  |           |                        |                |                |  |
| 370374-SG   | BV        | 3.0L (L) SureGrip™ Bag |                | 20             |  |
| 370374-NL   | BV        | 3.0L (NL)              |                | 20             |  |
| Elbow: Straight, GSP WYE: Swivel                                      |           |                        |                |                |  |
| 225-3703-810  | BV        | 3.0L (L)               |                | 20             |  |
| Elbow: Standard WYE: Straight   |           |                        |                |                |  |
| 225-3709-800  | BV        | 3.0L (L)               |                | 20             |  |
| Elbow: - WYE: Swivel with port  |           |                        |                |                |  |
| 373600  | BV        | 3.0L (L)               |                | 20             |  |
| Elbow: - WYE: Swivel  |           |                        |                |                |  |
| 225-3703-800  | BV        | 3.0L (L)               |                | 20             |  |

### ORDERING INFORMATION

| <b>2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit</b> |           |                    |               |                |                |
|---|-----------|--------------------|---------------|----------------|----------------|
| Elbow: GSP Wye: Straight  |           |                    |               |                |                |
| Item Number   | Filter(s) | Gas Sampling Line  | Breathing Bag | Breathing Mask | Units per Case |
| 370941  | BV        | GSL without filter | 3.0L (L)      | 5045           | 15             |
| 370941-NL   | BV        | GSL without filter | 3.0L (NL)     | 5045           | 15             |
| 370956D   | BV        | GSL without filter | 3.0L (L)      |                | 20             |
| 370956D-NL  | BV        | GSL without filter | 3.0L (NL)     |                | 20             |
| 370928  | BV        | GSL with filter    | 3.0L (L)      |                | 20             |
| 370928-NL   | BV        | GSL with filter    | 3.0L (NL)     |                | 20             |
| 371056  |           | GSL without filter | 3.0L (L)      |                | 20             |
| 371056-NL   |           | GSL without filter | 3.0L (NL)     |                | 20             |
| 225-3710-804  |           |                    | 3.0L (L)      |                | 20             |
| 671004  |           |                    | 3.0L (NL)     |                | 20             |
| 225-3710-807  |           |                    | 2.0L (L)      |                | 20             |

## ORDERING INFORMATION

### 2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Straight

|              |          |    |
|--------------|----------|----|
| 225-3710-802 | 2.0L (L) | 20 |
| 225-3710-800 | 3.0L (L) | 20 |

### 2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Swivel

|              |          |    |
|--------------|----------|----|
| 374000       | 3.0L (L) | 20 |
| 225-3704-800 | 3.0L (L) | 20 |

### 2x Adult 60" (152 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight

| Item Number     | Temp Port | Filter(s)     | Gas Sampling Line  | Breathing Bag | Breathing Mask | Units per Case |
|-----------------|-----------|---------------|--------------------|---------------|----------------|----------------|
| 370913-NL       | YES       | Breathing     | GSL without filter | 3.0L (NL)     | 5045           | 20             |
| 370916-NL       |           | Breathing     | GSL without filter | 2.0L (NL)     |                | 20             |
| 370923-NL       |           | Breathing (2) |                    | 3.0L (NL)     |                | 20             |
| 371804          |           | BV            |                    | 3.0L (L)      |                | 15             |
| 371804-NL       |           | BV            |                    | 3.0L (NL)     |                | 15             |
| 225-3718-806    | YES       | BV            |                    | 3.0L (L)      |                | 15             |
| 225-3718-806-NL | YES       | BV            |                    | 3.0L (NL)     |                | 15             |
| 225-3719-804    |           |               |                    | 3.0L (L)      |                | 15             |
| 225-3719-804-NL |           |               |                    | 3.0L (NL)     |                | 15             |

### Elbow: Standard Wye: Straight

|              |  |           |  |          |  |    |
|--------------|--|-----------|--|----------|--|----|
| 225-3718-800 |  | Breathing |  | 3.0L (L) |  | 15 |
| 225-3719-800 |  |           |  | 3.0L (L) |  | 15 |

### Elbow: Standard, Straight GS Connector Wye: Straight

|        |     |           |                    |          |  |    |
|--------|-----|-----------|--------------------|----------|--|----|
| 371821 | Yes | Breathing | GSL without filter | 3.0L (L) |  | 15 |
|--------|-----|-----------|--------------------|----------|--|----|





## ORDERING INFORMATION

| <b>2x Adult 72" (183 cm) Fixed Length Dual-Limb Breathing Circuit</b> |              |                    |               |                |                |
|---|--------------|--------------------|---------------|----------------|----------------|
| Elbow: GSP Wye: Straight  |              |                    |               |                |                |
| Item Number   | Filter(s)    | Gas Sampling Line  | Breathing Bag | Breathing Mask | Units per Case |
| 370902-NL   | BV (2)       |                    | 3.0L (NL)     |                | 20             |
| 370905-NL   | BV           | GSL without filter | 3.0L (NL)     | 5245           | 20             |
| 370908-NL   | BV Temp Port | GSL without filter | 3.0L (NL)     | 5045           | 20             |
| 370909-NL   | BV HME       | GSL without filter | 3.0L (L)      | 5045           | 20             |
| 225-3724-800  |              |                    | 3.0L (L)      |                | 15             |
| <b>Elbow: GSP Wye: Swivel</b>   |              |                    |               |                |                |
| 370300-NL   | BV           |                    | 3.0L (NL)     |                | 20             |

| <b>2x Adult 144" (366 cm) Fixed Length Dual-Limb Breathing Circuit</b> |           |                   |               |                |                |
|--|-----------|-------------------|---------------|----------------|----------------|
| Elbow: Standard Wye: Straight  |           |                   |               |                |                |
| Item Number  | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 225-3723-800   | BV        |                   | 3.0L (L)      |                | 15             |

| <b>2x Pediatric 36" (91 cm) Fixed Length Dual-Limb Breathing Circuit</b> |           |                   |               |                |                |
|--|-----------|-------------------|---------------|----------------|----------------|
| Elbow: Standard Wye: Straight  |           |                   |               |                |                |
| Item Number  | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 225-3720-800   | BV        |                   | 1.0L (L)      |                | 20             |
| 225-3708-802   |           |                   | 2.0L (L)      |                | 20             |
| 225-3720-802   | BV        |                   | 2.0L (L)      |                | 20             |
| <b>Elbow: GSP Wye: Straight</b>  |           |                   |               |                |                |
| 225-3708-804   |           |                   | 1.0L (L)      |                | 20             |
| 670804   |           |                   | 1.0L (NL)     |                | 20             |
| 672004   | BV        |                   | 1.0L (NL)     |                | 20             |
| 225-3720-804   | BV        |                   | 1.0L (L)      |                | 20             |
| 372036   | (2) BV    |                   | 1.0L (L)      |                | 20             |
| 370807   |           |                   | 2.0L (L)      |                | 20             |
| 372007   | BV        |                   | 2.0L (L)      |                | 20             |

| <b>2x Pediatric 54" (137 cm) Fixed Length Dual-Limb Breathing Circuit</b> |           |                   |               |                |                |
|---|-----------|-------------------|---------------|----------------|----------------|
| Elbow: Standard Wye: Straight   |           |                   |               |                |                |
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 372900  |           |                   | 1.0L (L)      |                | 15             |
| <b>Elbow: GSP Wye: Straight</b>   |           |                   |               |                |                |
| 372904  |           |                   | 1.0L (L)      |                | 15             |
| 375204  | BV        |                   | 1.0L (L)      |                | 15             |



## Fixed Length Single-Limb Circuits

### ORDERING INFORMATION

| <b>2x Adult 40" (101 cm) Fixed Length Single-Limb Breathing Circuit</b><br>Elbow: GSP |           |                   |               |                |                |
|---|-----------|-------------------|---------------|----------------|----------------|
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 384000  | Insp.     |                   | 3.0L (L)      |                | 20             |
| 384000-NL   | Insp.     |                   | 3.0L (NL)     |                | 20             |
| 384001-NL   | Insp.     |                   | 2.0L (NL)     |                | 20             |
| 384004  | Insp.     | Yes               | 3.0L (L)      | 5445           | 20             |
| 384004-NL   | Insp.     | Yes               | 3.0L (NL)     | 5445           | 20             |
| 384006  | Insp.     | Yes               | 3.0L (L)      | 5465           | 20             |
| 384006-NL   | Insp.     | Yes               | 3.0L (NL)     | 5465           | 20             |

| <b>2x Adult 60" (152 cm) Fixed Length Single-Limb Breathing Circuit</b><br>Elbow: GSP |           |                   |                        |                |                |
|---|-----------|-------------------|------------------------|----------------|----------------|
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag          | Breathing Mask | Units per Case |
| 386000  | Insp.     |                   | 3.0L (L)               |                | 15             |
| 386000-NL   | Insp.     |                   | 3.0L (NL)              |                | 15             |
| 386001  | Insp.     |                   | 3.0L (L) SureGrip™ Bag |                | 15             |
| 386003  | Insp.     |                   | 3.0L (L) SureGrip™ Bag | 5445           | 15             |

| <b>2x Adult 72" (183 cm) Fixed Length Single-Limb Breathing Circuit</b><br>Elbow: GSP |           |                   |               |                |                |
|---|-----------|-------------------|---------------|----------------|----------------|
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 387901  |           |                   |               |                | 15             |
| 387000-NL   | Insp.     |                   | 3.0L (NL)     |                | 15             |

| <b>2x Pediatric 40" (101 cm) Fixed Length Single-Limb Breathing Circuit</b><br>Elbow: GSP |           |                   |               |                |                |
|---|-----------|-------------------|---------------|----------------|----------------|
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 364001-NL   | Insp.     |                   | 1.0L (NL)     |                | 20             |
| 364002-NL   | Insp.     |                   | 2.0L (NL)     |                | 20             |

| <b>2x Pediatric 60" (152 cm) Fixed Length Single-Limb Breathing Circuit</b><br>Elbow: GSP |           |                   |               |                |                |
|---|-----------|-------------------|---------------|----------------|----------------|
| Item Number   | Filter(s) | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 366001-NL   | Insp.     |                   | 1.0L (NL)     |                | 15             |
| 366002-NL   | Insp.     |                   | 2.0L (NL)     |                | 15             |

| <b>Manual Ventilator Transport Kits</b> |   |                   |               |                |                |
|---|---|-------------------|---------------|----------------|----------------|
| Item Number                             | Tubing  | Gas Sampling Line | Breathing Bag | Breathing Mask | Units per Case |
| 385000                                  | Manual Ventilator Transport Kit with Dial Bleed Valve, Manometer Port with Tethered Cap without Breathing Bag   |                   |               |                | 25             |
| 385100                                  | Manual Ventilator Transport Kit with Dial Bleed Valve, Manometer Port with Tethered Cap with 3.0L Breathing Bag |                   |               |                | 20             |

# Custom Circuits

## CREATING YOUR OWN

Create your own custom order to package breathing circuit accessories and components for your specific procedures and protocols.

Choose from a wide variety of components to help streamline the ordering process (circuit tubing, wye connectors, patient elbows, breathing filters, HMEs, gas sampling lines, breathing masks, breathing bags, etc.)

Ordering custom kits, in some cases, can be more cost-effective than purchasing the components separately.

- Time savings. All the items from set-up to procedural application are in one package.
- Simple ordering. A single order for each circuit contains all the components you need, improving inventory management.
- Reduced wastage. The exact number of components is provided according to patient need and department protocol.
- Reduced stock. Custom circuits can be held by your local distributor according to forecast demand.



| Components         | Options  |
|--------------------|--|
| Breathing Circuits | Z-wave tubing, Square wave tubing, Expandable tubing, Adult and Pediatric single limb tubing |
| Breathing Masks    | Premium™, Premium Plus™, Premium Soft Plus™ disposable anesthesia face masks                 |
| Elbow Connectors   | 90° Elbow, 90° Elbow with gas sampling port and Straight connectors                          |
| HMEs and Filters   | Thermovent® HEPA Filters, HMEs   |
| Wye Connector      | Standard, Surefit, Straight, Parallel and Swivel connectors                                  |
| Gas Sampling       | Internal, External, Filters, Adapters  |
| Caps               | Dust, Female and Male Luer lock caps   |
| Breathing Bags     | Latex, Non-Latex, SureGrip™ Breathing Bags   |
| Ventilator Tubes   | With / Without Filters   |

## VALUE PACK Custom Kits

### CREATING YOUR OWN

- Clinicians can operate faster and more efficiently with a comprehensive kit of components.
- Operating Rooms can be set up faster while lowering turnover times.
- Value Pack means fewer items to stock, helping to make inventory more manageable, reduce re-ordering, and improve overall profitability for your healthcare facility.

Here is a list of optional components for a value Pack:

| Components              | Options  |
|-------------------------|--|
| Breathing Circuits      | Standard or Custom   |
| Breathing Masks         | Premium™, Premium Plus™, Premium Soft Plus™ disposable anesthesia face masks                   |
| Breathing Bags          | Latex, Non-Latex, SureGrip™ Breathing Bags   |
| HMEs and Filters        | Thermovent® HEPA Filters, HMEs   |
| Ventilator Tubes        | With / Without Filters   |
| ECG Pads                | 3 Pads / Pack, 5 Pads / Pack, Back Pad   |
| Suction Devices         | Yankauer Bulb Tip, Suction Catheter with Valve, Suction Catheter without Valve, Suction Tubing |
| Endotracheal Devices    | Sizes 7.0-8.0 mm<br>10 mL Syringes<br>Intubation Stylets                                       |
| Laryngoscope Blades     | Disposables in Macintosh or Miller styles  |
| Oral Airways            | Nasopharyngeal Airways I.D. 6.0-9.0 mm<br>Berman Sizes 6.0-10.0 cm<br>Guedel Sizes 6.0-10.0 cm |
| Esophageal Stethoscopes | 400 Series, 12-24 french<br>700 Series, 12-24 french<br>700 Series Sub, 12-24 french           |
| Temperature Circuits    | Rectal Probe, Skin Probe, ACCU-BAR™ Temperature Strip  |
| Blood Pressure Cuffs    | Single or Double Tube  |
| Syringes                | 10 mL or 20 mL   |
| Supplies                | Transpore Tape, Alcohol Prep Pads, Tongue Blades, Gloves, 3" x 3" Gauze Pads                   |



## Face Masks

### ORDERING INFORMATION

| Portex® Premium Soft Plus™ Face Masks with Inflation Valve |                |            |                |
|--|----------------|------------|----------------|
| Item Number  | Face Mask Size | Hook Ring  | Units per Case |
| 005471   | 1              |            | 30             |
| 005452   | 2              |            | 30             |
| 005448   | 3              | Clear      | 30             |
| 005455   | 4              | Yellow     | 30             |
| 005445   | 5              | Light Blue | 30             |
| 005465   | 6              | Dark Blue  | 30             |



| Portex® Disposable Face Masks without Inflation Valve |                |           |                |
|---|----------------|-----------|----------------|
| Item Number   | Face Mask Size | Hook Ring | Units per Case |
| 15071   | 1              | White     | 50             |
| 15052   | 2              | Pink      | 50             |
| 15070   | 3              | Yellow    | 50             |
| 15048   | 4              | Green     | 50             |
| 15055   | 5              | Red       | 50             |
| 15045   | 6              | Dark Blue | 50             |
| 15266   | 7              | Fuchsia   | 50             |



| Portex® Disposable Face Masks with Inflation Valve |                |           |                |
|--|----------------|-----------|----------------|
| Item Number  | Face Mask Size | Hook Ring | Units per Case |
| 15271  | 1              | White     | 50             |
| 15252  | 2              | Pink      | 50             |
| 15270  | 3              | Yellow    | 50             |
| 15248  | 4              | Green     | 50             |
| 15255  | 5              | Red       | 50             |
| 15245  | 6              | Dark Blue | 50             |
| 15265  | 7              | Fuchsia   | 50             |



## Breathing Bags

### ORDERING INFORMATION

| Breathing Bags, Latex |  |      |         |                |  |
|-----------------------|--|------|---------|----------------|--|
| Item Number           | Description                                    | Size | Bushing | Units per Case |  |
| 321525                | Anesthesia Breathing Bag, Natural Rubber Latex | 0.5L | 22.0 mm | 25             |  |
| 321625                | Anesthesia Breathing Bag, Natural Rubber Latex | 1.0L | 22.0 mm | 25             |  |
| 321725                | Anesthesia Breathing Bag, Natural Rubber Latex | 2.0L | 22.0 mm | 25             |  |
| 321825                | Anesthesia Breathing Bag, Natural Rubber Latex | 3.0L | 22.0 mm | 25             |  |



| Breathing Bags, Non-Latex |   |      |         |                |  |
|---------------------------|---|------|---------|----------------|--|
| Item Number               | Description                                   | Size | Bushing | Units per Case |  |
| 670012                    | Anesthesia Breathing Bag, Neoprene, Non-Latex | 0.5L | 22.0 mm | 25             |  |
| 670001                    | Anesthesia Breathing Bag, Neoprene, Non-Latex | 1.0L | 22.0 mm | 25             |  |
| 670002                    | Anesthesia Breathing Bag, Neoprene, Non-Latex | 2.0L | 22.0 mm | 25             |  |
| 670003                    | Anesthesia Breathing Bag, Neoprene, Non-Latex | 3.0L | 22.0 mm | 25             |  |



## Bacterial / Viral Filter

### ORDERING INFORMATION

| Bacterial / Viral Filter |                  |           |   |  |                |             |               |               |                       |                |
|--------------------------|------------------|-----------|---|--|----------------|-------------|---------------|---------------|-----------------------|----------------|
| Item Number              | Description      | Tidal VOL | Resistance (mg H <sub>2</sub> O at L/min) | Moisture Output (mg H <sub>2</sub> O / L@Vt) | Deadspace (mL) | Weight (gm) | Media         | Filter        | Filter Efficiency (%) | Units per Case |
| 002862                   | Breathing Filter | n/a       | 1.0@30<br>2.0@60<br>3.5@90                | NA   | 32             | 20          | Polypropylene | Electrostatic | >99.9                 | 40             |
| 100/585/000N             | Thermovent® HEPA | 150-1200  | 1.5@30<br>3.3@60<br>5.0@90                | 18.6@500<br>15.6@1000                        | 45             | 29          | Paper         | Hydrophobic   | >99.99999             | 50             |
| 002873                   | HEPA Filter      | 200-1500  | 1.5@30                                    | 25 @500                                      | 55             | 40          | Paper         | Hydrophobic   | 99.9999               | 50             |



## Heat and Moisture Exchangers (HMEs)

## ORDERING INFORMATION

| Heat and moisture Exchangers (HMEs) |  |                |   |  |                |             |                  |                |
|-------------------------------------|--|----------------|---|--|----------------|-------------|------------------|----------------|
| Item Number                         | Description  | Tidal VOL (mL) | Resistance (mg H <sub>2</sub> O at L/min) | Moisture Output (mg H <sub>2</sub> O / L@Vt) | Deadspace (mL) | Weight (gm) | Media            | Units per Case |
| 002815                              | HME, Neonatal, Straight                                  | >10            | 0.9@5                                     | 32@25  | 3              | 3           | Hygroscopic Foam | 50             |
| 002814P                             | HME with Gas Port, Pediatric, Angled                     | >45            | 0.2@10                                    | 31@50  | 13             | 9           | Hygroscopic Foam | 50             |
| 580011                              | Thermovent® 600 HME, Pediatric, Straight                 | 100-600        | 1.0@30<br>2.9@60<br>5.7@90                | 27@500                                       | 11             | 10          | Hygroscopic Foam | 50             |
| 002811P                             | HME with Gas Port, Rectangular, Adult, Straight          | 90-1500        | 0.98@30                                   | 32@500                                       | 76             | 38          | Hygroscopic Foam | 50             |
| 100/582/000                         | Thermovent® 1200 HME, Adult, Straight                    | 600-1200       | 0.4@30<br>1.2@60<br>2.4@90                | 24@1000                                      | 32             | 18.6        | Hygroscopic Foam | 20             |
| 002813                              | HME with Gas Port, Adult, Small, Straight                | 90-1500        | 1.0@15                                    | 31@250                                       | 26             | 18          | Hygroscopic Foam | 50             |
| 002817P                             | HME with Gas Port, Adult, Large, Straight                | 150-1500       | 1.8@60                                    | 31@500                                       | 53             | 27          | Hygroscopic Foam | 50             |
| 002818P                             | HME with Gas Port, Adult, Small, Angled                  | 90-1500        | 0.15@15                                   | 31@250                                       | 31             | 19          | Hygroscopic Foam | 50             |
| 002812                              | HME with Gas Port, Adult, Med., Angled                   | 120-1500       | 0.22@15                                   | >31@500                                      | 42             | 21          | Hygroscopic Foam | 50             |
| 002816                              | HME with Gas Port, Adult, Large, Angled                  | 200-1500       | 0.3@30                                    | 31@500                                       | 66             | 30          | Hygroscopic Foam | 50             |
| 002841                              | HME with Flex Tube, Adult, Small, Straight               | 400-1000       | 1.0@30<br>3.0@60<br>5.0@90                | 30@500<br>27@1000                            | 70             | 27          | Hygroscopic Foam | 20             |
| 002841P                             | HME with Gas Port, and Flex Tube, Adult, Small, Straight | 400-1000       | 1.0@30<br>3.0@60<br>5.0@90                | 30@500<br>27@1000                            | 70             | 27          | Hygroscopic Foam | 20             |





## Filtered Heat and Moisture Exchangers (HMEFs)

### ORDERING INFORMATION

| filtered Heat and moisture Exchangers (HMEFs) |   |                |   |  |                |             |                         |                       |                |
|---|---|----------------|---|--|----------------|-------------|-------------------------|-----------------------|----------------|
| Item Number                                   | Description   | Tidal VOL (mL) | Resistance (mg H <sub>2</sub> O at L/min) | Moisture Output (mg H <sub>2</sub> O / L@Vt) | Deadspace (mL) | Weight (gm) | Media                   | Filter Efficiency (%) | Units per Case |
| 002822  | Filtered HME with Gas Port, Rectangular, Adult, Straight          | 90-1500        | 0.98@30                                   | 32@500                                       | 76             | 38          | Hygroscopic Foam        | Electrostatic         | 50             |
| 002874  | HEPA Filtered HME with Gas Port, Adult, Straight                  | 150-1500       | 1.48@30                                   | 31@500                                       | 55             | 40          | Hygroscopic Foam, Paper | Hydrophobic           | 50             |
| 002823P                                       | Filtered HME with Gas Port, Small Adult, Straight                 | 120-1500       | 1.3@15                                    | >31@500                                      | 33             | 21          | Hygroscopic Foam        | Electrostatic         | 50             |
| 002825  | Filtered HME with Gas Port, Small Adult, Angled                   | 90-1500        | 1.1@15                                    | 31@250                                       | 32             | 18          | Hygroscopic Foam        | Electrostatic         | 50             |
| 002821  | Filtered HME with Gas Port, Large Adult, Angled                   | 200-1500       | 1.0@30                                    | 31@500                                       | 66             | 31          | Hygroscopic Foam        | Electrostatic         | 50             |
| 002865  | Filtered HME, Adult, Straight                                     | 400-1000       | 1.6@30<br>3.8@60<br>6.5@90                | 27@500<br>24@1000                            | 32             | 20          | Hygroscopic Foam        | Electrostatic         | 40             |
| 002866  | Filtered HME with Flex Tube, Adult, Straight                      | 400-1000       | 1.6@30<br>3.8@60<br>6.5@90                | 27@500<br>24@1000                            | 67             | 34          | Hygroscopic Foam        | Electrostatic         | 20             |
| 002850P                                       | Filtered HME with Gas Sampling Port, Adult, Straight              | 400-1000       | 1.5@30<br>3.0@60<br>5.0@90                | 28@500<br>26@1000                            | 68             | 38          | Hygroscopic Foam        | Electrostatic         | 20             |
| 002851P                                       | Filtered HME with Gas Sampling Port and Flex Tube Adult, Straight | 400-1000       | 1.5@30<br>3.0@60<br>5.0@90                | 28@500<br>26@1000                            | 103            | 52          | Hygroscopic Foam        | Electrostatic         | 20             |

| Speciality Filters |                                      |                |   |  |                |             |                    |                     |                       |                |
|--------------------|--------------------------------------|----------------|---|--|----------------|-------------|--------------------|---------------------|-----------------------|----------------|
| Item Number        | Description                          | Tidal VOL (mL) | Resistance (mg H <sub>2</sub> O at L/min) | Moisture Output (mg H <sub>2</sub> O / L@Vt) | Deadspace (mL) | Weight (gm) | Media              | Filter              | Filter Efficiency (%) | Units per Case |
| 002291             | 0.5µ Hydrophobic Disc Filter (Clear) | n/a            | n/a                                       | n/a  | 0.75           | 3.8         | Gore-Tex® Membrane | Direct Interception | >99.9                 | 200            |



# Accessories

## ORDERING INFORMATION

| Gas Sampling Lines |                                 |                 |            |                |
|--------------------|---------------------------------|-----------------|------------|----------------|
| Item Number        | Description                     | Luer Connection | Filter     | Units per Case |
| 225-3422-800       | Gas Sampling Line, 10'x.050" ID | Male / Male     |            | 50             |
| 225-3450-800       | Gas Sampling Line, 10'x.050" ID | Male / Female   |            | 50             |
| 225-3421-800       | Gas Sampling Line, 10'x.050" ID | Male / Male     | 0.5 micron | 50             |
| 225-3412-800       | Gas Sampling Line, 10'x.050" ID | Male / Female   | 0.5 micron | 50             |
| 225-3415-800       | Gas Sampling Line, 10'x.050" ID | Male / Female   | 5.0 micron | 50             |
| 225-3416-800       | Gas Sampling Line, 15'x.050" ID | Male / Female   | 5.0 micron | 50             |
| 225-3492-300       | Gas Sampling Line, 10'x.047" ID | Male / Male     |            | 100            |
| 225-3445-300       | Gas Sampling Line, 10'x.047" ID | Male / Female   |            | 100            |



| Patient Connectors (Elbows) |  |                |
|-----------------------------|--|----------------|
| Item Number                 | Description  | Units per Case |
| 225-3507-800                | (A) Patient Connector, Standard Elbow, Clear                         | 50             |
| 225-3523-804                | (B) Patient Connector, Gas Sampling Elbow with Capped Port, Clear    | 50             |
| 225-3524-804                | (C) Patient Connector, Gas Sampling Straight with Capped Port, Clear | 50             |



| Ventilator Tubing |  |                |
|-------------------|--|----------------|
| Item Number       | Description  | Units per Case |
| 225-3714-800      | Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter | 30             |
| 225-3717-800      | Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter | 30             |
| 225-3715-800      | Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter | 20             |
| 225-3716-800      | Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter | 20             |



| CO <sub>2</sub> Absorbents |  |                |
|----------------------------|--|----------------|
| Item Number                | Description  | Units per Case |
| 008860                     | Sodasorb® CO <sub>2</sub> Absorbent, Canister Pak (Refill Bag) (1.64 L)    | 12             |
| 008870                     | Sodasorb® CO <sub>2</sub> Absorbent, Pre Pak (Canister Discs) (1.3 L)      | 12             |
| 008880                     | Sodasorb® CO <sub>2</sub> Absorbent, Pail (5 Gallon) (19 L)                | 1              |
| 008865                     | Sodasorb® LF CO <sub>2</sub> Absorbent, Canister Pak (Refill Bag) (1.64 L) | 12             |
| 008875                     | Sodasorb® LF CO <sub>2</sub> Absorbent, Pre Pak (Canister Discs) (1.3 L)   | 12             |
| 008855                     | Sodasorb® LF CO <sub>2</sub> Absorbent, Pail (1 Gallon) (4.3L or 4.2 kg)   | 3              |



# Index

## ITEM NUMBERS

|                       |    |               |    |                  |    |
|-----------------------|----|---------------|----|------------------|----|
| 100/582/000 .....     | 16 | 002813 .....  | 16 | 15270 .....      | 13 |
| 100/585/ 000N .....   | 15 | 002814P ..... | 16 | 15271 .....      | 13 |
| 225-3412-800 .....    | 18 | 002815 .....  | 16 | 321525 .....     | 14 |
| 225-3415-800 .....    | 18 | 002816 .....  | 16 | 321625 .....     | 14 |
| 225-3416-800 .....    | 18 | 002817P ..... | 16 | 321725 .....     | 14 |
| 225-3421-800 .....    | 18 | 002818P ..... | 16 | 321825 .....     | 14 |
| 225-3422-800 .....    | 18 | 002821 .....  | 17 | 364001-NL .....  | 10 |
| 225-3445-300 .....    | 18 | 002822 .....  | 17 | 364002-NL .....  | 10 |
| 225-3450-800 .....    | 18 | 002823P ..... | 17 | 366001-NL .....  | 10 |
| 225-3492-300 .....    | 18 | 002825 .....  | 17 | 366002-NL .....  | 10 |
| 225-3507-500 .....    | 11 | 002841 .....  | 16 | 370300-NL .....  | 9  |
| 225-3507-800 .....    | 18 | 002841P ..... | 16 | 370374-NL .....  | 7  |
| 225-3523-804 .....    | 18 | 002850P ..... | 17 | 370374-SG .....  | 7  |
| 225-3524-804 .....    | 18 | 002851P ..... | 17 | 370807 .....     | 9  |
| 225-3703-800 .....    | 7  | 002862 .....  | 15 | 370900-NL .....  | 4  |
| 225-3703-810 .....    | 7  | 002865 .....  | 17 | 370902-NL .....  | 9  |
| 225-3704-800 .....    | 8  | 002866 .....  | 17 | 370905-NL .....  | 9  |
| 225-3708-802 .....    | 9  | 002873 .....  | 15 | 370908-NL .....  | 9  |
| 225-3708-804 .....    | 9  | 002874 .....  | 17 | 370909-NL .....  | 9  |
| 225-3709-800 .....    | 7  | 005445 .....  | 13 | 370913-NL .....  | 8  |
| 225-3709-804 .....    | 7  | 005448 .....  | 13 | 370916-NL .....  | 8  |
| 225-3710-800 .....    | 8  | 005452 .....  | 13 | 370923-NL .....  | 8  |
| 225-3710-802 .....    | 8  | 005455 .....  | 13 | 370924 .....     | 7  |
| 225-3710-804 .....    | 7  | 005465 .....  | 13 | 370928 .....     | 7  |
| 225-3710-807 .....    | 7  | 005471 .....  | 13 | 370928-NL .....  | 7  |
| 225-3714-800 .....    | 18 | 008855 .....  | 18 | 370936 .....     | 7  |
| 225-3715-800 .....    | 18 | 008860 .....  | 18 | 370936-NL .....  | 7  |
| 225-3716-800 .....    | 18 | 008865 .....  | 18 | 370941 .....     | 7  |
| 225-3717-800 .....    | 18 | 008870 .....  | 18 | 370941-NL .....  | 7  |
| 225-3718-800 .....    | 8  | 008875 .....  | 18 | 370956D .....    | 7  |
| 225-3718-806 .....    | 8  | 008880 .....  | 18 | 370956D-NL ..... | 7  |
| 225-3718-806-NL ..... | 8  | 15045 .....   | 13 | 371001-NL .....  | 4  |
| 225-3719-800 .....    | 8  | 15048 .....   | 13 | 371002 .....     | 4  |
| 225-3719-804 .....    | 8  | 15052 .....   | 13 | 371056 .....     | 7  |
| 225-3719-804-NL ..... | 8  | 15055 .....   | 13 | 371056-NL .....  | 7  |
| 225-3720-800 .....    | 9  | 15070 .....   | 13 | 371804 .....     | 8  |
| 225-3720-802 .....    | 9  | 15071 .....   | 13 | 371804-NL .....  | 8  |
| 225-3720-804 .....    | 9  | 15245 .....   | 13 | 371821 .....     | 8  |
| 225-3723-800 .....    | 9  | 15248 .....   | 13 | 372007 .....     | 9  |
| 225-3724-800 .....    | 9  | 15252 .....   | 13 | 372036 .....     | 9  |
| 002291 .....          | 17 | 15255 .....   | 13 | 372900 .....     | 9  |
| 002811P .....         | 16 | 15265 .....   | 13 | 372904 .....     | 9  |
| 002812 .....          | 16 | 15266 .....   | 13 | 373600 .....     | 7  |

## ITEM NUMBERS - CONT'D

|                 |    |                 |    |
|-----------------|----|-----------------|----|
| 374000 .....    | 8  | 453600-NL ..... | 5  |
| 375204 .....    | 9  | 454000 .....    | 5  |
| 384000 .....    | 10 | 454204 .....    | 5  |
| 384000-NL ..... | 10 | 455300 .....    | 5  |
| 384001-NL ..... | 10 | 455417D .....   | 5  |
| 384004 .....    | 10 | 455477D .....   | 5  |
| 384004-NL ..... | 10 | 480804 .....    | 6  |
| 384006 .....    | 10 | 482802 .....    | 6  |
| 384006-NL ..... | 10 | 482802-NL ..... | 6  |
| 385000 .....    | 10 | 482804 .....    | 6  |
| 385100 .....    | 10 | 482804-NL ..... | 6  |
| 386000 .....    | 10 | 487000 .....    | 4  |
| 386000-NL ..... | 10 | 490804 .....    | 6  |
| 386001 .....    | 10 | 490804-NL ..... | 6  |
| 386003 .....    | 10 | 490900-NL ..... | 6  |
| 387000-NL ..... | 10 | 490901-NL ..... | 6  |
| 387901 .....    | 10 | 490902-NL ..... | 6  |
| 450904 .....    | 5  | 490903-NL ..... | 6  |
| 450904-SG ..... | 5  | 491000-NL ..... | 6  |
| 450906 .....    | 5  | 492804 .....    | 6  |
| 450906-NL ..... | 5  | 492804-NL ..... | 6  |
| 450907 .....    | 5  | 580011 .....    | 16 |
| 450907-NL ..... | 5  | 650904 .....    | 5  |
| 450950-NL ..... | 4  | 670001 .....    | 14 |
| 450951-NL ..... | 4  | 670002 .....    | 14 |
| 450952-NL ..... | 4  | 670003 .....    | 14 |
| 450953-NL ..... | 4  | 670012 .....    | 14 |
| 450954-NL ..... | 4  | 670804 .....    | 9  |
| 450955-NL ..... | 4  | 670904 .....    | 7  |
| 450956-NL ..... | 4  | 671004 .....    | 7  |
| 450960-NL ..... | 5  | 672004 .....    | 9  |
| 450961-NL ..... | 5  | 4509901D .....  | 5  |
| 450962-NL ..... | 5  |                 |    |
| 450963-NL ..... | 5  |                 |    |
| 450964-NL ..... | 5  |                 |    |
| 450965-NL ..... | 5  |                 |    |
| 451004 .....    | 5  |                 |    |
| 451004-NL ..... | 5  |                 |    |
| 451040-NL ..... | 4  |                 |    |
| 451041-NL ..... | 4  |                 |    |
| 451050-NL ..... | 5  |                 |    |
| 451051-NL ..... | 5  |                 |    |
| 453600 .....    | 5  |                 |    |



## Notes:



For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical International**  
1500 Eureka Park, Lower Pemberton  
Ashford, Kent, TN25 4BF  
Tel: +44 (0)845 850 0445  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical ASD, Inc.**  
6000 Nathan Lane North  
Minneapolis, MN 55442, USA  
Tel: 1-614-210-7300  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

MHYTCA-1204

**Rx** **CE** **CE**  
**ONLY** 0123 2797  
**smiths medical**

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Products referenced are CE-marked. Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex®, and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. SM196462GB-122019