

smiths medical  
portex®



## Seldinger Chest Drainage Kit

A positive alternative to large bore drains



# A POSITIVE ALTERNATIVE TO LARGE BORE DRAINS

The Portex® Seldinger Chest Drainage Kit has been specifically designed to combine the benefits of the Seldinger approach with the advantages of small bore catheters. It is particularly suited to:

- Pneumothorax – spontaneous or iatrogenic
- Effusions – benign or malignant, but excluding those with a very high protein content
- Low viscosity empyema



- 1 Scalpel N°11 to make initial skin incision
- 2 10 mL syringe to connect to blunt Tuohy introducer needle to identify when the pleural cavity has been entered
- 3 Guidewire (50 cm long) to maintain track into pleural cavity and facilitate the placement of the pleural catheter
- 4 Dilator (14F/4.66 mm) to slightly enlarge track in sub-dermal tissues and facilitate passage of the pleural catheter
- 5 Pleural catheter (12F, 3.96 mm, 30 cm long) to permit drainage
- 6 Male connector to permit connection to a chest drain system
- 7 4-way stopcock to interface between the catheter and chest drainage bottle/bag to allow intermittent sampling



## Seldinger Chest Drainage Kit

Product Code	Description	Qty Per Pack
200/910/120	Seldinger Chest Drainage Kit - Supplied Sterile	5

For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical International**  
 1500 Eureka Park, Lower Pemberton  
 Ashford Kent, TN25 4BF  
 Tel: +44 (0)845 850 0445  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced are CE marked. ©2018 Smiths Medical. All rights reserved. CT194810GB-022018

MHYTCA-1268



smiths medical