



## aerFree™ Airway Management System (AMS)

Advancing Patient Safety Through Non-Invasive Airway Support





# RESOLVING AIRWAY ISSUES STARTS WITH AVOIDING THEM

## aerFree™ AMS

### ADVANCING PATIENT SAFETY THROUGH NON-INVASIVE AIRWAY SUPPORT

aer+™ Technology is the application of negative pressure to a patient's neck above the upper airway. It simply and safely supports the patency of the upper airway during medical and surgical procedures requiring mild to moderate sedation.

- May improve patient comfort during sedation
- Allows staff to maintain greater distance from radiation field
- Effective: Reduces the number of apnea events and respiratory impairment events<sup>1</sup>
- May save clinicians time
- Easy to use: Applied to patient in seconds
- Simple: Minimal training required and can be used by all experience levels of healthcare personnel
- Hands-free chin lift: Decreases airway interventions which can disrupt workflow
- Non-invasive external airway support
- Designed for mild to moderate sedation

### ORDERING INFORMATION

#### aerFree™ Airway Management System

Product Code	Description	Units per case
30-5001	aerFree™ Airway Management System with tubing	24

### REFERENCE:

1. Kais Susan S, et al. Continuous Negative External Pressure [cNEP] Diminishes Obstructive Apneas and Improves Oxygenation During Colonoscopy. Endoscopy 2016, 48(06): 584-587

Please refer to Instructions for Use provided with the device for a full listing of the indications and precautions. aer+ and aerFree™ are trademarks of Sommetrics, Inc.

### aerFree™ Airway Management System (AMS)

is a simple to use, non-invasive external airway device that effectively and safely reduces the occurrence of upper airway obstruction during sedation.



For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Distributed by:**  
**Smiths Medical ASD, Inc.**  
 6000 Nathan Lane North  
 Minneapolis, MN 55442, USA  
 Tel: 1-614-210-7300  
 Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Manufactured by:**  
 Sommetrics, Inc.  
 2384 La Mirada Drive  
 Vista, CA 92081  
 Tel: 1-760-842-7605

MPAUC-1783

**Rx  
ONLY**

smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2017 Smiths Medical. All rights reserved. GA194070EN-102017