

CADD® Ambulatory Infusion Pump Pouches



Enhancing patient mobility during infusion therapy



CADDTM

CADD® Ambulatory Infusion Pump Pouches

For over 25 years CADD® ambulatory infusion pumps and pouches have set the standard for convenient, reliable and safe ambulatory infusion therapy by minimizing the impact that treatments have on patients' daily lives.

Form-fitted

CADD® pump pouches are designed from nylon materials to accommodate CADD® pumps, dedicated disposables and medication bags or cassette reservoirs in a form-fitted single-use or reusable option.

Variety of Sizes

CADD® pump pouches are available in several sizes to accommodate a variety of infusion therapies and help promote patient ambulation in the hospital and at home.

Each CADD® pump pouch contains:

- Hook and loop straps on the inside to hold a CADD® pump, tubing and a medication bag or cassette reservoir in place
- A buckled strap on the outside for flexibility to wear the pouch over the shoulder or around the waist
- A clear plastic panel to view the pump's screen



Order Information

50/100 mL CADD® pump pouch	21-2165-25
50/100 mL CADD® pump single-use, black pouches, 10 pack	21-2170-25
250/500 mL CADD® pump pouch with waist belt	21-2167-25
250/500 mL CADD® pump pouch with waist belt and shoulder strap	21-2168-25
250/500 mL CADD® pump single-use, black pouches, 10 pack	21-2171-25
1-liter CADD® pump backpack	21-2166-24
3-liter CADD® pump backpack	21-2169-24

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.
St. Paul, MN 55112, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Product(s) described may not be licensed or available for sale in Canada or other countries outside of the United States.

CADD and the Smiths Medical and CADD design marks are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners.

©2010 Smiths Medical family of companies. All rights reserved. 2/10 IN193000

Smiths Medical Canada Ltd.
Markham, Ontario, Canada, L3R 4Y8
Phone: 905-477-2000
Toll-Free: 1-800-387-4346

Customer and Clinical Services
U.S. 1-800-258-5361
Canada 1-800-387-4346
www.smiths-medical.com

smiths medical