

PharmGuard® Server Software with Medfusion® 4000 Syringe Pumps v1.5 / v1.6

smiths medical
medfusion™

Version Comparisons

	PharmGuard® Server Software V1.1	PharmGuard® Server Software V2.0	PharmGuard® Server Software V2.1	PharmGuard® Server Software V2.2
Medfusion® 4000 Pump Features				
One server for multiple sites / facilities	No	Yes	Yes	Yes
Multiple libraries per server (via grouping)	No	Yes	Yes	Yes
LDAP – Active directory integration of user logins	No	Yes	Yes	Yes
Number of Medfusion® 4000 pump device specific reports	26	38	38	38
Number of system / administrative reports	13	13	13	13
Capability to drill-down within reports	No	Yes	Yes	Yes
Scheduling and e-mail distribution of customized reports	No	Yes	Yes	Yes
Multi-onsite reporting (via grouping)	No	Yes	Yes	Yes
Architecture for near real time data / backbone for PharmGuard® Interoperability software	No	Yes	Yes	Yes
Microsoft® Windows® Server 2012 support	No	Yes	Yes	Yes
Microsoft® SQL Server 2012 support	No	Yes	Yes	Yes
Compatible with Microsoft® Internet Explorer® 9 or 10	No	Yes	Yes	Yes
Compatible with Microsoft® Internet Explorer® 11	No	No	Yes	Yes
Supported Medfusion® 4000 pump firmware	V1.1.2	V1.5.0	V1.5.1	V1.6
Syringe Vol. Near Empty alarm default (may be toggled at pump)	Audible	Muted	Muted	Muted
Enhanced network connectivity – auto-reconnection following drop-offs	No	No	Yes	Yes
Supported Medfusion® 4000 pump v1.5	No	Yes	Yes	Yes
Supported Medfusion® 4000 pump v1.5.1	No	Yes	Yes	Yes
Supported Medfusion® 4000 pump v1.6 with smart pump programming enabled	No	No	No	Yes
Supported Medfusion® 4000 pump v1.6 without smart pump programming enabled	No	Yes	Yes	Yes
Supported PharmGuard® Interoperability software – IHE-Interface (PharmGuard® Interoperability software)	No	V1.0	V1.0	V1.1
Supported auto-documentation in EMR via PharmGuard® Interoperability software	No	Yes	Yes	Yes
Supported smart pump programming via EMR / PharmGuard® Interoperability software	No	No	No	Yes

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Phone: 1-614-210-7300
www.smiths-medical.com

MPAUCA-1608

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. PharmGuard, Medfusion and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2017 Smiths Medical. All rights reserved. IN193738EN-072017