

PharmGuard® Administrator Medication Safety Software v4.2 and CADD® Administrator Medication Safety Software v3.3

smiths medical
medfusion™

Specifications

The PharmGuard® Administrator and CADD® Administrator medication safety software is used to create and manage medication libraries that are transferred to the CADD®-Solis pump.

- PharmGuard® Administrator software is used with CADD®-Solis v4 pumps
- CADD® Administrator software is used with CADD®-Solis v3 and VIP pumps

PharmGuard® Administrator 4.2 Software Specifications

System Requirements

Operating System	PC Operating System 32 or 64-bit <ul style="list-style-type: none">• Microsoft® Windows® 7 Professional (SP1)• Microsoft® Windows® 8 Professional
Supported Languages	English, Danish, Spanish, French, German, Italian, Japanese, Dutch, Portuguese, Swedish, Norwegian, and Finnish
RAM	1 GB for the application software and associated database
Net Framework	.Net 4.5.2
Hard disk space	100 MB (application) 350 MB (database)
Ports	<ul style="list-style-type: none">• USB port on the pc• USB A to Mini-USB cable for connection to a CADD®-Solis or CADD®-Solis VIP pump
Reports	Adobe® Reader® 9.4 or newer for viewing reports in PDF

Additional requirements for SQL Server Installations

SQL Server	Microsoft® SQL Server® 2008, 2012, or 2014 if installing the PharmGuard® Administrator software database onto your Microsoft® SQL Server® instance
------------	--

CADD® Administrator 3.3 Software Specifications

System Requirements

Operating System	PC Operating System 32 or 64-bit <ul style="list-style-type: none">• Microsoft® Windows® 7 Professional (SP1)
Supported Languages	English, Danish, Spanish, French, German, Italian, Japanese, Dutch, Portuguese, Swedish, Norwegian, and Finnish
RAM	1 GB for the application software and associated database
Net Framework	.Net 4.5.2
Hard disk space	100 MB (application) 350 MB (database)
Ports	<ul style="list-style-type: none">• USB port on the pc• USB A to Mini-USB cable for connection to a CADD®-Solis or CADD®-Solis VIP pump
Reports	Adobe® Reader® 9.4 or newer for viewing reports in PDF

Additional requirements for SQL Server Installations

SQL Server	Microsoft® SQL Server® 2008, 2012, or 2014 if installing the CADD® Administrator software database onto your Microsoft® SQL Server® instance
------------	--

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Phone: 1-614-210-7300
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. PharmGuard, CADD and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2018 Smiths Medical. All rights reserved. IN193757EN-102018

MMSPCA-0406

**Rx
ONLY**

smiths medical