



smiths medical
medex™



Mangum™ Micro Fluid Delivery System

The innovative alternative
for accurate fluid delivery



THE MANGUM™ FLUID DELIVERY SYSTEM IS NON-DEHP FORMULATION, LATEX FREE AND OFFERS THESE PRIMARY BENEFITS:

ACCURACY

- Allows reliable, accurate, and predictable fluid delivery in conjunction with a Medfusion® series syringe pump
- Small bore tubing enables low volume infusions

FREE FLOW PROTECTION


- Controlled fluid flow from the bag to the syringe or the syringe to the patient
- 90° turning radius on the 2-way stopcock is designed to help eliminate free flow to the patient

STANDARDIZATION

- Allows for standardization to syringe pumps
- Reduces cost of nurse training on multiple pump devices
- Fewer programming errors associated with use of multiple pumps

COST SAVINGS

- Allows for accurate infusions at a fraction of the cost of burette sets



Mangum™ Micro Fluid Delivery System

Order #	Description
MX448FMG	74" (188cm) Small bore tubing, with 0.2-micron filter, 6" (15cm) detachable extension, anti free-flow stopcock (approx. priming volume 3.0 mL)
MX448HLMG	69" (175cm) Small bore tubing, 6" (15cm) detachable extension, anti free-flow stopcock, non PVC (approx. priming volume 2.5 mL)
MX448F78MG	87" (221cm) Small bore tubing, with 0.2 micron filter, bonded non-removable 6" (15cm) extension, anti free-flow stopcock (approx. priming volume 2.5 mL)
MX448L60MG	67" (170cm) Small bore tubing, 6" (15cm) detachable extension, anti free-flow stopcock (approx. priming volume 2.6 mL)
MX5096020MG	62" (157cm) Microbore bore tubing, small body anti free-flow stopcock, for use with 1 cc syringe (approx. priming volume 0.5 mL)
MX536020MG	68" (173cm) Microbore tubing, 6" (15cm) detachable extension, anti free-flow stopcock, (approx. priming volume 1.3 mL)
MX536020MGIL	60" (152cm) Microbore tubing, anti free-flow stopcock (approx. priming volume 0.5 mL)



For more information visit www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
 6000 Nathan Lane North
 Minneapolis, MN 55442, USA
 Tel: 1-614-210-7300
 Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Medfusion, Mangum and the Smiths Medical and Medex design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2016 Smiths Medical. All rights reserved. IN193857EN-062016

MDBLA-1024

**Rx
ONLY**

smiths medical