

smiths medical
medex™



C-Fusor® Pressure Infusor

A clear solution for
pressure infusion,
transfusion and flushing



C-FUSOR® PRESSURE INFUSOR

The Medex® C-Fusor® reusable pressure infusor is designed for rapid set-up and bag replacement. Straps made of VELCRO® brand material hold the side panels of the infusor to the IV container and can be opened quickly when the next bag is needed.

Maximum fluid delivery

Pressure is evenly applied to the bag surface, helping send maximum volume to the IV set.

Secure hold

A reinforced loop secures the C-Fusor® to the ram's horn while a robust hanger holds the bag in place.

Unobstructed visibility

Clear panels permit checking bag volume at a glance.

Puncture resistant

Durable materials and robust construction withstand the rigors of frequent use and potential hazards.

Easy to read

Pressure gauge with simple display and clear markings can be read at a distance.


Easy to clean

Surfaces are smooth and panels lay flat for easy cleaning.

The C-Fusor® reusable pressure infusor is available in two sizes – 500 mL and 1000 mL – and accommodates 2 and 3 litre fluid bags when fitted with an extension.



Reinforced Hangers



Clear Material



Large Pressure Gauge

Ordering Information

Infusors and Extension

Code	Description
MX4805P1CZ	C-Fusor® Pressure Infusor - 500 mL
MX4810P1CZ	C-Fusor® Pressure Infusor - 1000 mL
MX4830P1CZ	C-Fusor® extension piece for 3000 mL (to be used with MX4810P1CZ)

Replacements

Code	Description
MX1821-GCZ	Pressure gauge, 0-760 mm Hg
MX1821-BCZ	Squeeze bulb, male Luerlock and 4-way stopcock

Note: The C-Fusor® pressure infusor (bag only) is warranted against defects from materials and workmanship for one year from date of sale.

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical Australasia Pty. Ltd.
Suite 2.03, 97 Waterloo Road
Macquarie Park NSW 2113
Ph AU: 1800 654 949
www.smiths-medical.com

Distributed by:
USL Medical
494 Rosebank Road
Avondale, 1026 New Zealand
Ph NZ: 0800 658 814

MDBLCA-1120

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. C-Fusor, Medex and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries ©2018 Smiths Medical. All rights reserved. MX195010GBAUS-112018

smiths medical