

Mini-Torr Plus® Non-Invasive Blood Pressure Monitor



The BCI® Mini-Torr Plus® Non-Invasive Blood Pressure Monitor provides a reliable platform for your NIBP monitoring needs. It easily adapts to meet your departmental NIBP needs, with options ranging from a lightweight roll stand or IV pole mount, to a built-in printer, pulse oximetry and Welch Allyn™ SureTemp® temperature. Its combination of reliability and flexibility makes the Mini-Torr Plus® a practical starting point for your long-term NIBP monitoring needs.

For more information visit smiths-medical.com.

PATIENT MONITORING



Mini-Torr Plus® Non-Invasive Blood Pressure Monitor

Technical Specifications

DISPLAYS, INDICATORS AND KEYS

- **SpO₂ (Optional)**
3-digit LED display, 0.40 inches (10.16 mm) high
- **Pulse Rate**
3-digit LED display, 0.30 inches (7.62 mm) high
- **Pulse Strength**
Logarithmically scaled 10-segment LED bar graph
- **NIBP**
3-digit LED display
Systolic and Diastolic Pressures – 0.40 inches (10.16 mm) high
Mean Arterial Pressure – 0.30 inches (7.62 mm) high
- **Temperature (Optional)**
3-digit LED display – 0.40 inches (10.16 mm) high

NIBP

- **Range**
Adult: 20 to 250 mmHg
Neonate: 20 to 135 mmHg
- **NIBP Accuracy**
Meets SP10-1992 AAMI standard for NIBP accuracy
- **Cuff Pressure Accuracy**
±3 mmHg or 2% (whichever is greater)
- **Determination Time**
30 to 50 seconds typical (max. 120 seconds in adult mode, 90 seconds in neonate mode)
- **Interval Times**
1, 2, 2.5, 5, 10, 15, 20, 30, 45, 60, 90 mins
- **Calibration**
Factory calibrated
- **Inflation Pressure**
Adult: 50 to 285 mmHg
Neonate: 50 to 140 mmHg
- **Default Inflation Pressure**
Adult: 175 mmHg
Neonate: 100 mmHg
- **Alarm Range – Adult**
High: 15 to 285 mmHg and OFF (1 mmHg steps)
Low: 10 to 280 mmHg and OFF (1 mmHg steps)

- **Alarm Range – Neonate**
High: 15 to 135 mmHg and OFF (1 mmHg steps)
Low: 10 to 130 mmHg and OFF (1 mmHg steps)
- **Alarm Default – Adult**
Systolic: High – 200 mmHg; Low – 90 mmHg
Diastolic: High – 105 mmHg; Low – 40 mmHg
Mean Arterial: High – 140 mmHg; Low – 50 mmHg
- **Alarm Default – Neonate**
Systolic: High – 120 mmHg; Low – 60 mmHg
Diastolic: High – 100 mmHg; Low – 40 mmHg
Mean Arterial: High – 110 mmHg; Low – 50 mmHg

PULSE RATE

- **Range**
Adult: 30 to 180 bpm (NIBP)
Neonate: 40 to 180 bpm (NIBP)
30 to 254 bpm (SpO₂)
- **Accuracy**
±2%
- **Alarm Range – Adult and Neonate**
High 30 to 250 bpm and OFF (1 bpm steps)
Low: 30 to 250 bpm and OFF (1 bpm steps)
- **Alarm Default – Adult**
High: 150 bpm
Low: 45 bpm
- **Alarm Default – Neonate**
High: 150 bpm
Low: 50 bpm

SpO₂ (OPTIONAL)

- **Range**
0 to 100% SpO₂ (functional)
- **Accuracy**
±2% counts at 70 to 100%; <70% unspecified
- **Averaging**
8 beats
- **Pulse Tone**
Volume adjustable or OFF (pitch corresponds to SpO₂ value)
- **Alarm Range – Adult and Neonate**
High: 50 - 100% and OFF (1% steps)
Low: 50 - 99% and OFF (1% steps)
- **Alarm Default – Adult**
High: OFF
Low: 85%
- **Alarm Default – Neonate**
High: OFF
Low: 80%

TEMPERATURE (OPTIONAL)

- **Range**
84°F to 108°F (28.9°C to 42.2°C)
- **Accuracy**
±0.2°F, ASTM E1112-86 (±0.1°C)

- **Calibration**
Factory calibrated

DIMENSIONS

- **Base Unit**
3.24" H x 8.5" W x 5.5" D
8.2 cm H x 21.6 cm W x 14.0 cm D
- **With Temperature**
3.24" H x 10" W x 5.5" D
8.2 cm H x 25.4 cm W x 14.0 cm D

WEIGHT

- **Base Unit**
3.5 pounds (1.6 kg)
- **With Temperature**
4.2 pounds (1.9 kg)

AC CHARGER POWER REQUIREMENTS

105 to 125 VAC, 60 Hz
90 to 110 VAC, 50 Hz (optional);
230 VAC, 50/60 Hz (optional)
Output of 24VDC @ 1.5A with 4kV isolation

BATTERY LIFE

Nickel Cadmium battery, 6 VDC, internal rechargeable, fully charged continuous-life 6 hours approximate, charge time of 4 hours

SERIAL OUTPUT

- **RS232C**
- **Data Format**
19200 baud, 1 start bit, 8 data bits, 1 stop bit, no parity

INTEGRATED PRINTER (OPTIONAL)

Text only, no graphics. Patient data log or trend tables.

ENVIRONMENTAL

- **Operating Temperature**
Base Unit: 32°F to 131°F (0°C to 55°C)
- **With Temperature**
60.8°F to 104°F (16°C to 40°C)
- **Storage Temperature**
-40°F to 158°F (-40°C to 70°C)
- **Relative Humidity** (non-condensing)
10-95% (storage)
15-95% (operating)

STANDARDS

EN60601-1 (IEC60601-1), EN60601-1-1 (IEC60601-1-1), EN60601-1-2, (IEC60601-1-2)

REGULATORY COMPLIANCE

ISO 9001/ISO 13485

CATALOG NUMBER

- **6004**
Mini-Torr Plus® NIBP Monitor

Information subject to change without notice. Please read the instructions for use supplied with the product for detailed instructions, warnings and cautions.

Q-Line to U.S.A. customers call 1-800-348-1842 for emergency and clinical application information 24 hours a day, 365 days a year.

For more information, please call Smiths Medical PM, Inc., at 262-542-3100 or 800-558-2345, or your Smiths Medical distributor.

Smiths Medical PM, Inc.

Patient Monitoring and Ventilation

N7 W22025 Johnson Drive, Waukesha, Wisconsin 53186 USA

Phone: 262-542-3100 Fax: 262-542-0718 Toll-Free: 800-558-2345

www.smiths-medical.com

Smiths Medical – a part of Smiths Group plc

The Smiths Medical Design Mark, BCI, and Mini-Torr Plus are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. ©2007 Smiths Medical family of companies. All rights reserved. BCI6004 - PM195615EN Rev 07 9/07

smiths medical