

CADD® Ambulatory Infusion Systems

Ordering Information



CADD® Tubing Sets with Flow Stop Free-Flow Protection



CADD® Extension Sets



CADD® Pumps and PharmGuard® Software



CADD® Pump Accessories













CADD® Ambulatory Tubing Sets with Flow Stop Free-Flow Protection

CADD® ambulatory tubing sets with Flow Stop protect against unregulated gravity fluid infusion by automatically activating to close the fluid path when the administration set or reservoir is removed from the pump. The system not only protects against unintentional gravity infusion but also allows gravity priming when the set is not attached to the pump.



CADD® Standard Administration Sets with Flow Stop Free-Flow Protection


- Maximum delivery rate up to 250 mL/hr with CADD®-Solis model 2110 v3 and above and CADD®-Solis VIP pumps, all models
- Maximum delivery rate up to 175 mL/hr with CADD®-Solis model 2100 (all versions) and model 2110 v1.0 and v2.0 pumps
- Maximum delivery rate up to 125 mL/hr with CADD-Legacy® pumps

	Product Code	Description	Quantity
	21-7321-24	CADD® Administration Set with female Luer, flow stop, clamp, and one-way checkvalve with male Luer (for use with a syringe), not manufactured with DEHP Length: 102 in. (259 cm) Priming Volume: 2.7 mL	12/box
	21-7322-24	CADD® Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 78 in. (198 cm) Priming Volume: 2.3 mL	12/box
	21-7323-24	CADD® Blue Striped Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 78 in. (198 cm) Priming Volume: 2.3 mL	12/box
	21-7324-24	CADD® Yellow Striped Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 123 in. (312 cm) Priming Volume: 3.2 mL	12/box
	21-7359-24	CADD® Administration Set with male Luer, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 69 in. (175 cm) Priming Volume: 2.0 mL	12/box
	21-7390-24	CADD® Administration Set with female Luer, flow stop, y-extension, one-way checkvalve with female Luer, one-way checkvalve with male Luer, and clamps (for use with a syringe), not manufactured with DEHP Length: 102 in. (259 cm) with 3 in. (7.6 cm) pigtail Priming Volume: 4.0 mL	12/box
	21-7391-24	CADD® Administration Set with bag spike, flow stop, y-extension, one-way checkvalve with female Luer, one-way checkvalve with male Luer, and clamps, not manufactured with DEHP Length: 108 in. (274 cm) with 3 in. (7.6 cm) pigtail Priming Volume: 4.2 mL	12/box
	21-7624-24	CADD® Yellow Striped Administration Set with male NRFit™ Connector with one-way checkvalve, bag spike, flow stop, and clamp, not manufactured with DEHP Length: 123 in. (312 cm) Priming Volume: 3.2 mL	12/box



CADD® Standard Administration Sets with Air-Eliminating Filters and Flow Stop Free-Flow Protection

- Maximum delivery rate up to 250 mL/hr with CADD®-Solis model 2110 v3 and above and CADD®-Solis VIP pumps, all models
- Maximum delivery rate up to 175 mL/hr with CADD®-Solis model 2100 (all versions) and model 2110 v1.0 and v2.0 pumps
- Maximum delivery rate up to 125 mL/hr with CADD-Legacy® pumps



Product Code	Description	Quantity
21-7336-24	CADD® Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 94 in. (239 cm) Priming Volume: 10.0 mL	15/box
21-7383-24	CADD® Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 93 in. (236 cm) Priming Volume: 8.0 mL	15/box
21-7339-24	CADD® Yellow Striped Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 130 in. (330 cm) Priming Volume: 6.2 mL	12/box
21-7394-24	CADD® Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 108 in. (274 cm) Priming Volume: 3.3 mL	12/box
21-7395-24	CADD® Administration Set with female Luer, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 102 in. (259 cm) Priming Volume: 3.1 mL	12/box
21-7333-24	CADD® Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, needleless injection site, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 114 in. (290 cm) Priming Volume: 15.5 mL	15/box
21-7639-24	CADD® Yellow Striped Administration Set with male NRFit™ Connector with one-way checkvalve, bag spike, flow stop, 0.2-micron air-eliminating filter, and clamp, not manufactured with DEHP Length: 130 in. (330 cm) Priming Volume: 6.2 mL	12/box



CADD® High Volume Administration Sets with Flow Stop Free-Flow Protection

- Maximum delivery rate up to 500 mL/hr with CADD®-Solis model 2110 v3 and above or CADD®-Solis VIP pumps, all models
- **NOT COMPATIBLE** with CADD-Legacy®, CADD®-Solis model 2100 (all versions) and model 2110 v1.0 and v2.0 pumps

Product Code	Description	Quantity
21-7355-24	CADD® High Volume Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 120 in. (305 cm) Priming Volume: 20 mL	15/box
21-7357-24	CADD® High Volume Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 92 in. (234 cm) Priming Volume: 11 mL	15/box
21-7375-24	CADD® High Volume Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, Luer-activated needleless injection site, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 118 in. (300 cm) Priming Volume: 18 mL	15/box
21-7381-24	CADD® High Volume Administration Set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 118 in. (300 cm) Priming Volume: 18 mL	15/box
21-7384-24	CADD® Yellow Striped High Volume Administration Set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 123 in. (312 cm) Priming Volume: 16 mL	15/box
21-7385-24	CADD® Yellow Striped High Volume Administration Set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer, not manufactured with DEHP Length: 136 in. (345 cm) Priming Volume: 23 mL	15/box
21-7684-24	CADD® Yellow Striped High Volume Administration Set with male NRFit™ Connector with one-way checkvalve, bag spike, flow stop, and clamp, not manufactured with DEHP Length: 123 in. (312 cm) Priming Volume: 16 mL	15/box
21-7685-24	CADD® Yellow Striped High Volume Administration Set with male NRFit™ Connector with one-way checkvalve, bag spike, flow stop, 0.2-micron air-eliminating filter, and clamp, not manufactured with DEHP Length: 136 in. (345 cm) Priming Volume: 23 mL	15/box



CADD™ Medication Cassette Reservoirs with Flow Stop Free-Flow Protection

Enhance patient mobility and drug security with self-contained medication cassettes. Available to be pre-filled or empty for filling to match each prescription.

- Maximum delivery rate up to 250 mL/hr with CADD®-Solis model 2110 v3 and above and CADD®-Solis VIP pumps, all models
- Maximum delivery rate up to 175 mL/hr with CADD®-Solis model 2100 (all versions) and model 2110 v1.0 and v2.0 pumps
- Maximum delivery rate up to 125 mL/hr with CADD-Legacy® pumps


Must be used with CADD® extension sets

Product Code	Description	Quantity
 21-7301-24	50-mL CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7300-24	100-mL Yellow CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7302-24	100-mL CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7308-24	250-mL Clear CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, not manufactured with DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7309-24	250-mL Yellow CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, not manufactured with DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7310-24	250-mL Blue CADD™ Medication Cassette Reservoir with flow stop, female Luer, clamp, and non-vented stopper, not manufactured with DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7600-24	100-mL Yellow CADD™ Medication Cassette Reservoir with female NRFit™ Connector, flow stop, clamp, and non-vented stopper, DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
 21-7609-24	250-mL Yellow CADD™ Medication Cassette Reservoir with female NRFit™ Connector, flow stop, clamp, and non-vented stopper, not manufactured with DEHP Length: 8 in. (20 cm) Priming Volume: 0.2 mL	12/box
24-1009-64	Cap, female, yellow	50/box
24-1004-64	Cap, male, non-vented, yellow	50/box




CADD® Extension Sets


For use with CADD™ medication cassette reservoirs to ensure system level accuracy. All integrated with anti-siphon valves for extra protection against free-flow of fluids.



Product Code	Description	Quantity
21-7040	CADD® Extension Set with male/male Luer, 0.2 micron air-eliminating filter, and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 55 in. (140 cm) Priming volume: 4.6 mL	25/box



21-7045	CADD® Extension Set with male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP	50/box
21-7060	CADD® Extension Set with male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 30 in. (76 cm) Priming Volume: 1.1 mL	12/box




21-7046	CADD® Extension Set with male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP	50/box
21-7061	CADD® Extension Set with male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 45 in. (114 cm) Priming Volume: 1.4 mL	12/box




21-7047	CADD® Extension Set With male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP	50/box
21-7062	CADD® Extension Set With male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 60 in. (152 cm) Priming Volume: 1.7 mL	12/box



21-7052	CADD® Microbore Extension Set with male Luer, 0.2 micron air-eliminating filter, clamp, and integral anti-siphon valve with male Luer, not manufactured with DEHP Length: 60 in. (152 cm) Priming Volume: 1.1 mL	50/box
---------	---	--------




21-7106	CADD® Extension Set with male/male Luer, 0.2 micron air-eliminating filter, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 60 in. (152 cm) Priming Volume: 1.8 mL	50/box
---------	---	--------




21-7092	CADD® Checkvalve Extension Set with male/male Luer, y-extension with one-way checkvalve with female Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 68 in. (173 cm) with 3 in. (7.6 cm) pigtail Priming Volume: 2.9 mL	12/box
---------	---	--------



21-7105	CADD® Yellow Striped Extension Set with male/male Luer, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 90 in. (229 cm) Priming Volume: 2.3 mL	12/box
---------	---	--------




21-7107	CADD® Yellow Striped Extension Set with male/male Luer, 0.2 micron air-eliminating filter, clamp, and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 100 in. (254 cm) Priming volume: 5.1 mL	12/box
---------	--	--------




21-7108	CADD® Yellow Striped Extension Set with male/male Luer, 0.2 micron air-eliminating filter and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 90 in. (229 cm) Priming Volume: 2.4 mL	50/box
---------	---	--------




21-7109	CADD® Blue Striped Extension Set with male/male Luer, 0.2 micron air-eliminating filter and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 90 in. (229 cm) Priming Volume: 2.4 mL	50/box
---------	---	--------



21-7115	CADD® Checkvalve Extension Set with male/male Luer, y-extension with one-way checkvalve with female Luer, 0.2 micron air-eliminating filter and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 69 in. (175 cm) with 3 in. (7.6 cm) pigtail Priming Volume: 3.2 mL	12/box
---------	---	--------



21-7605-24	CADD® Yellow Striped Extension Set with male NRFit™ Connectors, integral anti-siphon valve, and clamp, not manufactured with DEHP Length: 90 in. (229 cm) Priming Volume: 2.3 mL	12/box
------------	---	--------



21-7608-24	CADD® Yellow Striped Extension Set with male NRFit™ Connectors, 0.2-micron air-eliminating filter and integral anti-siphon valve with clamp, not manufactured with DEHP Length: 90 in. (229 cm) Priming Volume: 2.4 mL	50/box
------------	---	--------

CADD® Ambulatory Infusion Pumps and Programming Software



CADD®-Solis Ambulatory Infusion Systems

Product Code	Description
21-2111-0401*	CADD®-Solis v4.1 PCA Ambulatory Pump, English Only Supports the Wireless Communication Module
21-2112-0401*	CADD®-Solis v4.1 Ambulatory Pump with yellow keypad, English Only Supports the Wireless Communication Module



CADD®-Solis VIP Ambulatory Infusion Pump

Product Code	Description
21-2120-0105*	CADD®-Solis VIP Pump with PharmGuard® Administrator Medication Safety Software (BLUE)
21-2125-0105*	CADD®-Solis VIP Pump with Select Protocols (GREEN)
21-2127-0105*	CADD®-Solis VIP Pump with Standard Settings (BLACK)



CADD-Legacy® Ambulatory Infusion Pumps

Product Code	Description
21-6300	CADD-Legacy® PCA Pump Supports Continuous, Demand Dose, Clinician Bolus
21-6400	CADD-Legacy® 1 Pump Supports Continuous
21-6500	CADD-Legacy® PLUS Pump Supports Continuous, Intermittent



Programming Software for CADD®-Solis Ambulatory Pumps

Product Code	Description
21-2194-0402-25	PharmGuard® Administrator v4.2 Medication Safety Software (CD) For CADD®-Solis Pumps
21-2193-0402-25	PharmGuard® Point of Care v4.2 Medication Safety Software (CD) For CADD®-Solis Pumps
21-2136-0200-25	PharmGuard® SureLink® Remote Support Software - PC Direct Connect
67-2470-040100-51*	Electronic Software Update to CADD®-Solis v4.1 Pump Firmware, English Only
21-2194-0303-25	CADD®-Solis Administrator v3.3 Medication Safety Software (CD) For CADD®-Solis and CADD®-Solis VIP Pumps
21-2193-0303-25	CADD®-Solis Point of Care v3.3 Medication Safety Software (CD) For CADD®-Solis Pumps
21-2194-0301-25	CADD®-Solis Administrator v3.1 Medication Safety Software (CD) For CADD®-Solis and CADD®-Solis VIP Pumps

*Language specific

CADD® Pump Accessories



CADD®-Solis Power Sources

Power to the pump is provided by either 4 AA batteries, an AC adapter with power cord, rechargeable battery pack, or wireless communication module.

Product Code	Description
21-2160	Battery Pack, Rechargeable
21-0270	AC Adapter - Wallmount For use with CADD®-Solis pumps
21-2130	CADD®-Solis Wireless Communication Module for use with CADD®-Solis v4 pumps

Replacement Parts

Product Code	Description
21-2185-51	CADD® Pump Key For use with all CADD® pumps
21-2184-51	CADD®-Solis Battery Door
21-6216-64	CADD-Legacy® Battery Door
67-2130-51	LI-ION, 5200 MAH Replacement Battery, CADD®-Solis Wireless Communication Module
67-2474-51	Replacement Polemount Capture, CADD®-Solis Wireless Communication Module

Remote Dose Cord

Product Code	Description
21-2186-25	Remote Dose Cord (CADD®-Solis Pump)
21-6220-64	Remote Dose Cord (CADD-Legacy® Pump)



CADD®-Solis Pump Lockboxes with Full Keypad Access

Product Code	Description
21-2188-25	Large LockBox with CADD® Key Lock - Clear Compatible with CADD®-Solis pump and standard IV bags up to 500 mL
21-2189 -25	Large LockBox with CADD® Key Lock - Yellow Compatible with CADD®-Solis pump and standard IV bags up to 500 mL
21-2172-25	Large Lockbox with Key Lock - Clear Compatible with CADD®-Solis pump and standard IV bags up to 500 mL
21-2174-25	Large Lockbox with Key Lock - Yellow Compatible with CADD®-Solis pump and standard IV bags up to 500 mL



Polemount Accessories

Product Code	Description
21-6120-24	Lockable Polemount Bracket Slides into recess on the polemount bracket adapter or lockbox to allow pole-mounting of pump or pump and lockbox
21-2135-25	Polemount Bracket Adapter For use with CADD®-Solis pump and polemount bracket or lockable polemount bracket. Not for use when the Communication Module is attached to the pump
21-2183-25	Polemount Swivel For use with CADD®-Solis pump and polemount bracket adapter (and move below bracket adapter)
21-6212-24	Security Shell Adapter For use with CADD-Legacy® pump
21-6210-24	Polemount Bracket Adapter For use with CADD-Legacy® pump



Pouches

For use with all CADD® pumps

Reusable Pouches

Product Code	Description
21-2169-24	Backpack (3 L capacity) (A)
21-2166-24	Backpack (1 L capacity) (B)
21-2168-25	250-/500-mL Pump Pouch with Shoulder Strap (C)
21-2165-64	50-/100-mL Pump Pouch (D)
21-2167-25	250-/500-mL Waist Pouch (E)

Single Use Pouches

Product Code	Description
21-2170-64	Black, 50-/100-mL Pump Pouches (F)
21-2171-25	Black, 250-/500-mL Pump Pouches (G)





For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford, Kent, TN25 4BF
Tel: +44 (0)845 850 0445
www.smiths-medical.com

Smiths Medical ASD, Inc.
6000 Nathan Lane North Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. CADD, PharmGuard, CADD-Legacy and the Smiths Medical design mark are trademarks of Smiths Medical. NRFit is a registered trademark of GEDSA used with their permission. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. IN193501GB-112019

MMSPCA-0491

Rx ONLY CE 0086 CE 2797

smiths medical