

smiths medical

bringing technology to life



Cleo® 90 Infusion Set

All-in-one subcutaneous
infusion set



THE ADVANTAGES OF THE CLEO® 90 INFUSION SET ARE CLEAR



Enclosed needle safety

The enclosed needle retracts automatically into the housing after insertion, reducing the risk of accidental needlestick injuries.

All-in-one design

Its unique patented housing includes both an introducer and a cannula. The pre-assembled Cleo® 90 infusion set makes insertions simple and allows you to control the speed of insertion.

Secure

The Cleo® 90 infusion set is a system with its own needle protection device. After insertion of the cannula, the needle retracts automatically into its housing in order to reduce the risk of needlestick incidents.

Hypoallergenic dressing

The small, transparent adhesive dressing breathes yet keeps water out – and stays in place with less tape. Its patented release tab makes removal easy, too.

Flexible 360° connection

Ability to rotate the connection of the tubing 360°, which offers you even more freedom and flexibility. The Flex Attach™ 360° site connection provides on-site rotation offering greater flexibility.



A COMPLETE RANGE OF SUBCUTANEOUS INFUSION SETS ADAPTED TO EACH PERSON'S NEEDS

Cleo® 90 Infusion Set - Straight angle catheter

Item reference	Cannula length (mm)	Tubing length
21-7220-24	6	24" (61 cm)
21-7221-24	6	31" (79 cm)
21-7222-24	6	42" (107 cm)
21-7230-24	9	24" (61 cm)
21-7231-24	9	31" (79 cm)
21-7232-24	9	42" (107 cm)

Supplied sterile in cases of 10 units

Variety of sizes.
Available in needle lengths 6 and 9 mm.

Standard luer lock connection.



Twist cap to remove.



Insert the cannula. Press device firmly until you hear a click, then maintain for 5 seconds. Pull straight off.



Connect the tubing set.

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford Kent, TN25 4BF
Tel: +44 (0)845 850 0445

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Cleo, Flex Attach and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked.
©2019 Smiths Medical. All rights reserved. IN193873GB-092019

MPAUC-1585

Rx ONLY **CE 2797**

smiths medical