


acapella® duet Vibratory PEP Therapy System Improving Quality of Life





Two therapies together in one convenient device

The Portex® acapella® duet is a smart, convenient addition to the popular acapella® family. The benefits of vibratory PEP therapy have been well documented, and now medication can be delivered simultaneously, saving time and effort for both clinicians and patients. Engineered with a built-in port for aerosolised medication via a small volume nebuliser, acapella® duet allows convenient administration of bronchodilator or mucolytic treatments concurrently with vibratory PEP therapy, without the need for a tee adapter.

Convenient versatile and smart

Having two therapies — vibratory PEP and nebulised medication — together in one lightweight, compact device helps to further establish Smiths Medical as a leader in providing high-quality, innovative and cost-effective devices for respiratory care.

- Treatment time is reduced to one combined session instead of two separate sessions
- SVN port with tethered cap helps reduce risk of contamination when in resting position
- Versatile design allows stand-alone acapella® therapy, free of SVN attachment
- Meets both low and high-flow expiratory requirements
- Lightweight and streamlined body makes acapella® duet easy for patients to hold

acapella® duet features:

- Vibratory PEP therapy and nebulisation in one device saves time, effort and cost
- Easily disassembled and reassembled for effortless and thorough cleaning by boiling, autoclaving, or in the dishwasher (top rack only)
- Can be used by the patient in any position — reclining, standing or sitting

ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION	UNITS/CASE
27-9000	acapella® duet kit, includes acapella® duet, mouthpiece, Portex® SVN, oxygen tubing, collapsible flex tubing	10
27-9001	acapella® duet and mouthpiece only	10
007760	Portex® updraft medication nebuliser	50

acapella® family



For further information on the acapella® family, please visit:
www.smiths-medical.com/respiratorycare

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical International Ltd.
 Ashford, UK TN25 4BF
 Phone: +44 (0) 845 8500445

Smiths Medical ASD, Inc.
 6000 Nathan Lane North, Minneapolis, MN 55442, USA
 Tel: 1-614-210-7300
 Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex, acapella and the Smiths Medical and Portex design marks are registered trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. RE1942086B-122019

Rx ONLY smiths medical

MROCCA-0014