

General Anesthesia Airway

Intubation Product Catalogue



Magill Tubes



Reinforced Tubes



Polar Preformed Tubes



Speciality Tubes



Laryngoscopes



Accessories





Contents

Cuffed Magill Tubes	6
Portex® SACETT® (Oral)	6
Portex® Cuffed Endotracheal Tube (Oral / Nasal)	6
Portex® Cuffed Siliconised Endotracheal Tube (Oral / Nasal)	6
Portex® Cuffed Ivory Endotracheal Tube (Nasal)	7
Portex® Choice Cuffed Endotracheal Tube (Oral / Nasal)	7
Uncuffed Magill Tubes	8
Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal)	8
Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal)	8
Portex® Uncuffed Ivory Endotracheal Tube (Oral / Nasal)	9
Reinforced Tubes	10
Portex® Reinforced Endotracheal Tube (Oral / Nasal)	10
Portex® Reinforced Silicone Endotracheal Tube (Oral / Nasal)	10
Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral / Nasal)	11
Portex® Reinforced Endotracheal Tube (Oral / Nasal)	11
Portex® Silicone Reinforced Endotracheal Tube (Oral / Nasal)	11
Cuffed Polar Preformed Tubes	12
Portex® North Polar Endotracheal Tube (Oral)	12
Portex® North Polar Ivory Endotracheal Tube (Nasal)	12
Portex® South Polar Endotracheal Tube (Oral)	12
Uncuffed Polar Preformed Tubes	13
Portex® North Polar Endotracheal Tube (Oral)	13
Portex® South Polar Endotracheal Tube (Oral)	13
Portex® North Polar Polyurethane Endotracheal Tube (Nasal)	13
Speciality Tubes	14
Portex® Extra-long Croup Endotracheal Tube (Oral)	14
Portex® Cole-style Endotracheal Tube (Oral)	14
Portex® Microlaryngeal Ivory Endotracheal Tube (Oral / Nasal)	14
Portex® Neonatal and Paediatric Endotracheal Tube (Oral / Nasal)	14
Portex® Neonatal and Paediatric Endotracheal Tube (Oral / Nasal)	15
Fiber Optic Laryngoscopes	16
Callisto Timesco™ Laryngoscope Handle	16
XLED Timesco™ Laryngoscope Handle	16
Xenon Timesco™ Laryngoscope Handle	16
Callisto Timesco™ Macintosh Laryngoscope Blade	16

Callisto Timesco™ Miller Laryngoscope Blade.....	16
Callisto Timesco™ Robertshaw Laryngoscope Blade.....	17
Callisto Timesco™ Seward Laryngoscope Blade.....	17
Callisto Timesco™ Eclipse Laryngoscope Blade.....	17
Freeway Timesco™ Macintosh Laryngoscope Blade.....	17
Freeway Timesco™ Miller Laryngoscope Blade.....	17
Optima CLX Timesco™ Macintosh Laryngoscope Blade.....	17
Optima CLX Timesco™ Miller Laryngoscope Blade.....	18
Optima CLX Timesco™ Robertshaw Laryngoscope Blade.....	18
Optima CLX Timesco™ Seward Laryngoscope Blade.....	18
Optima CLX Timesco™ Light Guides.....	18
Optima Eclipse Timesco™ Eclipse Laryngoscope Blade.....	18
Sirius XL Timesco™ Eclipse Laryngoscope Blade.....	18
Sirius XL Timesco™ Macintosh Laryngoscope Blade.....	19
Sirius XL Timesco™ Miller Laryngoscope Blade.....	19
Callisto Timesco™ Laryngoscope Set.....	19
Timesco™ EMS Laryngoscope Set.....	19
Standard Light Laryngoscopes.....	20
Europa Timesco™ Laryngoscope Handle.....	20
Orion Timesco™ Laryngoscope Handle.....	20
Europa Timesco™ Macintosh Laryngoscope Blade.....	20
Europa Timesco™ Miller Laryngoscope Blade.....	20
Europa Timesco™ Robertshaw Laryngoscope Blade.....	20
Europa LED Timesco™ Macintosh Laryngoscope Blade.....	21
Europa LED Timesco™ Miller Laryngoscope Blade.....	21
Thoracic Airways.....	22
Portex® Endobronchial Tube.....	22
Laryngeal Masks.....	23
Portex® Silicone Laryngeal Mask.....	23
Accessories.....	24
Portex® Laryngotracheal Anaesthesia Kit.....	24
Portex® Guedel Airway.....	24
Portex® Guedel Airway.....	24
Portex® Berman Airway.....	24
Portex® Nasopharyngeal Airway.....	25
Portex® Intubation Stylet.....	25
Portex® Bougie.....	25

Portex® Exchange Guide.....	25
PressureEasy® Cuff Pressure Indicator	26
Portex® Cuff Inflator Pressure Gauge	26
Portex® Paediatric Endotracheal Tube Holder.....	26
Portex® RSP® Paediatric Endotracheal Tube Holder.....	26
Portex® Adult Endotracheal Tube Holder.....	26
Connectors.....	27
Index.....	28

Cuffed Magill Tubes

ORDERING INFORMATION



Portex® SACETT® (Oral) Tracheal tube with additional lumen for improved drainage of the subglottic region, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/189/060	6.0	9.1	290	23	10
100/189/065	6.5	9.7	300	23	10
100/189/070	7.0	10.4	310	30	10
100/189/075	7.5	11.3	320	30	10
100/189/080	8.0	11.8	330	30	10
100/189/085	8.5	12.3	330	30	10
100/189/090	9.0	13.0	330	30	10



Portex® Cuffed Endotracheal Tube (Oral / Nasal) Tracheal tube for long-term intubation and short-term ICU stays, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/199/050	5.0	6.9	250	17	20
100/199/055	5.5	7.6	280	17	20
100/199/060	6.0	8.2	290	23	20
100/199/065	6.5	8.9	300	23	20
100/199/070	7.0	9.6	310	30	20
100/199/075	7.5	10.3	320	30	20
100/199/080	8.0	10.9	330	30	20
100/199/085	8.5	11.6	330	30	20
100/199/090	9.0	12.3	330	30	20
100/199/095	9.5	13.0	330	30	20
100/199/100	10.0	13.7	330	34	20



Portex® Cuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/166/050	5.0	6.9	250	17	10
100/166/055	5.5	7.6	280	17	10
100/166/060	6.0	8.2	290	23	10
100/166/065	6.5	8.9	300	23	10
100/166/070	7.0	9.6	310	30	10
100/166/075	7.5	10.3	320	30	10
100/166/080	8.0	10.9	330	30	10
100/166/085	8.5	11.6	330	30	10
100/166/090	9.0	12.3	330	30	10
100/166/095	9.5	13.0	330	30	10
100/166/100	10.0	13.7	330	34	10

ORDERING INFORMATION

**Portex® Cuffed Ivory Endotracheal Tube (Nasal)**

Tracheal tube made from softer material particularly suitable for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/179/050	5.0	7.3	250	17	10
100/179/055	5.5	8.0	280	17	10
100/179/060	6.0	8.8	290	23	10
100/179/065	6.5	9.5	300	23	10
100/179/070	7.0	10.2	310	30	10
100/179/075	7.5	10.9	320	30	10
100/179/080	8.0	11.6	330	30	10
100/179/085	8.5	12.4	330	30	10
100/179/090	9.0	13.1	330	30	10

**Portex® Choice Cuffed Endotracheal Tube (Oral / Nasal)**

Cuffed tracheal tube for short-term intubation, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/150/050	5.0	6.7	265	22	10
100/150/055	5.5	7.3	295	22	10
100/150/060	6.0	8.0	300	24	10
100/150/065	6.5	8.5	310	24	10
100/150/070	7.0	9.2	320	26	10
100/150/075	7.5	10.0	325	26	10
100/150/080	8.0	10.7	325	28	10
100/150/085	8.5	11.3	325	28	10
100/150/090	9.0	12.0	325	28	10
100/150/095	9.5	12.7	325	28	10

Uncuffed Magill Tubes

ORDERING INFORMATION



Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, magill tip				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/111/020	2.0	2.9	140	10
100/111/025	2.5	3.4	150	10
100/111/030	3.0	4.2	170	10
100/111/035	3.5	4.8	190	10
100/111/040	4.0	5.4	210	10
100/111/045	4.5	6.2	230	10
100/111/050	5.0	6.9	250	10
100/111/055	5.5	7.6	280	10
100/111/060	6.0	8.2	290	10
100/111/065	6.5	8.9	300	10
100/111/070	7.0	9.6	310	10
100/111/075	7.5	10.3	320	10
100/111/080	8.0	10.9	330	10
100/111/085	8.5	11.6	330	10
100/111/090	9.0	12.3	330	10



Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/141/025	2.5	3.4	150	10
100/141/030	3.0	4.2	170	10
100/141/035	3.5	4.8	190	10
100/141/040	4.0	5.4	210	10
100/141/045	4.5	6.2	230	10
100/141/050	5.0	6.9	250	10
100/141/055	5.5	7.6	280	10
100/141/060	6.0	8.2	290	10
100/141/065	6.5	8.9	300	10
100/141/070	7.0	9.6	310	10
100/141/075	7.5	10.3	320	10
100/141/080	8.0	10.9	330	10
100/141/085	8.5	11.6	330	10
100/141/090	9.0	12.3	330	10
100/141/095	9.5	13.0	330	10

ORDERING INFORMATION



Portex® Uncuffed Ivory Endotracheal Tube (Oral / Nasal)

Tracheal tube made from softer material particularly suitable for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/105/025	2.5	3.7	150	10
100/105/030	3.0	4.4	170	10
100/105/035	3.5	5.1	190	10
100/105/040	4.0	5.9	210	10
100/105/045	4.5	6.6	230	10
100/105/050	5.0	7.3	250	10
100/105/055	5.5	8.0	280	10
100/105/060	6.0	8.8	290	10
100/105/065	6.5	9.5	300	10
100/105/070	7.0	10.2	310	10
100/105/075	7.5	10.9	320	10
100/105/080	8.0	11.6	330	10
100/105/085	8.5	12.4	330	10
100/105/090	9.0	13.1	330	10

Reinforced Tubes

ORDERING INFORMATION



Portex® Reinforced Endotracheal Tube (Oral / Nasal) Cuffed single-use tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/110/050	5.0	8.2	259	17	5
100/110/055	5.5	8.6	287	17	5
100/110/060	6.0	9.1	297	23	5
100/110/065	6.5	9.8	307	23	5
100/110/070	7.0	10.4	318	30	5
100/110/075	7.5	11.1	328	30	5
100/110/080	8.0	11.9	338	30	5
100/110/085	8.5	12.4	338	30	5
100/110/090	9.0	13.1	338	30	5
100/110/095	9.5	13.7	338	30	5



Portex® Reinforced Silicone Endotracheal Tube (Oral / Nasal) Cuffed reusable tracheal tube, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/109/050	5.0	7.3	275	2
100/109/055	5.5	8.0	275	2
100/109/060	6.0	8.7	285	2
100/109/065	6.5	9.4	295	2
100/109/070	7.0	10.0	305	2
100/109/075	7.5	10.4	315	2
100/109/080	8.0	11.0	325	2
100/109/085	8.5	11.8	325	2
100/109/090	9.0	12.4	325	2

ORDERING INFORMATION

**Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral / Nasal)**

Reusable tracheal tube with TTS™ or small diameter Aire-Cuf®, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
TTS™					
35W025	2.5	4.0	140		1
35W030	3.0	4.7	160		1
35W035	3.5	5.3	180		1
35W040	4.0	6.0	200		1
35W045	4.5	6.7	220		1
35W050	5.0	7.3	240		1
35W055	5.5	8.0	260		1
Small Diameter Aire-Cuf®, magill tip					
35W060	6.0	8.7	280	11	1
35W065	6.5	9.4	280	11	1
35W070	7.0	10.0	300	13	1
35W075	7.5	10.4	300	13	1
35W080	8.0	11.0	310	14	1
35W085	8.5	11.8	310	14	1
35W090	9.0	12.3	310	16	1
35W095	9.5	13.3	310	16	1

**Portex® Reinforced Endotracheal Tube (Oral / Nasal)**

Uncuffed single-use tracheal tube, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/113/030	3.0	4.5	175	5
100/113/035	3.5	5.3	195	5
100/113/040	4.0	6.0	215	5
100/113/045	4.5	6.6	245	5
100/113/050	5.0	8.2	259	5

**Portex® Silicone Reinforced Endotracheal Tube (Oral / Nasal)**

Uncuffed reusable tracheal tube, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/107/030	3.0	4.7	160	2
100/107/035	3.5	5.4	180	2
100/107/040	4.0	6.0	200	2
100/107/045	4.5	6.6	220	2
100/107/050	5.0	7.3	240	2

Cuffed Polar Preformed Tubes

ORDERING INFORMATION



Portex® North Polar Endotracheal Tube (Oral) Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/131/060	6.0	8.2	290	23	10
100/131/065	6.5	8.9	300	23	10
100/131/070	7.0	9.6	310	30	10
100/131/075	7.5	10.3	320	30	10
100/131/080	8.0	10.9	330	30	10
100/131/085	8.5	11.6	330	30	10
100/131/090	9.0	12.3	330	30	10



Portex® North Polar Ivory Endotracheal Tube (Nasal) Cuffed preformed tracheal tube, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/133/060	6.0	8.8	270	23	10
100/133/065	6.5	9.4	275	23	10
100/133/070	7.0	10.2	285	30	10
100/133/075	7.5	10.8	295	30	10
100/133/080	8.0	11.6	305	30	10



Portex® South Polar Endotracheal Tube (Oral) Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/136/050	5.0	6.9	240	17	10
100/136/055	5.5	7.6	252	17	10
100/136/060	6.0	8.2	265	23	10
100/136/065	6.5	8.9	280	23	10
100/136/070	7.0	9.6	300	30	10
100/136/075	7.5	10.3	315	30	10
100/136/080	8.0	10.9	330	30	10
100/136/085	8.5	11.6	345	30	10
100/136/090	9.0	12.3	360	30	10

Uncuffed Polar Preformed Tubes

ORDERING INFORMATION



Portex® North Polar Endotracheal Tube (Oral)

Uncuffed preformed tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/132/030	3.0	4.2	205	10
100/132/035	3.5	4.8	225	10
100/132/040	4.0	5.4	245	10
100/132/045	4.5	6.2	270	10
100/132/050	5.0	6.8	285	10
100/132/055	5.5	7.4	300	10
100/132/060	6.0	8.2	315	10
100/132/065	6.5	8.8	330	10
100/132/070	7.0	9.6	347	10



Portex® South Polar Endotracheal Tube (Oral)

Uncuffed preformed tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/134/030	3.0	4.2	170	10
100/134/035	3.5	4.8	175	10
100/134/040	4.0	5.4	195	10
100/134/045	4.5	6.2	215	10
100/134/050	5.0	6.8	220	10
100/134/055	5.5	7.5	230	10
100/134/060	6.0	8.2	240	10
100/134/065	6.5	8.8	250	10
100/134/070	7.0	9.6	275	10
100/134/075	7.5	10.2	290	10



Portex® North Polar Polyurethane Endotracheal Tube (Nasal)

Uncuffed preformed tracheal tube, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/130/060	6.0	8.2	290	10
100/130/065	6.5	8.9	300	10
100/130/070	7.0	9.6	310	10
100/130/075	7.5	10.3	320	10
100/130/080	8.0	10.9	330	10

Speciality Tubes

ORDERING INFORMATION



Portex® Extra-long Croup Endotracheal Tube (Oral)				
Tracheal tube to help intubation where there is inflammation of the larynx and trachea, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/112/025	2.5	3.5	220	10
100/112/030	3.0	4.2	250	10
100/112/045	4.5	6.2	330	10



Portex® Cole-style Endotracheal Tube (Oral)					
Tracheal tube with stepped design for ease of placement and an effective seal, magill tip					
Item Number	I.D. (mm)	Proximal O.D. (mm)	Distal O.D. (mm)	Length (mm)	Units per Case
100/430/150	1.5	5.5	3.0	130	10
100/430/200	2.0	5.8	3.5	130	10
100/430/250	2.5	6.5	4.0	135	10
100/430/300	3.0	6.5	4.5	135	10
100/430/350	3.5	7.5	5.0	148	10
100/430/400	4.0	7.5	5.5	148	10



Portex® Microlaryngeal Ivory Endotracheal Tube (Oral / Nasal)					
Tracheal tube for use in microlaryngeal procedures, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/190/050	5.0	7.2	315	25	10



Portex® Neonatal and Paediatric Endotracheal Tube (Oral / Nasal)				
Tracheal tube with 15mm side port connector, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Tracheal tube				
H4021S	2.5	3.7	155	10
H4022S	3.0	4.4	175	10
H4023S	3.5	4.8	196	10
H4024S	4.0	5.5	216	10
H4025S	4.5	6.1	236	10
H4026S	5.0	6.9	254	10
Tracheal tube with tube holder				
H4003S	3.5	4.8	196	10
H4004S	4.0	5.5	216	10
H4005S	4.5	6.1	236	10
H4006S	5.0	6.9	254	10

ORDERING INFORMATION

**Portex® Neonatal and Paediatric Endotracheal Tube (Oral / Nasal)**

Tracheal tube with low dead space 15 mm connector and tube holder, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Tracheal tube				
H4021PS	2.5	3.7	155	10
H4022PS	3.0	4.4	175	10
H4023PS	3.5	4.8	196	10
Tracheal tube with tube holder				
H4003PS	3.5	4.8	196	10

Fiber Optic Laryngoscopes

ORDERING INFORMATION



Callisto Timesco™ Laryngoscope Handle Single-use LED fiber optic handle with or without preloaded batteries				
Item Number	Item Description	Width (mm)	Height (mm)	Units per Case
DS.3940.470.10	Medium C Dry Cell Handle	33.0	139	25
DS.3940.480.10	Medium Preloaded Handle	33.0	139	25
DS.3940.460.10	Compact Preloaded Handle	29.0	120	25



XLED Timesco™ Laryngoscope Handle Reusable LED fiber optic handle with knurled design				
Item Number	Item Description	Width (mm)	Height (mm)	Units per Case
3000.500.10	Medium C Dry Cell Handle	29.0	146	1
3000.500.05	Small AA Dry Cell Handle	19.0	146	1
3000.508.10	Stubby AA Dry Cell Handle	32.0	111	1



Xenon Timesco™ Laryngoscope Handle Reusable Xenon fiber optic handle with knurled design				
Item Number	Item Description	Width (mm)	Height (mm)	Units per Case
3000.350.10	Medium C Dry Cell Handle	29.0	146	1
3000.350.05	Small AA Dry Cell Handle	19.0	146	1
3000.358.10	Stubby AA Dry Cell Handle	32.0	111	1



Callisto Timesco™ Macintosh Laryngoscope Blade Single-use metal fiber optic blade				
Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.3940.150.05	0	13.0	80	10
DS.3940.150.10	1	18.0	90	10
DS.3940.150.15	2	18.0	100	10
DS.3940.150.20	3	18.0	130	10
DS.3940.150.25	4	17.5	153	10
DS.3940.150.30	5	16.2	174	10



Callisto Timesco™ Miller Laryngoscope Blade Single-use metal fiber optic blade				
Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.3940.185.03	00	12.0	65	10
DS.3940.185.05	0	12.0	79	10
DS.3940.185.10	1	13.0	105	10
DS.3940.185.15	2	15.0	157	10
DS.3940.185.20	3	15.5	197	10
DS.3940.185.25	4	18.0	208	10

ORDERING INFORMATION



Callisto Timesco™ Robertshaw Laryngoscope Blade

Single-use metal fiber optic blade

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.3940.235.05	0	9.5	94	10
DS.3940.235.10	1	11.5	112	10



Callisto Timesco™ Seward Laryngoscope Blade

Single-use metal fiber optic blade

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.3940.240.10	1	12.5	107	10



Callisto Timesco™ Eclipse Laryngoscope Blade

Single-use fiber optic blade with tilting tip for difficult intubations

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.3940.175.20	3	17.5	130	10
DS.3940.175.25	4	17.6	150	10



Freeway Timesco™ Macintosh Laryngoscope Blade

Single-use plastic fiber optic blade

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2965.150.10	1	30.0	111	10
2965.150.15	2	22.0	118	10
2965.150.20	3	18.0	135	10
2965.150.25	4	18.0	160	10



Freeway Timesco™ Miller Laryngoscope Blade

Single-use plastic fiber optic blade

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2965.185.05	0	16.0	104	10
2965.185.10	1	22.0	118	10
2965.185.15	2	18.0	135	10
2965.185.20	3	18.0	160	10



Optima CLX Timesco™ Macintosh Laryngoscope Blade

Reusable metal fiber optic blade with replaceable light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2970.150.05	0	15.5	79	1
2970.150.10	1	22.0	92	1
2970.150.15	2	18.0	100	1
2970.150.20	3	18.0	131	1
2970.150.25	4	25.2	155	1
2970.150.30	5	46/50	182	1

ORDERING INFORMATION



Optima CLX Timesco™ Miller Laryngoscope Blade

Reusable metal fiber optic blade with replaceable light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2970.185.03	00	11.4	66	1
2970.185.05	0	11.4	77	1
2970.185.10	1	11.7	102	1
2970.185.15	2	14.8	153	1
2970.185.20	3	14.8	197	1
2970.185.25	4	17.0	205	1



Optima CLX Timesco™ Robertshaw Laryngoscope Blade

Reusable metal fiber optic blade with replaceable light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2970.185.05	0	11.5	94	1
2970.185.10	1	13.0	114	1



Optima CLX Timesco™ Seward Laryngoscope Blade

Reusable metal fiber optic blade with replaceable light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2970.185.10	1	9.5	105	1
2970.185.15	2	15.5	130	1

Optima CLX Timesco™ Light Guides

Laryngoscope blade replacement parts

Macintosh	Miller	Robertshaw	Seward	Blade Size	Units per Case
-	2971.185.03	-	-	00	1
2971.150.05	2971.185.05	2971.235.05	2971.240.10	0	1
2971.150.10	2971.185.10	2971.235.10	2971.240.15	1	1
2971.150.15	2971.185.15	-	-	2	1
2971.150.20	2971.185.20	-	-	3	1
2971.150.25	2971.185.25	-	-	4	1



Optima Eclipse Timesco™ Eclipse Laryngoscope Blade

Reusable metal fiber optic laryngoscope blade with tilting tip for difficult intubations

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2960.175.15	2	31.0	100	1
2960.175.20	3	33.0	130	1
2960.175.25	4	38.0	158	1



Sirius XL Timesco™ Eclipse Laryngoscope Blade

Reusable metal fiber optic laryngoscope blade with tilting tip for difficult intubations

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2955.175.20	3	24.0	135	1
2955.175.25	4	24.0	160	1



ORDERING INFORMATION



Sirius XL Timesco™ Macintosh Laryngoscope Blade
Reusable metal fiber optic blade with enclosed light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2955.150.05	0	20.0	83	1
2955.150.10	1	20.0	95	1
2955.150.15	2	24.0	115	1
2955.150.20	3	24.0	135	1
2955.150.25	4	24.0	160	1
2955.150.30	5	28.0	180	1



Sirius XL Timesco™ Miller Laryngoscope Blade
Reusable metal fiber optic blade with enclosed light guide

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
2955.185.03	00	18.0	65	1
2955.185.05	0	18.0	77	1
2955.185.10	1	18.0	102	1
2955.185.15	2	19.0	150	1
2955.185.20	3	19.0	195	1
2955.185.25	4	22.0	207	1



Callisto Timesco™ Laryngoscope Set
Single-use set supplied with a pre-assembled blade and handle for convenience

Item Number	Item Description	Units per Case
DS.3950.950.30	Callisto Blade Miller 0 & Medium Preloaded Handle	10
DS.3950.950.35	Callisto Blade Miller 1 & Medium Preloaded Handle	10
DS.3950.950.15	Callisto Blade Macintosh 2 & Medium Preloaded Handle	10
DS.3950.950.20	Callisto Blade Macintosh 3 & Medium Preloaded Handle	10
DS.3950.950.25	Callisto Blade Macintosh 4 & Medium Preloaded Handle	10
DS.3950.952.30	Callisto Blade Miller 0 & Compact Preloaded Handle	10
DS.3950.952.35	Callisto Blade Miller 1 & Compact Preloaded Handle	10
DS.3950.952.15	Callisto Blade Macintosh 2 & Compact Preloaded Handle	10
DS.3950.952.20	Callisto Blade Macintosh 3 & Compact Preloaded Handle	10
DS.3950.952.25	Callisto Blade Macintosh 4 & Compact Preloaded Handle	10



Timesco™ EMS Laryngoscope Set
Set with three Macintosh blades and handle stored in an impact resistant case

Item Number	Item Description	Units per Case
2955.970.10	Sirius Macintosh 2,3 & 4 Blades & XLED Medium Handle	1
2955.970.20	Sirius Macintosh 2,3 & 4 Blades & Xenon Medium Handle	1
DS.3940.940.20	Callisto Macintosh 2,3 & 4 blades & Callisto Medium Preloaded Handle	1
DS.3940.940.25	Callisto Macintosh 2,3 & 4 blades & Callisto Compact Preloaded Handle	1
DS.3940.940.50	Callisto Macintosh 2,3 & 4 blades & XLED Medium Handle	1
DS.3940.940.60	Callisto Macintosh 2,3 & 4 blades & Xenon Medium Handle	1

Standard Light Laryngoscopes

ORDERING INFORMATION



Europa Timesco™ Laryngoscope Handle				
Single-use standard light handle				
Item Number	Item Description	Width (mm)	Height (mm)	Units per Case
DS.2940.350.10	Medium C Dry Cell Handle	34.0	138	10



Orion Timesco™ Laryngoscope Handle				
Reusable standard light handle with knurled design				
Item Number	Item Description	Width (mm)	Height (mm)	Units per Case
2905.300.10	Medium C Dry Cell Handle	29.0	147	1
2905.300.05	Small AA Dry Cell Handle	19.0	147	1
2905.308.10	Stubby AA Dry Cell Handle	32.0	112	1



Europa Timesco™ Macintosh Laryngoscope Blade				
Single-use metal standard light laryngoscope blade				
Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.2940.150.05	0	12.6	80	10
DS.2940.150.10	1	17.6	90	10
DS.2940.150.15	2	17.6	100	10
DS.2940.150.20	3	17.1	130	10
DS.2940.150.25	4	17.1	153	10
DS.2940.150.30	5	15.8	174	10



Europa Timesco™ Miller Laryngoscope Blade				
Single-use metal standard light blade				
Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.2940.185.03	00	11.6	65	10
DS.2940.185.05	0	11.6	79	10
DS.2940.185.10	1	12.6	105	10
DS.2940.185.15	2	14.6	157	10
DS.2940.185.20	3	15.1	197	10
DS.2940.185.25	4	17.6	208	10



Europa Timesco™ Robertshaw Laryngoscope Blade				
Single-use metal standard light blade				
Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.2940.235.10	1	11.1	112	10

ORDERING INFORMATION

**Europa LED Timesco™ Macintosh Laryngoscope Blade**

Single-use metal standard light blade with bright white illumination

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.2950.150.05	0	15.0	80	10
DS.2950.150.10	1	20.0	90	10
DS.2950.150.15	2	20.0	100	10
DS.2950.150.20	3	19.5	130	10
DS.2950.150.25	4	19.5	153	10
DS.2950.150.30	5	18.2	174	10

**Europa LED Timesco™ Miller Laryngoscope Blade**

Single-use metal standard light blade with bright white illumination

Item Number	Blade Size	Height (mm)	Length (mm)	Units per Case
DS.2950.185.03	0	14.0	65	10
DS.2950.185.05	0	14.0	79	10
DS.2950.185.10	1	15.0	105	10
DS.2950.185.15	2	15.0	157	10
DS.2950.185.20	3	17.0	196	10
DS.2950.185.25	4	20.0	208	10

Thoracic Airways

ORDERING INFORMATION



Portex® Endobronchial Tube						
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Tracheal Cuff Ø	Bronchial Cuff Ø	Units per Case
Left-sided endobronchial tube						
198-28L	4.6	10.3	390	23	12	1
198-32L	5.4	11.6	390	26	13	1
198-35L	6.0	12.8	400	27	18	1
198-37L	6.3	13.3	400	28	18	1
198-39L	6.6	14.2	400	29	23	1
198-41L	7.0	14.8	400	31	23	1
Right-sided endobronchial tube						
197-28R	4.6	10.3	390	23	13	1
197-32R	5.4	11.6	390	26	13	1
197-35R	6.0	12.8	400	27	21	1
197-37R	6.3	13.3	400	28	21	1
197-39R	6.6	14.2	400	29	21	1
197-41R	7.0	14.8	400	31	21	1

Laryngeal Masks

ORDERING INFORMATION



Portex® Silicone Laryngeal Mask

Single piece laryngeal mask with aperture bars

Item Number	Size (mm)	I.D. Tube (mm)	O.D. Tube (mm)	Max Cuff Volume (ml)	Units per Case
100/222/200	2.0	7.0	11.5	10	10
100/222/250	2.5	8.3	13.3	14	10
100/222/300	3.0	10.0	15.0	20	10
100/222/400	4.0	10.0	15.0	30	10
100/222/500	5.0	11.5	16.5	40	10

Accessories

ORDERING INFORMATION

Portex® Laryngotracheal Anaesthesia Kit

For patients with difficult airways

Item Number	Syringe	Catheter	Units per Case
F031-61	5cc	with Luer Lock	20

Portex® Guedel Airway

Sterile and disposable without colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/320/200	00	50	NA	50
100/320/100	0	60	NA	50
100/320/010	1	70	NA	50
100/320/020	2	90	NA	50
100/320/030	3	100	NA	50
100/320/040	4	110	NA	50

Portex® Guedel Airway

Sterile and disposable with colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/322/040	000	40	Pink	10
100/322/050	00	50	Dark Blue	10
100/322/060	0	60	Black	10
100/322/070	1	70	White	10
100/322/080	2	80	Green	10
100/322/090	3	90	Yellow	10
100/322/100	4	100	Red	10
100/322/110	5	110	Light Blue	10
100/322/120	6	120	Orange	10

Portex® Berman Airway

Sterile and disposable with two channels for suctioning



Item Number	Size	Length (mm)	Units per Case
100/324/040	000	40	10
100/324/050	00	50	10
100/324/060	0	60	10
100/324/070	1	70	10
100/324/080	2	80	10
100/324/090	3	90	10
100/324/100	4	100	10
100/324/110	5	110	10

ORDERING INFORMATION




Portex® Nasopharyngeal Airway Soft and kink resistant				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/210/060	6.0	8.8	135	10
100/210/070	7.0	10.2	155	10
100/210/080	8.0	11.6	175	10
100/210/090	9.0	13.1	195	10



Portex® Intubation Stylet Standard and extra length				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Standard length stylets				
100/120/100	2.5 – 4.5	2.0	255	10
100/120/200	5.0 – 8.0	4.0	335	10
100/120/300	8.5 – 11.0	5.0	365	10
Extra length stylets				
100/121/200	5.0 – 8.0	4.0	673	10
100/121/300	8.5 – 11.0	5.0	693	10



Portex® Bougie Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable bougies				
14-504-17	6.0 – 11.0	5.0	600	1
Single-use bougies				
100/123/515	6.0 – 11.0	5.0	700	10



Portex® Exchange Guide Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable exchange guides				
14-504-76	4.0 – 5.5	3.3	700	1
14-504-68	6.0 – 11.0	5.0	700	1
Single-use exchange guides				
100/125/005	2.0 – 3.5	1.7	500	10
100/125/010	4.0 – 5.5	3.3	700	10
100/123/015	6.0 – 11.0	5.0	700	10

ORDERING INFORMATION



PressureEasy® Cuff Pressure Indicator
Continuously monitor tracheal cuff pressure

Item Number	Item Description	Units per Case
292003	PressureEasy® Cuff Pressure Monitor	10
292004	PressureEasy® Cuff Pressure Monitor with feedback line	10



Portex® Cuff Inflator Pressure Gauge
Controlled inflation and monitoring

Item Number	Item Description	Units per Case
100/568/000	Cuff Inflator Pressure Gauge with Connecting Tube	1
100/569/000	Separate Connecting Tube	10



Portex® Paediatric Endotracheal Tube Holder

Secure fixation to avoid risk of accidental extubation, disconnection or complete displacement of the tube

Item Number	Recommended tube I.D. (mm)	Units per Case
Paediatric Tracheal Tube Holder		
100/603/025	2.5	10
100/603/030	3.0	10
100/603/035	3.5	10
100/603/040	4.0	10
100/603/045	4.5	10
Siliconized Tracheal Tube with Holder		
100/611/020	2.0	10
100/611/025	2.5	10
100/611/030	3.0	10
100/611/035	3.5	10
100/611/040	4.0	10
100/611/045	4.5	10



Portex® RSP® Paediatric Endotracheal Tube Holder

Hypoallergenic, adhesive-backed base to reduce skin irritation and allow easy fixation

Item Number	Item Description	Units per Case
H4053	Holder for 3.5 mm endotracheal tube	10
H4054	Holder for 4.0 mm endotracheal tube	10
H4055	Holder for 4.5 mm endotracheal tube	10
H4056	Holder for 5.0 mm endotracheal tube	10



Portex® Adult Endotracheal Tube Holder

Hypoallergenic, adhesive-backed base to reduce skin irritation and allow easy fixation

Item Number	Item Description	Units per Case
K1000	Holder for 4.5 – 6.5 mm endotracheal tubes	30
K2000	Holder for 7.0 – 8.5 mm endotracheal tubes	30
K3000	Holder for 9.0 – 10.0 mm endotracheal tubes	30

ORDERING INFORMATION

**Connectors**

Durable and Robust

Item Number	Item Description	Units per Case
100/250/001	Reusable swivel 15mm connector	10
100/251/001	Catheter mount adaptor female 15mm	10
100/255/150	Double swivel connector	10
100/257/000	Fibreoptic bronchoscope swivel connector	10
100/590/000	Catheter Mount without double swivel connector	20
100/594/000	Catheter Mount with double swivel connector	10

Index

NUMERICS

14-504-17	25	100/109/085	10	100/130/065	13	100/141/030	8
14-504-68	25	100/109/090	10	100/130/070	13	100/141/035	8
14-504-76	25	100/110/050	10	100/130/075	13	100/141/040	8
35W025	11	100/110/055	10	100/130/080	13	100/141/045	8
35W030	11	100/110/060	10	100/131/060	12	100/141/050	8
35W035	11	100/110/065	10	100/131/065	12	100/141/055	8
35W040	11	100/110/070	10	100/131/070	12	100/141/060	8
35W045	11	100/110/075	10	100/131/075	12	100/141/065	8
35W050	11	100/110/080	10	100/131/080	12	100/141/070	8
35W055	11	100/110/085	10	100/131/085	12	100/141/075	8
35W060	11	100/110/090	10	100/131/090	12	100/141/080	8
35W065	11	100/110/095	10	100/132/030	13	100/141/085	8
35W070	11	100/111/020	8	100/132/035	13	100/141/090	8
35W075	11	100/111/025	8	100/132/040	13	100/141/095	8
35W080	11	100/111/030	8	100/132/045	13	100/150/050	7
35W085	11	100/111/035	8	100/132/050	13	100/150/055	7
35W090	11	100/111/040	8	100/132/055	13	100/150/060	7
35W095	11	100/111/045	8	100/132/060	13	100/150/065	7
100/105/025	9	100/111/050	8	100/132/065	13	100/150/070	7
100/105/030	9	100/111/055	8	100/132/070	13	100/150/075	7
100/105/035	9	100/111/060	8	100/133/060	12	100/150/080	7
100/105/040	9	100/111/065	8	100/133/065	12	100/150/085	7
100/105/045	9	100/111/070	8	100/133/070	12	100/150/090	7
100/105/050	9	100/111/075	8	100/133/075	12	100/150/095	7
100/105/055	9	100/111/080	8	100/133/080	12	100/166/050	6
100/105/060	9	100/111/085	8	100/134/030	13	100/166/055	6
100/105/065	9	100/111/090	8	100/134/035	13	100/166/060	6
100/105/070	9	100/112/025	14	100/134/040	13	100/166/065	6
100/105/075	9	100/112/030	14	100/134/045	13	100/166/070	6
100/105/080	9	100/112/045	14	100/134/050	13	100/166/075	6
100/105/085	9	100/113/030	11	100/134/055	13	100/166/080	6
100/105/090	9	100/113/035	11	100/134/060	13	100/166/085	6
100/107/030	11	100/113/040	11	100/134/065	13	100/166/090	6
100/107/035	11	100/113/045	11	100/134/070	13	100/166/095	6
100/107/040	11	100/113/050	11	100/134/075	13	100/166/100	6
100/107/045	11	100/120/100	25	100/136/050	12	100/179/050	7
100/107/050	11	100/120/200	25	100/136/055	12	100/179/055	7
100/109/050	10	100/120/300	25	100/136/060	12	100/179/060	7
100/109/055	10	100/121/200	25	100/136/065	12	100/179/065	7
100/109/060	10	100/121/300	25	100/136/070	12	100/179/070	7
100/109/065	10	100/123/015	25	100/136/075	12	100/179/075	7
100/109/070	10	100/123/515	25	100/136/080	12	100/179/080	7
100/109/075	10	100/125/005	25	100/136/085	12	100/179/085	7
100/109/080	10	100/125/010	25	100/136/090	12	100/179/090	7
		100/130/060	13	100/141/025	8	100/189/060	6

100/189/065	6	100/322/120	24	2955.150.05	19	2971.185.25	18
100/189/070	6	100/324/040	24	2955.150.10	19	3000.350.05	16
100/189/075	6	100/324/050	24	2955.150.15	19	3000.350.10	16
100/189/080	6	100/324/060	24	2955.150.20	19	3000.358.10	16
100/189/085	6	100/324/070	24	2955.150.25	19	3000.500.05	16
100/189/090	6	100/324/080	24	2955.185.03	19	3000.500.10	16
100/190/050	14	100/324/090	24	2955.185.05	19	3000.508.10	16
100/199/050	6	100/324/100	24	2955.185.10	19	292003	26
100/199/055	6	100/324/110	24	2955.185.15	19	292004	26
100/199/060	6	100/430/150	14	2955.185.20	18 19		
100/199/065	6	100/430/200	14	2955.185.25	18 19		
100/199/070	6	100/430/250	14	2955.970.10	19	D	
100/199/075	6	100/430/300	14	2955.970.20	19	DS.2940.150.05	20
100/199/080	6	100/430/350	14	2965.150.10	17	DS.2940.150.10	20
100/199/085	6	100/430/400	14	2965.150.15	17	DS.2940.150.15	20
100/199/090	6	100/568/000	26	2965.150.20	17	DS.2940.150.20	20
100/199/095	6	100/569/000	26	2965.150.25	17	DS.2940.150.25	20
100/199/100	6	100/590/000	27	2965.185.05	17	DS.2940.150.30	20
100/210/060	25	100/594/000	27	2965.185.10	17	DS.2940.185.03	20
100/210/070	25	100/603/025	26	2965.185.15	17	DS.2940.185.05	20
100/210/080	25	100/603/030	26	2965.185.20	17	DS.2940.185.10	20
100/210/090	25	100/603/035	26	2970.150.03	18	DS.2940.185.15	20
100/222/200	23	100/603/040	26	2970.150.05	17 18	DS.2940.185.20	20
100/222/250	23	100/603/045	26	2970.150.10	18	DS.2940.185.25	20
100/222/300	23	100/611/020	26	2970.150.15	17 18	DS.2940.235.10	20
100/222/400	23	100/611/025	26	2970.150.20	17	DS.2940.350.10	20
100/222/500	23	100/611/030	26	2970.150.25	17	DS.2950.150.05	21
100/250/001	27	100/611/035	26	2970.150.25	17	DS.2950.150.10	21
100/251/001	27	100/611/040	26	2970.175.15	18	DS.2950.150.15	21
100/255/150	27	100/611/045	26	2970.175.20	18	DS.2950.150.20	21
100/257/000	27	197-28R	22	2970.175.25	18	DS.2950.150.25	21
100/320/010	24	197-32R	22	2970.185.05	18	DS.2950.150.30	21
100/320/020	24	197-35R	22	2970.185.10	18	DS.2950.185.03	21
100/320/030	24	197-37R	22	2970.185.15	18	DS.2950.185.05	21
100/320/040	24	197-39R	22	2970.185.20	18	DS.2950.185.10	21
100/320/090	24	197-41R	22	2970.185.25	18	DS.2950.185.15	21
100/320/100	24	198-28L	22	2970.185.25	18	DS.2950.185.20	21
100/320/200	24	198-32L	22	2971.150.05	18	DS.2950.185.25	21
100/322/040	24	198-35L	22	2971.150.10	18	DS.3940.150.05	16
100/322/050	24	198-37L	22	2971.150.15	18	DS.3940.150.10	16
100/322/060	24	198-39L	22	2971.150.20	18	DS.3940.150.15	16
100/322/070	24	198-41L	22	2971.150.25	18	DS.3940.150.20	16
100/322/080	24	2905.300.05	20	2971.185.03	18	DS.3940.150.25	16
100/322/100	24	2905.300.10	20	2971.185.05	18	DS.3940.150.30	16
100/322/110	24	2905.308.10	20	2971.185.10	18	DS.3940.175.20	17
				2971.185.15	18	DS.3940.175.25	17
				2971.185.20	18	DS.3940.185.03	16

DS.3940.185.05	16	H4055	26
DS.3940.185.10	16	H4056	26
DS.3940.185.15	16		
DS.3940.185.20	16	K	
DS.3940.185.25	16	K1000	26
DS.3940.235.05	17	K2000	26
DS.3940.235.10	17	K3000	26
DS.3940.240.10	17		
DS.3940.460.10	16		
DS.3940.470.10	16		
DS.3940.480.10	16		
DS.3940.940.20	19		
DS.3940.940.25	19		
DS.3940.940.50	19		
DS.3940.940.60	19		
DS.3950.950.15	19		
DS.3950.950.20	19		
DS.3950.950.25	19		
DS.3950.950.30	19		
DS.3950.950.35	19		
DS.3950.952.15	19		
DS.3950.952.20	19		
DS.3950.952.25	19		
DS.3950.952.30	19		
DS.3950.952.35	19		
F			
F031-61	24		
H			
H4003PS	15		
H4003S	14		
H4004S	14		
H4005S	14		
H4006S	14		
H4021PS	15		
H4021S	14		
H4022PS	15		
H4022S	14		
H4023PS	15		
H4023S	14		
H4024S	14		
H4025S	14		
H4026S	14		
H4053	26		
H4054	26		



For more information visit our website at www.smiths-medical.com
or call Customer Service +44 (0)845 850 0445

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford, Kent, TN25 4BF
Tel: +44 (0)845 850 0445

Smiths Medical ASD, Inc.
6000 Nathan Lane North, Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

MHYTCA-1274

Rx
ONLY

CE
2797

CE
0413

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. SACETT, AIRCARE, Aire-Cuf, Bivona, PressureEasy, SecureEasy and Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners.
©2019 Smiths Medical. All rights reserved. SM196408GB-122019