

smiths medical
portex®

Portex® Arterial Blood Sampling Kits

Safe Application. Simple Activation. Consistent Results.



PORTEX® ARTERIAL BLOOD SAMPLING KITS



SIMPLE ACTIVATION
OF SAFETY
DEVICE REDUCES
POSSIBILITY OF
NEEDLE STICK
INJURIES

Pro-Vent™ Arterial Blood Sampling Kit

Large venting filter media for quick filling to minimize patient discomfort

Needle Pro

- Engages the needle with a simple one-handed technique, freeing the other hand to attend to the site and patient
- A single safety activation mechanism to help prevent needlestick injuries

Advanced Formula Heparin

- Available with either 1 mL or 3 mL syringes and advanced formula dry lithium heparin or low heparin formula.
- Allows for, consistent and accurate analyzer readings, while reducing the potential for clots

Pro-Vent® Blood Sampling Kit Advantages

- Fully depressed syringe plunger position provides ease of aspiration in code situations
- Positive stop on 3 mL syringe eliminates plunger disengagement
- All Pro-Vent® kits with pre-attached needles are packaged in sterile pouches for use in a sterile field



More accurate analyzer readings

Safer sample transport

Minimizes patient discomfort



CONSISTENT
RESULTS WITH
UNIQUE HEPARIN
FORMULATION FOR
OVER 20 YEARS AS
INDUSTRY STANDARD



Portex® Line Draw Plus Syringes Filter-Pro® air bubble removal device

Needleless blood sampling directly from an arterial line

Dry Heparin Formulation

- Allows for consistent and accurate readings
- Available with either 1 mL Luer slip, 3 mL Luer lock, or 3 mL Luer slip syringes

Low Heparin Formulation

- Allows the full range of blood gas, electrolyte, and metabolite test results from very small blood samples without heparin interference
- Available with either 1 mL Luer slip, 3 mL Luer lock, or 3 mL Luer slip syringe

Filter-Pro® Device

- Safely expels air bubbles from the sample without exposure to blood or bloodborne pathogens
- Caps blood-filled Luer lock and Luer slip syringes for added safety during transport

Blood sampling directly from an arterial line

Easy connection to the stopcock

Ideal for use in the Labor and Delivery and NICU



Ordering Information



Pro-Vent® Plus Arterial Blood Sampling Kits					
Product Code	Description	Needle	Volume	Heparin	Units per case
4619PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	25g x 5/8"	1 mL	23.5 IU/mL	200
4629PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23g x 1.00"	1 mL	23.5 IU/mL	200
4645PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	22g x 1.00"	3 mL	23.5 IU/mL	200
4698PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23g x 1.00"	3 mL	23.5 IU/mL	200



Line Draw Plus Arterial Blood Sampling Kits					
Product Code	Description	Needle	Volume	Heparin	Units per case
4041E	Line Draw Plus Arterial Blood Sampling Kit, 1mL Luer Slip	-	1 mL	23.5 IU/mL	200
4043E	Line Draw Plus Arterial Blood Sampling Kit, 3mL Luer Slip	-	3 mL	23.5 IU/mL	200


For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford Kent, TN25 4BF
Tel: +44 (0)845 850 0445
www.smiths-medical.com

Smiths Medical ASD, Inc.
6000 Nathan Lane North Minneapolis,
MN 55442, USA
Tel: 1-614-210-7300 Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us



smiths medical