

Portex® Blue Line Ultra® Tracheostomy Tubes and Kits



EASY INSERTION AND COMFORT

PORTEX®

Portex® Blue Line Ultra® Tracheostomy Tubes

The Blue Line Ultra® Tracheostomy Tube is your one tube solution. Tubes are available separately for surgical placement or included in percutaneous tracheostomy kits and can be used in both an in-patient and out-patient setting. Simplify and standardise with the Portex® Blue Line Ultra® tracheostomy tube.



Why choose Portex® Blue Line Ultra® Tracheostomy Tubes?

- Thermosensitive material, with sufficient initial rigidity for insertion, conforms to the individual patient's upper respiratory tract at body temperature
- 15 mm connector is bonded to the tube, allowing attachment of a ventilator breathing circuit even when an inner cannula is not present
- The reusable fenestrated and non-fenestrated inner cannulae are colour coded for ease of identification
- Smooth arc shaped design with tapered tip to assist placement and removal of both tracheostomy tube and smooth inner cannula with minimal disruption to the patients body
- Clear, anatomically-shaped flange provides better access for stoma care
- Profile® Soft-Seal® cuff design with reduced bulk and smooth transition to tracheostomy tube to improve ease of insertion.

Portex® Blue Line Ultra® Suctionaid® Tracheostomy Tubes

The Blue Line Ultra Suctionaid® tube is ideal for improving patient well-being by maintaining a clean, hygienic and unobstructed airway and offering the ability to remove secretions from above the cuff. The Suctionaid® can be used as part of a VAP reduction plan.



If contaminated secretions can be removed from above the cuff of an artificial airway then it may be possible to reduce the incidence of VAP*

*Goodhill et al. Evaluation of the Blue Line Ultra "Suctionaid" tracheostomy tube. Care of the Critically Ill. 2004; 20: 72-75.




100/800

Blue Line Ultra® tube with Profile Soft-Seal Cuff

Tube Size	Pack Size	Product Code	Fenestrated Product Code
6.0mm	10 units	100/800/060	100/802/060
7.0mm	10 units	100/800/070	100/802/070
7.5mm	10 units	100/800/075	100/802/075
8.0mm	10 units	100/800/080	100/802/080
8.5mm	10 units	100/800/085	100/802/085
9.0mm	10 units	100/800/090	100/802/090
10.0mm	10 units	100/800/100	100/802/100

Inner cannula offered separately



100/810

Blue Line Ultra® Tube Kit with Inner Cannulae

Tube Size	Pack Size	Cuffed Product Code	Uncuffed Product Code	Cuffed Fenestrated Product Code	Uncuffed Fenestrated Product Code
6.0mm	1 unit	100/810/060	100/811/060	100/812/060	100/813/060
7.0mm	1 unit	100/810/070	100/811/070	100/812/070	100/813/070
7.5mm	1 unit	100/810/075	100/811/075	100/812/075	100/813/075
8.0mm	1 unit	100/810/080	100/811/080	100/812/080	100/813/080
8.5mm	1 unit	100/810/085	100/811/085	100/812/085	100/813/085
9.0mm	1 unit	100/810/090	100/811/090	100/812/090	100/813/090
10.0mm	1 unit	100/810/100	100/811/100	100/812/100	100/813/100



100/833

Blue Line Ultra® Tube Orator Speaking Valve Kit with Inner Cannulae

Tube Size	Pack Size	Cuffed Fenestrated Product Code	Uncuffed Fenestrated Product Code
7.0mm	1 unit	100/832/070	100/833/070
8.0mm	1 unit	100/832/080	100/833/080
9.0mm	1 unit	100/832/090	100/833/090

When a fenestrated kit is ordered, it will be provided with one fenestrated and one plain inner cannula



100/860

Blue Line Ultra® Suctionaid® Tube with Profile Soft-Seal Cuff

Tube Size	Pack Size	Product Code	Tube I.D.	Tube O.D.	Tube Length
6.0 mm	10 units	100/860/060	6.0 mm	9.2 mm	64.5 mm
7.0 mm	10 units	100/860/070	7.0 mm	10.5 mm	70.0 mm
7.5 mm	10 units	100/860/075	7.5 mm	11.3 mm	73.0 mm
8.0 mm	10 units	100/860/080	8.0 mm	11.9 mm	75.5 mm
8.5 mm	10 units	100/860/085	8.5 mm	12.6 mm	78.0 mm
9.0 mm	10 units	100/860/090	9.0 mm	13.3 mm	81.0mm
10.0mm	10 units	100/860/100	10.0 mm	14.0 mm	87.5 mm



100/870

Blue Line Ultra® Suctionaid® Tube with Soft-Seal Cuff and Inner Cannulae

Tube Size	Pack Size	Product Code	Tube I.D.	Tube O.D.	Tube Length
6.0mm	1 unit	100/870/060	6.0 mm	9.2 mm	64.5 mm
7.0mm	1 unit	100/870/070	7.0 mm	10.5 mm	70.0 mm
7.5mm	1 unit	100/870/075	7.5 mm	11.3 mm	73.0 mm
8.0mm	1 unit	100/870/080	8.0 mm	11.9 mm	75.5 mm
8.5mm	1 unit	100/870/085	8.5 mm	12.6 mm	78.0 mm
9.0mm	1 unit	100/870/090	9.0 mm	13.3 mm	81.0mm
10.0mm	1 unit	100/870/100	10.0 mm	14.0 mm	87.5 mm

100/850



100/851

Inner Cannulae for Blue Line Ultra® Tracheostomy Tubes

Tube Size	Pack Size	Plain Product Code	Fenestrated Product Code
6.0mm	2 units	100/850/060	100/851/060
7.0mm	2 units	100/850/070	100/851/070
7.5mm	2 units	100/850/075	100/851/075
8.0mm	2 units	100/850/080	100/851/080
8.5mm	2 units	100/850/085	100/851/085
9.0mm	2 units	100/850/090	100/851/090
10.0mm	2 units	100/850/100	100/851/100

Smiths Medical is a leading manufacturer of airway management devices. Our Portex® and Bivona® lines comprise one of healthcare's most complete tracheostomy product ranges:

- **Portex® PVC standard tubes for paediatric and adult patients**

- Blue Line® and Blue Line Ultra® Tracheostomy Tubes

- **Bivona® silicone tubes for neonates, paediatrics and adults**

- TTS™ Tubes
- Fome-Cuf® Tubes
- Aire-Cuf® Tubes

- **Portex® Percutaneous Kits**

- ULTRAPerc® percutaneous tubes and kits
- UniPerc® percutaneous and tubes

- **Customisation Services**

In addition to our wide selection of standard offerings, we provide Customisation Services for our Bivona® silicone tracheostomy tubes. Technicians in our Customisation Lab craft thousands of tubes each year to fit the unique requirements of patients, based on clinicians' orders.

A Smiths Medical representative can help you select the products best suited to the exact needs of your patients.



THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical International Ltd

1500 Eureka Park, Lower Pemberton, Ashford, Kent, TN25 4BF
Tel: +44 (0)1303 260551 Fax: +44 (0)1303 266761

www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Smiths Medical ASD, Inc.

6000 Nathan Lane North, Minneapolis, MN 55442, USA
Tel: 1-614-210-7300 Toll-Free USA: 1-800-258-5361

www.smiths-medical.com

The Portex and Smiths Medical design marks, Blue Line, Blue Line Ultra, Suctionaid, Soft-Seal, Bivona and UniPerc are trademarks of the Smiths Medical family of companies. The symbol ® indicates that the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are trade names, trademarks, and service marks of their respective owners.

©2019 Smiths Medical family of companies. All rights reserved. TR194415GL-122019

