

Thin Walled Inner Cannulae

For use with the Portex® Blue Line Ultra® Tracheostomy Tube Range



AIRWAY MANAGEMENT

New Thin Wall

Smiths Medical is pleased to supply the new Thin Walled inner cannulae option, for use with the Blue Line Ultra® range of tracheostomy tubes. The new inner cannulae, available alongside Smiths original inner cannulae, can now be purchased through the NHS supplies route.



Increased Airflow

The new Thin Walled inner cannulae provides increased airflow in sizes 7.0 to 9.0 of the Blue Line Ultra® tracheostomy tube range. In this size range the cross sectional area of the Thin Walled inner cannulae has increased, when compared to the original Blue Line Ultra® inner cannulae. For example, the cross sectional area of size 7.0 has increased by 19.38%. This increased air flow therefore reduces the work load of the patient.

Benefits

- New thin wall provides a greater internal diameter
- Increased air flow
- Compatible with the Blue Line Ultra® range of tracheostomy tubes
- Robust ring pull design maximises confidence for safe removal
- Reduced price per inner cannulae



Ordering Information

Smiths Part Number	NPC	Price	Case Size
100/856/070	FDF1417	£140.35 inc. VAT	20 units
100/856/075	FDF1418	£140.35 inc. VAT	20 units
100/856/080	FDF1419	£140.35 inc. VAT	20 units
100/856/085	FDF1420	£140.35 inc. VAT	20 units
100/856/090	FDF1421	£140.35 inc. VAT	20 units

For more details contact Smiths Medical Customer Services on +44 (0)1923 241411

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical International Ltd

1500 Eureka Park, Lower Pemberton,
Ashford, Kent, TN25 4BF
Tel: +44 (0)1233 722351
Fax: +44 (0)1233 722153

www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Blue Line Ultra, Portex and the Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S Patent and Trademark Office and certain other countries.

©2019 Smiths Medical family of companies. All rights reserved.

TR194449GB-102019

Smiths Medical ASD, Inc.

6000 Nathan Lane North
Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361

www.smiths-medical.com



smiths medical