


smiths medical  
deltec®



**Family of GRIPPER®  
Portal Access Safety Needles**

Needle options for your therapy needs



# NEEDLE OPTIONS FOR YOUR THERAPY NEEDS



By keeping clinicians and patients in mind, Smiths Medical's GRIPPER® needles offer the innovative features that have set the standard in portal access and are asked for by name.

## GRIPPER® Portal Access Safety Needles

### GRIPPER® Micro Blunt Cannula Safety Needles



A new standard in comfort and safety

- Small, low profile infusion site
- Non-sharp, blunt cannula for safe removal and disposal
- Trocar, non-coring needle that locks into place for sharps disposal

### GRIPPER PLUS® POWER P.A.C. Safety Huber Needles



For power injection

- Designed for power injections of contrast media during CT scans
- Blue tubing differentiates from non-power needles
- Clamp insert indicates maximum flow up to 5 mL/sec



**GRIPPER PLUS®**  
Safety Huber Needles



Confidence and flexibility

---

- Wide variety of needle styles to meet patient needs
- “Click” for confirmation of needle protection

FOR OVER  
20 YEARS,  
SMITHS MEDICAL  
HAS BEEN A  
PROVIDER OF  
PORTAL ACCESS  
NEEDLES

## Ordering Information



GRIPPER® Micro Blunt  
Cannula Safety Needles

GRIPPER® Micro Blunt Cannula Safety Needles				
Gauge	Length	Reorder Numbers		
		Non Y-site	NAC	
19	19 mm	21-3260-24	21-3280-24	
19	25 mm	21-3261-24	21-3281-24	
19	32 mm	21-3262-24	21-3282-24	
20	19 mm	21-3256-24	21-3276-24	
20	25 mm	21-3257-24	21-3277-24	
20	32 mm	21-3258-24	21-3278-24	
22	19 mm	21-3251-24	21-3271-24	
22	25 mm	21-3252-24	21-3272-24	
22	32 mm	21-3253-24	21-3273-24	



GRIPPER PLUS® POWER P.A.C.  
Safety Huber Needles

GRIPPER PLUS® POWER P.A.C. Safety Huber Needles				
Gauge	Length	Reorder Numbers		
		Without Y-site	NAC	
19	19 mm	21-3368-24	21-3468-24	
19	25 mm	21-3364-24	21-3469-24	
19	32 mm	21-3365-24	21-3470-24	
20	16 mm	21-3369-24	21-3464-24	
20	19 mm	21-3367-24	21-3465-24	
20	25 mm	21-3362-24	21-3466-24	
20	32 mm	21-3363-24	21-3467-24	



GRIPPER PLUS®  
Safety Huber Needles

GRIPPER PLUS® Safety Huber Needles				
Gauge	Length	Reorder Numbers		
		Without Y-site	Y-site	NAC
19	19 mm	21-2768-24	21-2868-24	21-2968-24
19	25 mm	21-2764-24	21-2869-24	21-2969-24
19	32 mm	21-2765-24	21-2870-24	21-2970-24
20	16 mm	21-2769-24	21-2864-24	21-2964-24
20	19 mm	21-2767-24	21-2865-24	21-2965-24
20	25 mm	21-2762-24	21-2866-24	21-2966-24
20	32 mm	21-2763-24	21-2867-24	21-2967-24
22	16 mm	21-2770-24	21-2860-24	21-2960-24
22	19 mm	21-2766-24	21-2861-24	21-2961-24
22	25 mm	21-2760-24	21-2862-24	21-2962-24
22	32 mm	21-2761-24	21-2863-24	21-2963-24

For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical International**  
1500 Eureka Park, Lower Pemberton  
Ashford Kent, TN25 4BF  
Tel: +44 (0)845 850 0445  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical ASD, Inc.**  
6000 Nathan Lane North  
Minneapolis, MN 55442, USA  
Tel: 1-614-210-7300  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

MMSPCA-0504



smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. GRIPPER PLUS, GRIPPER and the Smiths Medical and Deltec design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2019 Smiths Medical. All rights reserved. VA19448GB-102019