


DuraFlex® and DuraFlex® Plus Wire Reinforced Epidural Catheters



**CHOICES WITHOUT COMPROMISE
IN WIRE REINFORCED EPIDURAL CATHETERS**

DuraFlex® and DuraFlex® Plus Epidural Catheters


For clinicians who use wire reinforced catheters, DuraFlex® and DuraFlex® Plus catheters provide a choice in catheter softness and tip geometries. Multiple catheter styles ensure that you can perform procedures with the catheter you prefer.



DuraFlex®
Wire reinforced catheters with reduced stiffness



DuraFlex® Plus
Wire reinforced catheters with a stiffer base to prevent coiling and promote easy threading



Extra Strength and Flexibility

DuraFlex® and DuraFlex® Plus catheters combine a flexible wire core enclosed within a soft nylon jacket. Both catheters maintain a balance of strength and flexibility with the availability of open ended and multiport tip configurations to enable you to deliver medication effectively.

Extra Choice, Extra Features

- Available with a choice of stiffness to ease insertion and placement
- Available with multiport or open-ended tip configurations
- Rounded tips to soften placement
- Large “windows” to visualize flash when aspirating

ORDERING INFORMATION

DURAFLEX™ EPIDURAL CATHETERS

4917-19	DuraFlex® Open End Epidural Catheter	10/case
4917M-19	DuraFlex® Closed End Multiport Epidural Catheter	10/case

DURAFLEX™ PLUS EPIDURAL CATHETERS

4917S-17	DuraFlex® Plus Open End Epidural Catheter	10/case
4917SM-17	DuraFlex® Plus Closed End Multiport Epidural Catheter	10/case

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370
Phone: (800) 848-1757 Fax: (800) 621-2590
www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group.

DuraFlex, Portex and the Smiths Medical design marks are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2009 Smiths Medical family of companies.

All rights reserved. 05.2009 19349

smiths medical