

**smiths medical**

bringing technology to life



## **CADD® Infusion Systems Portfolio with the NRFit™ Connector**

Secure Connections from Pump to Patient



# CADD® PORTFOLIO WITH THE NRFIT™ CONNECTOR

## Complete Solutions For Secure Infusion Procedures

CADD®-Solis Infusion Pump

Drug protocol always displayed on Protocol Bar with selected therapy, qualifier and drug

Intuitive programming with step-by-step instructions on color-coded screen

Eliminate numeric keypress errors with scrolling soft keys


Help avoid medication errors with clearly displayed infusion delivery parameters with unit of measure

Convenient PCA dosing while ambulatory with PCA dose key or handheld remote dose cord

CADD™ 100 mL Medication Cassette Reservoir with NRFit™ Connector



CADD™ 250 mL Medication Cassette Reservoir with NRFit™ Connector



CADD® Administration Set with NRFit™ Connector

# Ordering Information



## CADD™ Medication Cassette Reservoirs with NRFit™ Connector

Description	Length	Priming Volume	Plasticizer	Product Code	Unit
100 mL Yellow CADD™ medication cassette reservoir with NRFit™ connector, Flow Stop, female Luer, clamp, and non-vented stopper	20 cm (8")	0.2 mL	DEHP	21-7600-24	12/box
250 mL Yellow CADD™ medication cassette reservoir with NRFit™ connector, Flow Stop, female Luer, clamp, and non-vented stopper	20 cm (8")	0.2 mL	TOTM	21-7609-24	12/box



## CADD® Extension Sets with NRFit™ Connector

Description	Length	Priming Volume	Plasticizer	Product Code	Unit
CADD® yellow-striped extension set with NRFit™ connector, bag spike, Flow Stop, clamp, and one-way checkvalve with male Luer	229 cm (90")	2.3 mL	TOTM	21-7605-24	12/box
CADD® yellow-striped extension set with NRFit™ connector, male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	229 cm (90")	2.4 mL	TOTM	21-7608-24	12/box



## CADD® Administration Sets with NRFit™ Connector

Description	Length	Priming Volume	Plasticizer	Product Code	Unit
CADD® yellow-striped administration set with NRFit™ connector, bag spike, Flow Stop, clamp, and one-way checkvalve with male Luer	312 cm (123")	3.2 mL	TOTM	21-7624-24	12/box
CADD® yellow-striped administration set with NRFit™ connector, bag spike, Flow Stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	330 cm (130")	6.2 mL	TOTM	21-7639-24	12/box
CADD® yellow-striped high volume administration set with NRFit™ connector, bag spike, Flow Stop, clamp, and one-way checkvalve with male Luer	312 cm (123")	16 mL	TOTM	21-7684-24*	15/box
CADD® yellow-striped high volume administration set with NRFit™ connector, bag spike, Flow Stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	345 cm (136")	23 mL	TOTM	21-7685-24*	15/box

\*Coming soon. Pending 510k clearance. Not available for sale in the United States.



## Caps

Cap Type	Connector	Product Code	Unit
Yellow female	NRFit™	24-1009-64	50/box
Yellow male non-vented	NRFit™	24-1004-64	50/box



## Adapters

Description	Product Code	Unit
Infusion adapter to lock a NRFit™ female connector on to a female Luer	24-1001-64	10/box
Filling adapter to convert female NRFit™ connector to a female Luer for filling purposes	21-7664-24	50/box



**smiths medical**

bringing technology to life

# SECURITY FROM PUMP TO PATIENT

## A complete neuraxial solution that promotes patient security

The CADD® Infusion Systems Portfolio with the NRFit™ connector is an ISO 80369-6 compliant connection system intended to enhance patient security through a reduction in the risk of misconnection, while minimizing any change to clinical technique.

**Pinpoint stability**  
**Connections in a snap**  
**Free-flow protection with a flip**  
**The feel of precision**  
**Worry-less workflow**

For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical ASD, Inc.**

6000 Nathan Lane North  
Minneapolis, MN 55442, USA  
Phone: 1-614-210-7300

[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. CADD and the Smiths Medical design mark are trademarks of Smiths Medical. NRFit is a registered trademark of GEDSA used with their permission. Products referenced are CE-marked. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2017 Smiths Medical. All rights reserved. RA194666EN-122017

**MMSPCA-0049**

**Rx  
ONLY**

**smiths medical**