

Vascular Access

Product Catalog



PORT-A-CATH®
PORT-A-CATH® II
POWER P.A.C.



PORT-A-CATH® II



P.A.S. PORT®
Implantable
Peripheral Venous
Access Systems



GRIPPER® Micro
Blunt Cannula
Non-coring Safety
Needles



GRIPPER PLUS
POWER P.A.C.®
Safety Huber Needles



Contents

Implantable Access Systems	4
PORT-A-CATH® II POWER P.A.C. Polysulfone/Titanium Implantable Venous Access Systems	4
PORT-A-CATH® II POWER P.A.C. Systems with SafeSheath® Valved Peel-Away Introducers	5
PORT-A-CATH® POWER P.A.C. Titanium Implantable Venous Access Systems	5
P.A.S. PORT® T2 POWER P.A.C. Titanium Implantable Peripheral Venous Access Systems	6
PORT-A-CATH® II Systems with SafeSheath® Valved Peel-Away Introducers	6
PORT-A-CATH® II Polysulfone/Titanium Implantable Venous Access Systems	7
PORT-A-CATH® Titanium Implantable Venous Access Systems	8
P.A.S. PORT® Implantable Peripheral Venous Access Systems	9
ProPort® Plastic Implantable Venous Access Systems	10
PORT-A-CATH® System Accessories	11
Portal Access Needles	12
GRIPPER PLUS® POWER P.A.C. Safety Huber Needles	12
GRIPPER® Micro Blunt Cannula, Non-coring Safety Needles	13
GRIPPER PLUS® Safety Huber Needles	14
GRIPPER® Huber Needles	15
PORT-A-CATH® Plastic Hub Huber Needles	16
Index	20

PORT-A-CATH® II POWER P.A.C. Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Polysulfone/Titanium Portals with PolyFlow® Polyurethane Catheters									
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit	
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4452-24 Kit	(01) 0 0610586 03267 0	each	
						21-4453-24 Tray	(01) 0 0610586 03268 7	each	
Single-lumen Portal Preassembled	6 Fr	Polyurethane	1.9	1.0	76	21-4463-24 Tray	(01) 0 0610586 03271 7	each	
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4454-24 Kit	(01) 0 0610586 03269 4	each	
						21-4455-24 Tray	(01) 0 0610586 03270 0	each	
Single-lumen Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4465-24 Tray	(01) 0 0610586 03272 4	each	
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4482-24 Kit	(01) 0 0610586 03276 2	each	
						21-4483-24 Tray	(01) 0 0610586 03277 9	each	
Single-lumen Low Profile™ Portal Preassembled	6 Fr	Polyurethane	1.9	1.0	76	21-4485-24 Tray	(01) 0 0610586 03279 3	each	
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4470-24 Kit	(01) 0 0610586 03273 1	each	
						21-4471-24 Tray	(01) 0 0610586 03274 8	each	
Single-lumen Low Profile™ Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4473-24 Tray	(01) 0 0610586 03275 5	each	
Dual-lumen Low Profile™ Portal	10 Fr	Polyurethane	3.2	1.4/1.4	76	21-8468-24 Tray	(01) 0 0610586 03285 4	each	
Dual-lumen Low Profile™ Portal	7 Fr	Polyurethane	2.2	1.0/1.0	76	21-8466-24 Tray	(01) 0 0610586 03283 0	each	

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle, vein pick, and Point-Lok® device.

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® POWER P.A.C. safety needle, 2 syringes, tunneling tool and Point-Lok® device(s).

All systems are unassembled, unless otherwise noted.

PORT-A-CATH® II POWER P.A.C. Systems with SafeSheath® Valved Peel-Away Introducers

Introducer sheath length – 13 cm Dilator sheath length – 18 cm



PORT-A-CATH® POWER P.A.C. Systems with SafeSheath® Valved Peel-Away Introducers								
	Introducer Size	Catheter Material	OD (mm)	Length (cm)	Connector	Reorder No.	UPC No.	Unit
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK® connector	21-4474-24 Tray	(01) 1 0610586 03996 6	each
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK® connector	21-4475-24 Tray	(01) 1 0610586 03997 3	each
Single-lumen Portal	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK® connector	21-4476-24 Tray	(01) 1 0610586 03998 0	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK® connector	21-4477-24 Tray	(01) 1 0610586 03999 7	each
Dual-lumen Low Profile™ Portal	7 Fr	Polyurethane	2.2	76	WING-LOCK® connector	21-8469-24 Tray	(01) 1 0610586 04000 9	each
Dual-lumen Low Profile™ Portal	10 Fr	Polyurethane	3.2	76	WING-LOCK® connector	21-8470-24 Tray	(01) 1 0610586 04001 6	each

SafeSheath® introducers are available in Smiths Medical PORT-A-CATH® system trays and also as separate sets.

PORT-A-CATH® POWER P.A.C. Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Titanium Contoured Portals with PolyFlow® Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4422-24 Kit	(01) 0 0610586 03261 8	each
						21-4423-24 Tray	(01) 0 0610586 03263 5	each
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4424-24 Kit	(01) 0 0610586 03263 2	each
						21-4425-24 Tray	(01) 0 0610586 03264 9	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4436-24 Kit	(01) 0 0610586 03265 6	each
						21-4437-24 Tray	(01) 0 0610586 03266 3	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle, vein pick, and Point-Lok® device.

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® POWER P.A.C. safety needle(s), 2 syringes, tunneling tool and Point-Lok® device(s).

All systems are unassembled, unless otherwise noted.

P.A.S. PORT® T2 POWER P.A.C. Titanium Implantable Peripheral Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Titanium Contoured Portals with PolyFlow® Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Peripheral Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4872-24 Kit	(01) 0 0610586 03280 9	each
						21-4873-24 Tray	(01) 0 0610586 03281 6	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle, vein pick, and Point-Lok® device.

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® POWER P.A.C. safety needle(s), 2 syringes, tunneling tool and Point-Lok® device(s).

All systems are unassembled, unless otherwise noted.

PORT-A-CATH® II Systems with SafeSheath® Valved Peel-Away Introducers



PORT-A-CATH® II Systems with SafeSheath® Valved Peel-Away Introducers								
Introducer sheath length – 13 cm Dilator sheath length – 18 cm								
	Introducer Size	Catheter Material	OD (mm)	Length (cm)	Connector	Reorder No.	UPC No.	Unit
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK® connector	21-4066-24 Tray	(01) 1 0610586 03990 4	each
Single-lumen Portal	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK® connector	21-4067-24 Tray	(01) 1 0610586 03991 1	each
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK® connector	21-4068-24 Tray	(01) 1 0610586 03992 8	each
Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	76	ULTRA-LOCK® connector	21-4069-24 Tray	(01)1 0610586 03993 5	each

SafeSheath® introducers are available in Smiths Medical PORT-A-CATH® system trays and also as separate sets.

PORT-A-CATH® II

Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla



Polysulfone/Titanium Portals with PolyFlow® Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4052-24 Kit	(01) 0 0610586 02352 4	each
						21-4053-24 Tray	(01) 0 0610586 02353 1	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4054-24 Kit	(01) 0 0610586 02354 8	each
						21-4055-24 Tray	(01) 0 0610586 02355 5	each
Single-lumen Portal Preassembled	6 Fr	Polyurethane	1.9	1.0	76	21-4063-24 Tray	(01) 0 0610586 02357 9	each
Single-lumen Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4065-24 Tray	(01) 0 0610586 02358 6	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4082-24 Kit	(01) 0 0610586 02385 2	each
						21-4083-24 Tray	(01) 0 0610586 02361 6	
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4070-24 Kit	(01) 0 0610586 01782 0	each
						21-4071-24 Tray	(01) 0 0610586 01784 4	each
Single-lumen Low Profile™ Portal Preassembled	6 Fr	Polyurethane	1.9	1.0	76	21-4084-24 Kit	(01) 0 0610586 02386 9	each
						21-4085-24 Tray	(01) 0 0610586 02362 3	each
Single-lumen Low Profile™ Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4073-24 Tray	(01) 0 0610586 01787 5	each
Dual-lumen Low Profile™ Portal	7 Fr	Polyurethane	2.2	1.0/1.0	76	21-8066-24 Tray	(01) 0 0610586 02070 7	each
Dual-lumen Low Profile™ Portal	10 Fr	Polyurethane	3.2	1.4/1.4	76	21-8068-24 Tray	(01) 0 0610586 02073 8	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick and Point-Lok® device(s).

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® safety needle, 2 syringes, tunneling tool and Point-Lok® device.

All systems are unassembled, unless otherwise noted.

Polysulfone/Titanium Portals with Silicone Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4050-24 Kit	(01) 0 0610586 02350 0	each
						21-4051-24 Tray	(01) 0 0610586 02351 7	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick and Point-Lok® device(s).

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® safety needle, 2 syringes, tunneling tool and Point-Lok® device.

All systems are unassembled, unless otherwise noted.

PORT-A-CATH® Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla



Titanium Contoured Portals with PolyFlow® Polyurethane Catheters

	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4024-24 Kit	(01) 0 0610586 02379 1	each
						21-4025-24 Tray	(01) 0 0610586 02380 7	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4022-24 Kit	(01) 0 0610586 02377 7	each
						21-4023-24 Tray	(01) 0 0610586 02378 4	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4036-24 Kit	(01) 0 0610586 02383 8	each
						21-4037-24 Tray	(01) 0 0610586 02384 5	each

Titanium Contoured Portals with Silicone Catheters

	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4008-24 Kit	(01) 0 0610586 02374 6	each
						21-4009-24 Tray	(01) 0 0610586 02375 3	each
Single-lumen Low Profile™ Portal	9 Fr	Silicone	2.8	1.0	76	21-4034-24 Kit	(01) 0 0610586 02381 4	each
						21-4035-24 Tray	(01) 0 0610586 02382 1	each
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4000-24 Kit	(01) 0 0610586 02371 5	each
						21-4003-24 Tray	(01) 0 0610586 02372 2	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, and Point-Lok® device(s).

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® safety needle(s), 2 syringes, tunneling tool and Point-Lok® device(s).

All systems are unassembled, unless otherwise noted.

P.A.S. PORT® Implantable Peripheral Venous Access Systems

MR Conditional up to 3-Tesla



Titanium Portal with Polyurethane Catheter								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
P.A.S. PORT® T2 Titanium Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4572-24 Kit	(01) 0 0610586 02516 0	each
						21-4573-24 Tray	(01) 0 0610586 02517 7	each



Plastic/Titanium Portal with Polyurethane Catheter								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
P.A.S. PORT® Elite Plastic/Titanium Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4590-24 Kit	(01) 0 0610586 02154 4	each
						21-4591-24 Tray	(01) 0 0610586 02156 8	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle, vein pick and Point-Lok® device.

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/ sheath introducer assembly, 90° PORT-A-CATH® needle, GRIPPER PLUS® safety needle, 2 syringes, tunneling tool and Point-Lok® device(s).

All systems are unassembled, unless otherwise noted.

ProPort® Plastic Implantable Venous Access Systems

MR Safe



Plastic Portals with PolyFlow® Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4152-24 Kit	(01) 0 0610586 01254 2	each
						21-4153-24 Tray	(01) 0 0610586 01256 6	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4154-24 Kit	(01) 0 0610586 01257 3	each
						21-4155-24 Tray	(01) 0 0610586 01259 7	each
Single-lumen Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4165-24 Tray	(01) 0 0610586 01268 9	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4182-24 Kit	(01) 0 0610586 01275 7	each
						21-4183-24 Tray	(01) 0 0610586 01277 1	each
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4170-24 Kit	(01) 0 0610586 01269 6	each
						21-4171-24 Tray	(01) 0 0610586 01271 9	each
Single-lumen Low Profile™ Portal Preassembled	8.5 Fr	Polyurethane	2.6	1.6	76	21-4172-24 Kit	(01) 0 0610586 01272 3	each
						21-4173-24 Tray	(01) 0 0610586 01274 0	each

Plastic Portals with Silicone Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	UPC No.	Unit
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4150-24 Kit	(01) 0 0610586 01251 8	each
						21-4151-24 Tray	(01) 0 0610586 01253 5	each
Single-lumen Low Profile™ Portal	9 Fr	Silicone	2.8	1.0	76	21-4186-24 Kit	(01) 0 0610586 01281 8	each
						21-4187-24 Tray	(01) 0 0610586 01283 2	each

Kits include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), and vein pick.

Trays include: Portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® straight needle, GRIPPER PLUS® safety needle, 2 syringes, and tunneling tool.

All systems are unassembled, unless otherwise noted.

PORT-A-CATH® System Accessories



Adelante® Non-Valved Peel-Away Introducers

	Sheath Length (cm)	Dilator Length (cm)	Reorder No.	UPC No.	Unit
6 French Includes two guidewires	16	19	21-2320-24	(01) 3 0610586 02511 6	5/box
7 French	16	19	21-2471-24*	(01) 3 0610586 04015 7	5/box
8.5 French	16	19	21-2470-24	(01) 3 0610586 02514 7	5/box
9 French	16	19	21-2340-24	(01) 3 0610586 02512 3	5/box
10 French	16	19	21-2472-24*	(01) 3 0610586 04016 4	5/box
11 French	16	19	21-2386-24	(01) 3 0610586 02513 0	5/box

Adelante® introducers are available in Smiths Medical PORT-A-CATH® system trays and also as separate sets.

*Subject to availability



Introducer Sets - SafeSheath® Valved Peel-Away Introducer Sets

	Sheath Length (cm)	Dilator Length (cm)	Reorder No.	UPC No.	Unit
6 French	13	18	21-2476-24*	(01) 3 0610586 04020 1	5/box
7 French	13	18	21-2477-24*	(01) 3 0610586 04021 8	5/box
8.5 French	13	18	21-2478-24*	(01) 3 0610586 04022 5	5/box
9 French	13	18	21-2479-24*	(01) 3 0610586 04023 2	5/box
10 French	13	18	21-2480-24*	(01) 3 0610586 04024 9	5/box
11 French	13	18	21-2481-24*	(01) 3 0610586 04025 6	5/box

*Subject to availability

Subcutaneous Tunneling Device

	Reorder No.	UPC No.	Unit
3.2 mm OD x 24 cm length, solid malleable stainless steel, disposable	21-0125-24	(01) 3 0610586 01319 9	5/box

Introducer sets include: 18G extra thin wall introducer needle, "J" guidewire, vessel dilator/sheath introducer assembly, 90° PORT-A-CATH® needle (except 21-2320), 2 syringes, and Point-Lok® device.

GRIPPER PLUS® POWER P.A.C. Safety Huber Needles

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



GRIPPER PLUS® POWER P.A.C. Safety Huber Needles				
Not made with natural made rubber latex, non-coring, non-siliconized, removable contoured grip, cushioned needle platform, 8" extension tubing* with or without Y-site				
Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	UPC No.	Unit
20 gauge	5/8 inch	21-3464-24	(01) 3 0610586 03363 6	12/box
	3/4 inch	21-3465-24	(01) 3 0610586 03304 3	12/box
	1 inch	21-3466-24	(01) 3 0610586 03305 0	12/box
	1 ¼ inch	21-3467-24	(01) 3 0610586 03306 7	12/box
19 gauge	3/4 inch	21-3468-24	(01) 3 0610586 03307 4	12/box
	1 inch	21-3469-24	(01) 3 0610586 03308 1	12/box
	1 ¼ inch	21-3470-24	(01) 3 0610586 03309 8	12/box
Without Y-site	Length	Reorder No.	UPC No.	Unit
20 gauge	5/8 inch	21-3369-24	(01) 3 0610586 03302 9	12/box
	3/4 inch	21-3367-24	(01) 3 0610586 03300 5	12/box
	1 inch	21-3362-24	(01) 3 0610586 03296 1	12/box
	1 ¼ inch	21-3363-24	(01) 3 0610586 03297 8	12/box
19 gauge	3/4 inch	21-3368-24	(01) 3 0610586 03301 2	12/box
	1 inch	21-3364-24	(01) 3 0610586 03298 5	12/box
	1 ¼ inch	21-3365-24	(01) 3 0610586 03299 2	12/box

*Manufactured with Non-DEHP plasticized tubing

GRIPPER® Micro Blunt Cannula, Non-coring Safety Needles

MR Conditional up to 3-Tesla



GRIPPER® Micro Blunt Cannula, Non-coring Safety Needles

Not made with natural rubber latex, non-coring, non-siliconized, cushioned infusion site, 8" extension tubing* with or without Y-site

Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	UPC No.	Unit
22 gauge	3/4 inch	21-3271-24	(01) 3 0610586 02988 6	12/box
	1 inch	21-3272-24	(01) 3 0610586 02989 3	12/box
	1¼ inch	21-3273-24	(01) 3 0610586 02990 9	12/box
20 gauge	3/4 inch	21-3276-24	(01) 3 0610586 02993 0	12/box
	1 inch	21-3277-24	(01) 3 0610586 02994 7	12/box
	1¼ inch	21-3278-24	(01) 3 0610586 02995 4	12/box
19 gauge	3/4 inch	21-3280-24	(01) 3 0610586 02997 8	12/box
	1 inch	21-3281-24	(01) 3 0610586 02998 5	12/box
	1¼ inch	21-3282-24	(01) 3 0610586 02999 2	12/box
Without Y-site	Length	Reorder No.	UPC No.	Unit
22 gauge	3/4 inch	21-3251-24	(01) 3 0610586 02974 9	12/box
	1 inch	21-3252-24	(01) 3 0610586 02975 6	12/box
	1¼ inch	21-3253-24	(01) 3 0610586 02976 3	12/box
20 gauge	3/4 inch	21-3256-24	(01) 3 0610586 02979 4	12/box
	1 inch	21-3257-24	(01) 3 0610586 02980 0	12/box
	1¼ inch	21-3258-24	(01) 3 0610586 02981 7	12/box
19 gauge	3/4 inch	21-3260-24	(01) 3 0610586 02983 1	12/box
	1 inch	21-3261-24	(01) 3 0610586 02984 8	12/box
	1¼ inch	21-3262-24	(01) 3 0610586 02985 5	12/box

*Manufactured with Non-DEHP plasticized tubing

GRIPPER PLUS® Safety Huber Needles

MR Conditional up to 3-Tesla



GRIPPER PLUS® Safety Huber Needles				
Not made with natural rubber latex, non-coring, non-siliconized, removable contoured grip, cushioned needle platform, 8" extension tubing* with or without Y-site				
Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	UPC No.	Unit
22 gauge	5/8 inch	21-2960-24	(01) 3 0610586 02768 4	12/box
	3/4 inch	21-2961-24	(01) 3 0610586 02769 1	12/box
	1 inch	21-2962-24	(01) 3 0610586 02770 7	12/box
	1 ¼ inch	21-2963-24	(01) 3 0610586 02771 4	12/box
20 gauge	5/8 inch	21-2964-24	(01) 3 0610586 02772 1	12/box
	3/4 inch	21-2965-24	(01) 3 0610586 02773 8	12/box
	1 inch	21-2966-24	(01) 3 0610586 02774 5	12/box
	1 ¼ inch	21-2967-24	(01) 3 0610586 02775 2	12/box
19 gauge	3/4 inch	21-2968-24	(01) 3 0610586 02776 9	12/box
	1 inch	21-2969-24	(01) 3 0610586 02777 6	12/box
	1 ¼ inch	21-2970-24	(01) 3 0610586 02778 3	12/box
Y-site Removable split septum injection cap on Luer-lock Y-site	Length	Reorder No.	UPC No.	Unit
22 gauge	5/8 inch	21-2860-24	(01) 3 0610586 02564 2	12/box
	3/4 inch	21-2861-24	(01) 3 0610586 02565 9	12/box
	1 inch	21-2862-24	(01) 3 0610586 02566 6	12/box
	1 ¼ inch	21-2863-24	(01) 3 0610586 02567 3	12/box
20 gauge	5/8 inch	21-2864-24	(01) 3 0610586 02568 0	12/box
	3/4 inch	21-2865-24	(01) 3 0610586 02569 7	12/box
	1 inch	21-2866-24	(01) 3 0610586 02570 3	12/box
	1 ¼ inch	21-2867-24	(01) 3 0610586 02571 0	12/box
19 gauge	3/4 inch	21-2868-24	(01) 3 0610586 02572 7	12/box
	1 inch	21-2869-24	(01) 3 0610586 02573 4	12/box
	1 ¼ inch	21-2870-24	(01) 3 0610586 02574 1	12/box
Without Y-site	Length	Reorder No.	UPC No.	Unit
22 gauge	5/8 inch	21-2770-24	(01) 3 0610586 02563 5	12/box
	3/4 inch	21-2766-24	(01) 3 0610586 02559 8	12/box
	1 inch	21-2760-24	(01) 3 0610586 02553 6	12/box
	1 ¼ inch	21-2761-24	(01) 3 0610586 02554 3	12/box
20 gauge	5/8 inch	21-2769-24	(01) 3 0610586 02562 8	12/box
	3/4 inch	21-2767-24	(01) 3 0610586 02560 4	12/box
	1 inch	21-2762-24	(01) 3 0610586 02555 0	12/box
	1 ¼ inch	21-2763-24	(01) 3 0610586 02556 7	12/box
19 gauge	3/4 inch	21-2768-24	(01) 3 0610586 02561 1	12/box
	1 inch	21-2764-24	(01) 3 0610586 02557 4	12/box
	1 ¼ inch	21-2765-24	(01) 3 0610586 02558 1	12/box

*Manufactured with Non-DEHP plasticized tubing

GRIPPER® Huber Needles

MR Conditional up to 3-Tesla



GRIPPER® Huber Needles

Not made with natural rubber latex, non-coring, non-siliconized, removable contoured grip, cushioned needle platform, 8" extension tubing* with or without Y-site

Y-site	Length	Reorder No.	UPC No.	Unit
Removable injection cap on Luer-lock Y-site				
22 gauge	5/8 inch	21-2938-24	(01) 3 0610586 02321 1	12/box
	3/4 inch	21-2939-24	(01) 3 0610586 01674 9	12/box
	1 inch	21-2940-24	(01) 3 0610586 01675 6	12/box
	1 ¼ inch	21-2941-24	(01) 3 0610586 01676 3	12/box
20 gauge	5/8 inch	21-2946-24	(01) 3 0610586 02323 5	12/box
	3/4 inch	21-2947-24	(01) 3 0610586 01677 0	12/box
	1 inch	21-2948-24	(01) 3 0610586 01678 7	12/box
	1 ¼ inch	21-2949-24	(01) 3 0610586 01679 4	12/box
19 gauge	3/4 inch	21-2955-24	(01) 3 0610586 01680 0	12/box
	1 inch	21-2956-24	(01) 3 0610586 01681 7	12/box
	1 ¼ inch	21-2957-24	(01) 3 0610586 01682 4	12/box
Without Y-site	Length	Reorder No.	UPC No.	Unit
22 gauge	5/8 inch	21-2737-24	(01) 3 0610586 02320 4	12/box
	3/4 inch	21-2733-24	(01) 3 0610586 01584 1	12/box
	1 inch	21-2714-24	(01) 3 0610586 01578 0	12/box
	1 ¼ inch	21-2715-24	(01) 3 0610586 01579 7	12/box
20 gauge	5/8 inch	21-2736-24	(01) 3 0610586 02322 8	12/box
	3/4 inch	21-2734-24	(01) 3 0610586 01585 8	12/box
	1 inch	21-2717-24	(01) 3 0610586 01580 3	12/box
	1 ¼ inch	21-2718-24	(01) 3 0610586 01581 0	12/box
19 gauge	3/4 inch	21-2735-24	(01) 3 0610586 01586 5	12/box
	1 inch	21-2720-24	(01) 3 0610586 01582 7	12/box
	1 ¼ inch	21-2721-24	(01) 3 0610586 01583 4	12/box

*Manufactured with Non-DEHP plasticized tubing

PORT-A-CATH® Plastic Hub Huber Needles

MR Conditional up to 3-Tesla



PORT-A-CATH® Plastic Hub Huber Needles				
Non-coring, non-siliconized				
Straight	Length	Reorder No.	UPC No.	Unit
22 gauge	1 inch	21-2005-24	(01) 3 0610586 02399 0	12/box
	1 ½ inch	21-2007-24	(01) 3 0610586 02401 0	12/box
20 gauge	1 inch	21-2013-24	(01) 3 0610586 02403 4	12/box
	1 ½ inch	21-2015-24	(01) 3 0610586 02405 8	12/box
19 gauge	1 inch	21-2021-24	(01) 3 0610586 02407 2	12/box
	2 inch	21-2027-24	(01) 3 0610586 02846 9	12/box
90° Bent	Length	Reorder No.	UPC No.	Unit
22 gauge	1/2 inch	21-2316-24	(01) 3 0610586 02411 9	12/box
	1 inch	21-2006-24	(01) 3 0610586 02400 3	12/box
	1 ½ inch	21-2008-24	(01) 3 0610586 02402 7	12/box
20 gauge	1/2 inch	21-2317-24	(01) 3 0610586 02412 6	12/box
	3/4 inch	21-2287-24	(01) 3 0610586 02408 9	12/box
	1 inch	21-2014-24	(01) 3 0610586 02404 1	12/box
	1 ½ inch	21-2016-24	(01) 3 0610586 02406 5	12/box
19 gauge	1 inch	21-2022-24	(01) 3 0610586 02805 6	12/box
	1 ½ inch	21-2026-24	(01) 3 0610586 02806 3	12/box

Index

Product Codes

21-0125-24	11	21-2761-24	14	21-2968-24	14
21-2005-24	16	21-2762-24	14	21-2969-24	14
21-2006-24	16	21-2763-24	14	21-2970-24	14
21-2007-24	16	21-2764-24	14	21-3251-24	13
21-2008-24	16	21-2765-24	14	21-3252-24	13
21-2013-24	16	21-2766-24	14	21-3253-24	13
21-2014-24	16	21-2767-24	14	21-3256-24	13
21-2015-24	16	21-2768-24	14	21-3257-24	13
21-2016-24	16	21-2769-24	14	21-3258-24	13
21-2021-24	16	21-2770-24	14	21-3260-24	13
21-2022-24	16	21-2860-24	14	21-3261-24	13
21-2026-24	16	21-2861-24	14	21-3262-24	13
21-2027-24	16	21-2862-24	14	21-3271-24	13
21-2287-24	16	21-2863-24	14	21-3272-24	13
21-2316-24	16	21-2864-24	14	21-3273-24	13
21-2317-24	16	21-2865-24	14	21-3276-24	13
21-2320-24	11	21-2866-24	14	21-3277-24	13
21-2340-24	11	21-2867-24	14	21-3278-24	13
21-2386-24	11	21-2868-24	14	21-3280-24	13
21-2470-24	11	21-2869-24	14	21-3281-24	13
21-2471-24	11	21-2870-24	14	21-3282-24	13
21-2472-24	11	21-2938-24	15	21-3362-24	12
21-2476-24	11	21-2939-24	15	21-3363-24	12
21-2477-24	11	21-2940-24	15	21-3364-24	12
21-2478-24	11	21-2941-24	15	21-3365-24	12
21-2479-24	11	21-2946-24	15	21-3367-24	12
21-2480-24	11	21-2947-24	15	21-3368-24	12
21-2481-24	11	21-2948-24	15	21-3369-24	12
21-2714-24	15	21-2949-24	15	21-3464-24	12
21-2715-24	15	21-2955-24	15	21-3465-24	12
21-2717-24	15	21-2956-24	15	21-3466-24	12
21-2718-24	15	21-2957-24	15	21-3467-24	12
21-2720-24	15	21-2960-24	14	21-3468-24	12
21-2721-24	15	21-2961-24	14	21-3469-24	12
21-2733-24	15	21-2962-24	14	21-3470-24	12
21-2734-24	15	21-2963-24	14	21-4000-24	8
21-2735-24	15	21-2964-24	14	21-4003-24	8
21-2736-24	15	21-2965-24	14	21-4008-24	8
21-2737-24	15	21-2966-24	14	21-4009-24	8
21-2760-24	14	21-2967-24	14	21-4022-24	8

21-4023-24	8	21-4422-24	5
21-4024-24	8	21-4423-24	5
21-4025-24	8	21-4424-24	5
21-4034-24	8	21-4425-24	5
21-4035-24	8	21-4436-24	5
21-4036-24	8	21-4437-24	5
21-4037-24	8	21-4452-24	4
21-4050-24	7	21-4453-24	4
21-4051-24	7	21-4454-24	4
21-4052-24	7	21-4455-24	4
21-4053-24	7	21-4463-24	4
21-4054-24	7	21-4465-24	4
21-4055-24	7	21-4470-24	4
21-4063-24	7	21-4471-24	4
21-4065-24	7	21-4473-24	4
21-4066-24	6	21-4474-24	5
21-4067-24	6	21-4475-24	5
21-4068-24	6	21-4476-24	5
21-4069-24	6	21-4477-24	5
21-4070-24	7	21-4482-24	4
21-4071-24	7	21-4483-24	4
21-4073-24	7	21-4485-24	4
21-4082-24	7	21-4572-24	9
21-4083-24	7	21-4573-24	9
21-4084-24	7	21-4590-24	9
21-4085-24	7	21-4591-24	9
21-4150-24	10	21-4872-24	6
21-4151-24	10	21-4873-24	6
21-4152-24	10	21-8066-24	7
21-4153-24	10	21-8068-24	7
21-4154-24	10	21-8466-24	4
21-4155-24	10	21-8468-24	4
21-4165-24	10	21-8469-24	5
21-4170-24	10	21-8470-24	5
21-4171-24	10		
21-4172-24	10		
21-4173-24	10		
21-4182-24	10		
21-4183-24	10		
21-4186-24	10		
21-4187-24	10		

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.

6000 Nathan Lane North
Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. PORT-A-CATH, P.A.S PORT, ProPort, GRIPPER, GRIPPER PLUS, Low Profile, Point-Lok, ULTRA-LOCK, WING-LOCK and the Smiths Medical and Deltec design marks are trademarks of Smiths Medical. SafeSheath is a registered trademark of Pressure Products Medical Supplies Inc. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks, or service marks of their respective owners ©2019 Smiths Medical. All rights reserved. SM190002EN-052019

MMSPCA-0614

**Rx
ONLY**

smiths medical