

General Anesthesia Airway

Intubation Product Catalogue



Magill Tubes



Reinforced Tubes



Polar Preformed Tubes



Speciality Tubes



Accessories



Contents

Cuffed Magill Tubes	6
Portex® SACETT® (Oral)	6
Portex® Cuffed Endotracheal Tube (Oral / Nasal)	6
Portex® Cuffed Siliconised Endotracheal Tube (Oral / Nasal)	6
Portex® Cuffed Ivory Endotracheal Tube (Nasal)	7
Portex® Choice Cuffed Endotracheal Tube (Oral / Nasal)	7
Uncuffed Magill Tubes	8
Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal)	8
Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal)	8
Portex® Uncuffed Ivory Endotracheal Tube (Oral / Nasal)	9
Reinforced Tubes	10
Portex® Reinforced Endotracheal Tube (Oral / Nasal)	10
Portex® Reinforced Silicone Endotracheal Tube (Oral / Nasal)	10
Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral / Nasal)	11
Portex® Silicone Reinforced Endotracheal Tube (Oral / Nasal)	11
Cuffed Polar Preformed Tubes	12
Portex® North Polar Endotracheal Tube (Oral)	12
Portex® North Polar Ivory Endotracheal Tube (Nasal)	12
Portex® South Polar Endotracheal Tube (Oral)	12
Uncuffed Polar Preformed Tubes	13
Portex® North Polar Endotracheal Tube (Oral)	13
Portex® South Polar Endotracheal Tube (Oral)	13
Portex® North Polar Polyurethane Endotracheal Tube (Nasal)	13
Speciality Tubes	14
Portex® Extra-long Croup Endotracheal Tube (Oral)	14
Portex® Cole-style Endotracheal Tube (Oral)	14
Portex® Microlaryngeal Ivory Endotracheal Tube (Oral / Nasal)	14

Thoracic Airways	15
Portex® Endobronchial Tube.....	15
Laryngeal Masks	16
Portex® Silicone Laryngeal Mask	16
Accessories	17
Portex® Laryngotracheal Anaesthesia Kit.....	17
Portex® Guedel Airway.....	17
Portex® Guedel Airway.....	17
Portex® Berman Airway.....	17
Portex® Nasopharyngeal Airway	18
Portex® Intubation Stylet	18
Portex® Bougie.....	18
Portex® Exchange Guide	18
PressureEasy® Cuff Pressure Indicator	19
Portex® Cuff Inflator Pressure Gauge	19
Portex® Paediatric Endotracheal Tube Holder.....	19
Connectors.....	19
Index	20

Cuffed Magill Tubes

ORDERING INFORMATION



Portex® SACETT® (Oral) Tracheal tube with additional lumen for improved drainage of the subglottic region, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/189/060	6.0	9.1	290	23	10
100/189/065	6.5	9.7	300	23	10
100/189/070	7.0	10.4	310	30	10
100/189/075	7.5	11.3	320	30	10
100/189/080	8.0	11.8	330	30	10
100/189/085	8.5	12.3	330	30	10
100/189/090	9.0	13.0	330	30	10



Portex® Cuffed Endotracheal Tube (Oral / Nasal) Tracheal tube for long-term intubation and short-term ICU stays, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/199/050	5.0	6.9	250	17	20
100/199/055	5.5	7.6	280	17	20
100/199/060	6.0	8.2	290	23	20
100/199/065	6.5	8.9	300	23	20
100/199/070	7.0	9.6	310	30	20
100/199/075	7.5	10.3	320	30	20
100/199/080	8.0	10.9	330	30	20
100/199/085	8.5	11.6	330	30	20
100/199/090	9.0	12.3	330	30	20
100/199/095	9.5	13.0	330	30	20
100/199/100	10.0	13.7	330	34	20



Portex® Cuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/166/050	5.0	6.9	250	17	10
100/166/055	5.5	7.6	280	17	10
100/166/060	6.0	8.2	290	23	10
100/166/065	6.5	8.9	300	23	10
100/166/070	7.0	9.6	310	30	10
100/166/075	7.5	10.3	320	30	10
100/166/080	8.0	10.9	330	30	10
100/166/085	8.5	11.6	330	30	10
100/166/090	9.0	12.3	330	30	10
100/166/095	9.5	13.0	330	30	10
100/166/100	10.0	13.7	330	34	10

ORDERING INFORMATION

**Portex® Cuffed Ivory Endotracheal Tube (Nasal)**

Tracheal tube made from softer material particularly suitable for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/179/050	5.0	7.3	250	17	10
100/179/055	5.5	8.0	280	17	10
100/179/060	6.0	8.8	290	23	10
100/179/065	6.5	9.5	300	23	10
100/179/070	7.0	10.2	310	30	10
100/179/075	7.5	10.9	320	30	10
100/179/080	8.0	11.6	330	30	10
100/179/085	8.5	12.4	330	30	10
100/179/090	9.0	13.1	330	30	10

**Portex® Choice Cuffed Endotracheal Tube (Oral / Nasal)**

Cuffed tracheal tube for short-term intubation, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/150/050	5.0	6.7	265	22	10
100/150/055	5.5	7.3	295	22	10
100/150/060	6.0	8.0	300	24	10
100/150/065	6.5	8.5	310	24	10
100/150/070	7.0	9.2	320	26	10
100/150/075	7.5	10.0	325	26	10
100/150/080	8.0	10.7	325	28	10
100/150/085	8.5	11.3	325	28	10
100/150/090	9.0	12.0	325	28	10
100/150/095	9.5	12.7	325	28	10

Uncuffed Magill Tubes

ORDERING INFORMATION



Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, magill tip				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/111/020	2.0	2.9	140	10
100/111/025	2.5	3.4	150	10
100/111/030	3.0	4.2	170	10
100/111/035	3.5	4.8	190	10
100/111/040	4.0	5.4	210	10
100/111/045	4.5	6.2	230	10
100/111/050	5.0	6.9	250	10
100/111/055	5.5	7.6	280	10
100/111/060	6.0	8.2	290	10
100/111/065	6.5	8.9	300	10
100/111/070	7.0	9.6	310	10
100/111/075	7.5	10.3	320	10
100/111/080	8.0	10.9	330	10
100/111/085	8.5	11.6	330	10
100/111/090	9.0	12.3	330	10



Portex® Uncuffed Siliconised Endotracheal Tube (Oral / Nasal) Tracheal tube coated in soft silicone, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/141/025	2.5	3.4	150	10
100/141/030	3.0	4.2	170	10
100/141/035	3.5	4.8	190	10
100/141/040	4.0	5.4	210	10
100/141/045	4.5	6.2	230	10
100/141/050	5.0	6.9	250	10
100/141/055	5.5	7.6	280	10
100/141/060	6.0	8.2	290	10
100/141/065	6.5	8.9	300	10
100/141/070	7.0	9.6	310	10
100/141/075	7.5	10.3	320	10
100/141/080	8.0	10.9	330	10
100/141/085	8.5	11.6	330	10
100/141/090	9.0	12.3	330	10
100/141/095	9.5	13.0	330	10

ORDERING INFORMATION



Portex® Uncuffed Ivory Endotracheal Tube (Oral / Nasal)

Tracheal tube made from softer material particularly suitable for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/105/025	2.5	3.7	150	10
100/105/030	3.0	4.4	170	10
100/105/035	3.5	5.1	190	10
100/105/040	4.0	5.9	210	10
100/105/045	4.5	6.6	230	10
100/105/050	5.0	7.3	250	10
100/105/055	5.5	8.0	280	10
100/105/060	6.0	8.8	290	10
100/105/065	6.5	9.5	300	10
100/105/070	7.0	10.2	310	10
100/105/075	7.5	10.9	320	10
100/105/080	8.0	11.6	330	10
100/105/085	8.5	12.4	330	10
100/105/090	9.0	13.1	330	10

Reinforced Tubes

ORDERING INFORMATION

**Portex® Reinforced Endotracheal Tube (Oral / Nasal)**

Cuffed single-use tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/110/050	5.0	8.2	259	17	5
100/110/055	5.5	8.6	287	17	5
100/110/060	6.0	9.1	297	23	5
100/110/065	6.5	9.8	307	23	5
100/110/070	7.0	10.4	318	30	5
100/110/075	7.5	11.1	328	30	5
100/110/080	8.0	11.9	338	30	5
100/110/085	8.5	12.4	338	30	5
100/110/090	9.0	13.1	338	30	5
100/110/095	9.5	13.7	338	30	5

**Portex® Reinforced Silicone Endotracheal Tube (Oral / Nasal)**

Cuffed reusable tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/109/055	5.5	8.0	275	2
100/109/060	6.0	8.7	285	2
100/109/065	6.5	9.4	295	2
100/109/070	7.0	10.0	305	2
100/109/075	7.5	10.4	315	2
100/109/080	8.0	11.0	325	2
100/109/085	8.5	11.8	325	2
100/109/090	9.0	12.4	325	2

ORDERING INFORMATION

**Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral / Nasal)**

Reusable tracheal tube with TTS™ or small diameter Aire-Cuf®, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
-------------	-----------	-----------	-------------	--------	----------------

TTS™

35W025	2.5	4.0	140		1
35W030	3.0	4.7	160		1
35W035	3.5	5.3	180		1
35W040	4.0	6.0	200		1
35W045	4.5	6.7	220		1
35W050	5.0	7.3	240		1
35W055	5.5	8.0	260		1

Small Diameter Aire-Cuf®, magill tip

35W060	6.0	8.7	280	11	1
35W065	6.5	9.4	280	11	1
35W070	7.0	10.0	300	13	1
35W075	7.5	10.4	300	13	1
35W080	8.0	11.0	310	14	1
35W085	8.5	11.8	310	14	1
35W090	9.0	12.3	310	16	1
35W095	9.5	13.3	310	16	1

**Portex® Silicone Reinforced Endotracheal Tube (Oral / Nasal)**

Uncuffed reusable tracheal tube, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
-------------	-----------	-----------	-------------	----------------

100/107/030	3.0	4.7	160	2
100/107/035	3.5	5.4	180	2
100/107/040	4.0	6.0	200	2
100/107/045	4.5	6.6	220	2
100/107/050	5.0	7.3	240	2

Cuffed Polar Preformed Tubes

ORDERING INFORMATION



Portex® North Polar Endotracheal Tube (Oral) Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/131/060	6.0	8.2	290	23	10
100/131/065	6.5	8.9	300	23	10
100/131/070	7.0	9.6	310	30	10
100/131/075	7.5	10.3	320	30	10
100/131/080	8.0	10.9	330	30	10
100/131/085	8.5	11.6	330	30	10
100/131/090	9.0	12.3	330	30	10



Portex® North Polar Ivory Endotracheal Tube (Nasal) Cuffed preformed tracheal tube, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/133/060	6.0	8.8	270	23	10
100/133/065	6.5	9.4	275	23	10
100/133/070	7.0	10.2	285	30	10
100/133/075	7.5	10.8	295	30	10
100/133/080	8.0	11.6	305	30	10



Portex® South Polar Endotracheal Tube (Oral) Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/136/050	5.0	6.9	240	17	10
100/136/055	5.5	7.6	252	17	10
100/136/060	6.0	8.2	265	23	10
100/136/065	6.5	8.9	280	23	10
100/136/070	7.0	9.6	300	30	10
100/136/075	7.5	10.3	315	30	10
100/136/080	8.0	10.9	330	30	10
100/136/085	8.5	11.6	345	30	10
100/136/090	9.0	12.3	360	30	10

Uncuffed Polar Preformed Tubes

ORDERING INFORMATION



Portex® North Polar Endotracheal Tube (Oral)

Uncuffed preformed tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/132/030	3.0	4.2	205	10
100/132/035	3.5	4.8	225	10
100/132/040	4.0	5.4	245	10
100/132/045	4.5	6.2	270	10
100/132/050	5.0	6.8	285	10
100/132/055	5.5	7.4	300	10
100/132/060	6.0	8.2	315	10
100/132/065	6.5	8.8	330	10
100/132/070	7.0	9.6	347	10



Portex® South Polar Endotracheal Tube (Oral)

Uncuffed preformed tracheal tube, murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/134/030	3.0	4.2	170	10
100/134/035	3.5	4.8	175	10
100/134/040	4.0	5.4	195	10
100/134/045	4.5	6.2	215	10
100/134/050	5.0	6.8	220	10
100/134/055	5.5	7.5	230	10
100/134/060	6.0	8.2	240	10
100/134/065	6.5	8.8	250	10
100/134/070	7.0	9.6	275	10
100/134/075	7.5	10.2	290	10



Portex® North Polar Polyurethane Endotracheal Tube (Nasal)

Uncuffed preformed tracheal tube, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/130/060	6.0	8.2	290	10
100/130/065	6.5	8.9	300	10
100/130/070	7.0	9.6	310	10
100/130/075	7.5	10.3	320	10
100/130/080	8.0	10.9	330	10

Speciality Tubes

ORDERING INFORMATION



Portex® Extra-long Croup Endotracheal Tube (Oral)				
Tracheal tube to help intubation where there is inflammation of the larynx and trachea, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/112/025	2.5	3.5	220	10
100/112/030	3.0	4.2	250	10
100/112/045	4.5	6.2	330	10



Portex® Cole-style Endotracheal Tube (Oral)					
Tracheal tube with stepped design for ease of placement and an effective seal, magill tip					
Item Number	I.D. (mm)	Proximal O.D. (mm)	Distal O.D. (mm)	Length (mm)	Units per Case
100/430/150	1.5	5.5	3.0	130	10
100/430/200	2.0	5.8	3.5	130	10
100/430/250	2.5	6.5	4.0	135	10
100/430/300	3.0	6.5	4.5	135	10
100/430/350	3.5	7.5	5.0	148	10
100/430/400	4.0	7.5	5.5	148	10



Portex® Microlaryngeal Ivory Endotracheal Tube (Oral / Nasal)					
Tracheal tube for use in microlaryngeal procedures, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/190/050	5.0	7.2	315	25	10

Thoracic Airways

ORDERING INFORMATION



Portex® Endobronchial Tube						
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Tracheal Cuff Ø	Bronchial Cuff Ø	Units per Case
Left-sided endobronchial tube						
198-28L	4.6	10.3	390	23	12	1
198-32L	5.4	11.6	390	26	13	1
198-35L	6.0	12.8	400	27	18	1
198-37L	6.3	13.3	400	28	18	1
198-39L	6.6	14.2	400	29	23	1
198-41L	7.0	14.8	400	31	23	1
Right-sided endobronchial tube						
197-28R	4.6	10.3	390	23	13	1
197-32R	5.4	11.6	390	26	13	1
197-35R	6.0	12.8	400	27	21	1
197-37R	6.3	13.3	400	28	21	1
197-39R	6.6	14.2	400	29	21	1
197-41R	7.0	14.8	400	31	21	1

Laryngeal Masks

ORDERING INFORMATION



Portex® Silicone Laryngeal Mask

Single piece laryngeal mask with aperture bars

Item Number	Size (mm)	I.D. Tube (mm)	O.D. Tube (mm)	Max Cuff Volume (ml)	Units per Case
100/222/200	2.0	7.0	11.5	10	10
100/222/250	2.5	8.3	13.3	14	10
100/222/300	3.0	10.0	15.0	20	10
100/222/400	4.0	10.0	15.0	30	10
100/222/500	5.0	11.5	16.5	40	10

Accessories

ORDERING INFORMATION

Portex® Laryngotracheal Anaesthesia Kit

For patients with difficult airways

Item Number	Syringe	Catheter	Units per Case
F031-61	5cc	with Luer Lock	20

Portex® Guedel Airway

Sterile and disposable without colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/320/200	00	50	NA	50
100/320/100	0	60	NA	50
100/320/010	1	70	NA	50
100/320/020	2	90	NA	50
100/320/030	3	100	NA	50
100/320/040	4	110	NA	50

Portex® Guedel Airway

Sterile and disposable with colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/322/040	000	40	Pink	10
100/322/050	00	50	Dark Blue	10
100/322/060	0	60	Black	10
100/322/070	1	70	White	10
100/322/080	2	80	Green	10
100/322/090	3	90	Yellow	10
100/322/100	4	100	Red	10
100/322/110	5	110	Light Blue	10
100/322/120	6	120	Orange	10

Portex® Berman Airway

Sterile and disposable with two channels for suctioning



Item Number	Size	Length (mm)	Units per Case
100/324/040	000	40	10
100/324/050	00	50	10
100/324/060	0	60	10
100/324/070	1	70	10
100/324/080	2	80	10
100/324/090	3	90	10
100/324/100	4	100	10
100/324/110	5	110	10

ORDERING INFORMATION



Portex® Nasopharyngeal Airway Soft and kink resistant				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/210/060	6.0	8.8	135	10
100/210/070	7.0	10.2	155	10
100/210/080	8.0	11.6	175	10
100/210/090	9.0	13.1	195	10



Portex® Intubation Stylet Standard and extra length				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Standard length stylets				
100/120/100	2.5 – 4.5	2.0	254	10
100/120/200	5.0 – 8.0	4.0	335	10
100/120/300	8.5 – 11.0	5.0	365	10
Extra length stylets				
100/121/200	5.0 – 8.0	4.0	673	10
100/121/300	8.5 – 11.0	5.0	693	10



Portex® Bougie Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable bougies				
14-504-17	6.0 – 11.0	5.0	600	1
Single-use bougies				
100/123/515	6.0 – 11.0	5.0	700	10



Portex® Exchange Guide Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable exchange guides				
14-504-76	4.0 – 5.5	3.3	700	1
14-504-68	6.0 – 11.0	5.0	700	1
Single-use exchange guides				
100/125/005	2.0 – 3.5	1.7	500	10
100/125/010	4.0 – 5.5	3.3	700	10
100/123/015	6.0 – 11.0	5.0	700	10

ORDERING INFORMATION

**PressureEasy® Cuff Pressure Indicator**

Continuously monitor tracheal cuff pressure

Item Number	Item Description	Units per Case
292003	PressureEasy® Cuff Pressure Monitor	10
292004	PressureEasy® Cuff Pressure Monitor with feedback line	10

**Portex® Cuff Inflator Pressure Gauge**

Controlled inflation and monitoring

Item Number	Item Description	Units per Case
100/568/000	Cuff Inflator Pressure Gauge with Connecting Tube	1

**Portex® Paediatric Endotracheal Tube Holder**

Secure fixation to avoid risk of accidental extubation, disconnection or complete displacement of the tube

Item Number	Recommended tube I.D. (mm)	Units per Case
Paediatric Tracheal Tube Holder		
100/603/025	2.5	10
100/603/030	3.0	10
100/603/035	3.5	10
100/603/040	4.0	10
100/603/045	4.5	10
Siliconized Tracheal Tube with Holder		
100/611/020	2.0	10
100/611/025	2.5	10
100/611/030	3.0	10
100/611/035	3.5	10
100/611/040	4.0	10
100/611/045	4.5	10

**Connectors**

Durable and Robust

Item Number	Item Description	Units per Case
100/250/001	Reusable swivel 15mm connector	10
100/251/001	Catheter mount adaptor female 15mm	10
100/255/150	Double swivel connector	10
100/257/000	Fibreoptic bronchoscope swivel connector	10
100/590/000	Catheter Mount without double swivel connector	20
100/594/000	Catheter Mount with double swivel connector	10

SYMBOLS

14-504-17 19	100/110/050 10	100/130/080 13	100/141/050 8
14-504-68 19	100/110/055 10	100/131/060 12	100/141/055 8
14-504-76 19	100/110/060 10	100/131/065 12	100/141/060 8
35W025 11	100/110/065 10	100/131/070 12	100/141/065 8
35W030 11	100/110/070 10	100/131/075 12	100/141/070 8
35W035 11	100/110/075 10	100/131/080 12	100/141/075 8
35W040 11	100/110/080 10	100/131/085 12	100/141/080 8
35W045 11	100/110/085 10	100/131/090 12	100/141/085 8
35W050 11	100/110/090 10	100/132/030 13	100/141/090 8
35W055 11	100/110/095 10	100/132/035 13	100/141/095 8
35W060 11	100/111/020 8	100/132/040 13	100/150/050 7
35W065 11	100/111/025 8	100/132/045 13	100/150/055 7
35W070 11	100/111/030 8	100/132/050 13	100/150/060 7
35W075 11	100/111/035 8	100/132/055 13	100/150/065 7
35W080 11	100/111/040 8	100/132/060 13	100/150/070 7
35W085 11	100/111/045 8	100/132/065 13	100/150/075 7
35W090 11	100/111/050 8	100/132/070 13	100/150/080 7
35W095 11	100/111/055 8	100/133/060 12	100/150/085 7
100/105/025 9	100/111/060 8	100/133/065 12	100/150/090 7
100/105/030 9	100/111/065 8	100/133/070 12	100/150/095 7
100/105/035 9	100/111/070 8	100/133/075 12	100/166/050 6
100/105/040 9	100/111/075 8	100/133/080 12	100/166/055 6
100/105/045 9	100/111/080 8	100/134/030 13	100/166/060 6
100/105/050 9	100/111/085 8	100/134/035 13	100/166/065 6
100/105/055 9	100/111/090 8	100/134/040 13	100/166/070 6
100/105/060 9	100/112/025 14	100/134/045 13	100/166/075 6
100/105/065 9	100/112/030 14	100/134/050 13	100/166/080 6
100/105/070 9	100/112/045 14	100/134/055 13	100/166/085 6
100/105/075 9	100/113/030 11	100/134/060 13	100/166/090 6
100/105/080 9	100/113/035 11	100/134/065 13	100/166/095 6
100/105/085 9	100/113/040 11	100/134/070 13	100/166/100 6
100/105/090 9	100/113/045 11	100/134/075 13	100/179/050 7
100/107/030 11	100/113/050 11	100/136/050 12	100/179/055 7
100/107/035 11	100/120/100 19	100/136/055 12	100/179/060 7
100/107/040 11	100/120/200 19	100/136/060 12	100/179/065 7
100/107/045 11	100/120/300 19	100/136/065 12	100/179/070 7
100/107/050 11	100/121/200 19	100/136/070 12	100/179/075 7
100/109/050 10	100/121/300 19	100/136/075 12	100/179/080 7
100/109/055 10	100/123/015 19	100/136/080 12	100/179/085 7
100/109/060 10	100/123/515 19	100/136/085 12	100/179/090 7
100/109/065 10	100/125/005 19	100/136/090 12	100/189/060 6
100/109/070 10	100/125/010 19	100/141/025 8	100/189/065 6
100/109/075 10	100/130/060 13	100/141/030 8	100/189/070 6
100/109/080 10	100/130/065 13	100/141/035 8	100/189/075 6
100/109/085 10	100/130/070 13	100/141/040 8	100/189/080 6
100/109/090 10	100/130/075 13	100/141/045 8	100/189/085 6

100/189/090	6	100/324/090	18	H4006S	14
100/190/050	14	100/324/100	18	H4021PS	15
100/199/050	6	100/324/110	18	H4021S	14
100/199/055	6	100/430/150	14	H4022PS	15
100/199/060	6	100/430/200	14	H4022S	14
100/199/065	6	100/430/250	14	H4023PS	15
100/199/070	6	100/430/300	14	H4023S	14
100/199/075	6	100/430/350	14	H4024S	14
100/199/080	6	100/430/400	14	H4025S	14
100/199/085	6	100/568/000	20	H4026S	14
100/199/090	6	100/569/000	20	H4053	20
100/199/095	6	100/590/000	21	H4054	20
100/199/100	6	100/594/000	21	H4055	20
100/210/060	19	100/603/025	20	H4056	20
100/210/070	19	100/603/030	20		
100/210/080	19	100/603/035	20	K	
100/210/090	19	100/603/040	20	K1000	20
100/222/200	17	100/603/045	20	K2000	20
100/222/250	17	100/611/020	20	K3000	20
100/222/300	17	100/611/025	20		
100/222/400	17	100/611/030	20		
100/222/500	17	100/611/035	20		
100/250/001	21	100/611/040	20		
100/251/001	21	100/611/045	20		
100/255/150	21	197-28R	16		
100/257/000	21	197-32R	16		
100/320/010	18	197-35R	16		
100/320/020	18	197-37R	16		
100/320/030	18	197-39R	16		
100/320/040	18	197-41R	16		
100/320/090	18	198-28L	16		
100/320/100	18	198-32L	16		
100/320/200	18	198-35L	16		
100/322/040	18	198-37L	16		
100/322/050	18	198-39L	16		
100/322/060	18	198-41L	16		
100/322/070	18	292003	20		
100/322/080	18	292004	20		
100/322/100	18				
100/322/110	18	F			
100/322/120	18	F031-61	18		
100/324/040	18	H			
100/324/050	18	H4003PS	15		
100/324/060	18	H4003S	14		
100/324/070	18	H4004S	14		
100/324/080	18	H4005S	14		

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical Australasia Pty. Ltd.
Suite 2.03, 97 Waterloo Road
Macquarie Park NSW 2113 Australia
Tel AU: 1800 654 949
www.smiths-medical.com

Distributed by:
USL Medical
494 Rosebank Road
Avondale, 1026 New Zealand
Tel NZ: 0800 658 814

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. SACETT, AIRCARE, Aire-Cuf, Bivona, PressureEasy, SecureEasy and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2018 Smiths Medical. All rights reserved. SM196408GBAUS-052018

MHYTCA-1267



smiths medical