

## General Anesthesia Airway


Intubation Systems & Accessories Catalogue



Magill Tubes



Reinforced  
Tubes



Polar Preformed  
Tubes



Speciality Tubes



Accessories





# Contents

- Cuffed Magill Tubes ..... 5**
  - SACETT® (Oral)..... 5
  - Cuffed Endotracheal Tube (Oral/Nasal) ..... 5
  - Cuffed Siliconised Endotracheal Tube (Oral/Nasal)..... 5
  - Cuffed Ivory Endotracheal Tube (Nasal)..... 6
  - Choice Cuffed Endotracheal Tube (Oral/Nasal) ..... 6
- Uncuffed Magill Tubes ..... 7**
  - Uncuffed Siliconised Endotracheal Tube (Oral/Nasal)..... 7
  - Uncuffed Siliconised Endotracheal Tube (Oral/Nasal)..... 7
  - Uncuffed Ivory Endotracheal Tube (Oral/Nasal) ..... 8
- Reinforced Tubes ..... 9**
  - Reinforced Endotracheal Tube (Oral/Nasal) ..... 9
  - Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral/Nasal) ..... 9
  - Bivona® Wire Reinforced Endotracheal Tubes, Uncuffed, Magill..... 9
- Cuffed Polar Preformed Tubes ..... 10**
  - North Polar Endotracheal Tube (Oral) ..... 10
  - North Polar Ivory Endotracheal Tube (Nasal) ..... 10
  - South Polar Endotracheal Tube (Oral) ..... 10
- Uncuffed Polar Preformed Tubes ..... 11**
  - North Polar Endotracheal Tube (Oral) ..... 11
  - South Polar Endotracheal Tube (Oral) ..... 11
  - North Polar Polyurethane Endotracheal Tube (Nasal) ..... 11
- Speciality Tubes ..... 12**
  - Extra-long Croup Endotracheal Tube (Oral)..... 12
  - Cole-style Endotracheal Tube (Oral) ..... 12
  - Microlaryngeal Ivory Endotracheal Tube (Oral/Nasal)..... 12
- Thoracic Airways..... 13**
  - Endobronchial Tube ..... 13
- Laryngeal Masks..... 14**
  - Silicone Laryngeal Mask..... 14
  - CO<sub>2</sub> Absorbents ..... 14

<b>Accessories</b> .....	15
Laryngotracheal Anaesthesia Kit .....	15
Guedel Airway .....	15
Guedel Airway .....	15
Berman Airway .....	15
Nasopharyngeal Airway.....	16
Intubation Stylet.....	16
Bougie .....	16
Exchange Guide .....	16
Cuff Inflator Pressure Gauge.....	17
Paediatric Endotracheal Tube Holder .....	17
Connectors.....	17

## Cuffed Magill Tubes

### ORDERING INFORMATION



<b>Blue Line® SACETT® (Oral)</b> Tracheal tube with additional lumen for drainage of the subglottic region, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/189/060	6.0	9.1	290	23	10
100/189/065	6.5	9.7	300	23	10
100/189/070	7.0	10.4	310	30	10
100/189/075	7.5	11.3	320	30	10
100/189/080	8.0	11.8	330	30	10
100/189/085	8.5	12.3	330	30	10
100/189/090	9.0	13.0	330	30	10




<b>Cuffed Endotracheal Tube (Oral/Nasal)</b> Tracheal tube for long-term intubation and short-term ICU stays, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/199/050	5.0	6.9	250	17	20
100/199/055	5.5	7.6	280	17	20
100/199/060	6.0	8.2	290	23	20
100/199/065	6.5	8.9	300	23	20
100/199/070	7.0	9.6	310	30	20
100/199/075	7.5	10.3	320	30	20
100/199/080	8.0	10.9	330	30	20
100/199/085	8.5	11.6	330	30	20
100/199/090	9.0	12.3	330	30	20
100/199/095	9.5	13.0	330	30	20
100/199/100	10.0	13.7	330	34	20



<b>Cuffed Siliconised Endotracheal Tube (Oral/Nasal)</b> Tracheal tube coated in silicone, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/166/050	5.0	6.9	250	17	10
100/166/055	5.5	7.6	280	17	10
100/166/060	6.0	8.2	290	23	10
100/166/065	6.5	8.9	300	23	10
100/166/070	7.0	9.6	310	30	10
100/166/075	7.5	10.3	320	30	10
100/166/080	8.0	10.9	330	30	10
100/166/085	8.5	11.6	330	30	10
100/166/090	9.0	12.3	330	30	10
100/166/095	9.5	13.0	330	30	10
100/166/100	10.0	13.7	330	34	10

## ORDERING INFORMATION



### Cuffed Ivory Endotracheal Tube (Nasal)

Tracheal tube for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/179/050	5.0	7.3	250	17	10
100/179/055	5.5	8.0	280	17	10
100/179/060	6.0	8.8	290	23	10
100/179/065	6.5	9.5	300	23	10
100/179/070	7.0	10.2	310	30	10
100/179/075	7.5	10.9	320	30	10
100/179/080	8.0	11.6	330	30	10
100/179/085	8.5	12.4	330	30	10
100/179/090	9.0	13.1	330	30	10



### Choice Cuffed Endotracheal Tube (Oral/Nasal)

Murphy eye

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/150/050	5.0	6.7	265	22	10
100/150/055	5.5	7.3	295	22	10
100/150/060	6.0	8.0	300	24	10
100/150/065	6.5	8.5	310	24	10
100/150/070	7.0	9.2	320	26	10
100/150/075	7.5	10.0	325	26	10
100/150/080	8.0	10.7	325	28	10
100/150/085	8.5	11.3	325	28	10
100/150/090	9.0	12.0	325	28	10
100/150/095	9.5	12.7	325	28	10

# Uncuffed Magill Tubes

## ORDERING INFORMATION



<b>Uncuffed Siliconised Endotracheal Tube (Oral/Nasal)</b> Tracheal tube coated in silicone, magill tip				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/111/020	2.0	2.9	140	10
100/111/025	2.5	3.4	150	10
100/111/030	3.0	4.2	170	10
100/111/035	3.5	4.8	190	10
100/111/040	4.0	5.4	210	10
100/111/045	4.5	6.2	230	10
100/111/050	5.0	6.9	250	10
100/111/055	5.5	7.6	280	10
100/111/060	6.0	8.2	290	10
100/111/065	6.5	8.9	300	10
100/111/070	7.0	9.6	310	10
100/111/075	7.5	10.3	320	10
100/111/080	8.0	10.9	330	10
100/111/085	8.5	11.6	330	10
100/111/090	9.0	12.3	330	10



<b>Uncuffed Siliconised Endotracheal Tube (Oral/Nasal)</b> Tracheal tube coated in silicone, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/141/025	2.5	3.4	150	10
100/141/030	3.0	4.2	170	10
100/141/035	3.5	4.8	190	10
100/141/040	4.0	5.4	210	10
100/141/045	4.5	6.2	230	10
100/141/050	5.0	6.9	250	10
100/141/055	5.5	7.6	280	10
100/141/060	6.0	8.2	290	10
100/141/065	6.5	8.9	300	10
100/141/070	7.0	9.6	310	10
100/141/075	7.5	10.3	320	10
100/141/080	8.0	10.9	330	10
100/141/085	8.5	11.6	330	10
100/141/090	9.0	12.3	330	10

## ORDERING INFORMATION



### Uncuffed Ivory Endotracheal Tube (Oral/Nasal)


Tracheal tube for nasal intubation, magill tip

Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/105/025	2.5	3.7	150	10
100/105/030	3.0	4.4	170	10
100/105/035	3.5	5.1	190	10
100/105/040	4.0	5.9	210	10
100/105/045	4.5	6.6	230	10
100/105/050	5.0	7.3	250	10
100/105/055	5.5	8.0	280	10
100/105/060	6.0	8.8	290	10
100/105/065	6.5	9.5	300	10
100/105/070	7.0	10.2	310	10
100/105/080	8.0	11.6	330	10



## Reinforced Tubes

### ORDERING INFORMATION



<b>Reinforced Endotracheal Tube (Oral/Nasal)</b> Cuffed single-use tracheal tube, Murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/110/050	5.0	8.2	259	17	5
100/110/055	5.5	8.6	287	17	5
100/110/060	6.0	9.1	297	23	5
100/110/065	6.5	9.8	307	23	5
100/110/070	7.0	10.4	318	30	5
100/110/075	7.5	11.1	328	30	5
100/110/080	8.0	11.9	338	30	5
100/110/085	8.5	12.4	338	30	5
100/110/090	9.0	13.1	338	30	5
100/110/095	9.5	13.7	338	30	5



<b>Bivona® Cuffed Wire Reinforced Endotracheal Tube (Oral/Nasal)</b> Reusable tracheal tube with TTS™ or small diameter Aire-Cuf®, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
<b>TTS™</b>					
35W040	4.0	6.0	200		1
35W045	4.5	6.7	220		1
35W050	5.0	7.3	240		1
35W055	5.5	8.0	260		1
<b>Small Diameter Aire-Cuf®, magill tip</b>					
35W060	6.0	8.7	280	11	1
35W065	6.5	9.4	280	11	1
35W070	7.0	10.0	300	13	1



<b>Bivona® Wire Reinforced Endotracheal Tubes, Uncuffed, Magill</b>				
Product Code	I.D. (mm)	O.D. (mm)	Length (mm)	Units Per Case
15W045	4.5	6.7	220	1
15W050	5.0	7.3	240	1
15W055	5.5	8.7	280	1

## Cuffed Polar Preformed Tubes

### ORDERING INFORMATION



<b>North Polar Endotracheal Tube (Oral)</b> Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/131/060	6.0	8.2	290	23	10
100/131/065	6.5	8.9	300	23	10
100/131/070	7.0	9.6	310	30	10
100/131/075	7.5	10.3	320	30	10
100/131/080	8.0	10.9	330	30	10
100/131/085	8.5	11.6	330	30	10
100/131/090	9.0	12.3	330	30	10



<b>North Polar Ivory Endotracheal Tube (Nasal)</b> Cuffed preformed tracheal tube, magill tip					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/133/060	6.0	8.8	270	23	10
100/133/065	6.5	9.4	275	23	10
100/133/070	7.0	10.2	285	30	10
100/133/075	7.5	10.8	295	30	10
100/133/080	8.0	11.6	305	30	10



<b>South Polar Endotracheal Tube (Oral)</b> Cuffed preformed tracheal tube, murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/136/050	5.0	6.9	240	17	10
100/136/055	5.5	7.6	252	17	10
100/136/060	6.0	8.2	265	23	10
100/136/065	6.5	8.9	280	23	10
100/136/070	7.0	9.6	300	30	10
100/136/075	7.5	10.3	315	30	10
100/136/080	8.0	10.9	330	30	10
100/136/085	8.5	11.6	345	30	10
100/136/090	9.0	12.3	360	30	10

## Uncuffed Polar Preformed Tubes

### ORDERING INFORMATION



<b>North Polar Endotracheal Tube (Oral)</b> Uncuffed preformed tracheal tube, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/132/030	3.0	4.2	205	10
100/132/035	3.5	4.8	225	10
100/132/040	4.0	5.4	245	10
100/132/045	4.5	6.2	270	10
100/132/050	5.0	6.8	285	10
100/132/055	5.5	7.4	300	10
100/132/060	6.0	8.2	315	10
100/132/065	6.5	8.8	330	10
100/132/070	7.0	9.6	347	10



<b>South Polar Endotracheal Tube (Oral)</b> Uncuffed preformed tracheal tube, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/134/030	3.0	4.2	170	10
100/134/035	3.5	4.8	175	10
100/134/040	4.0	5.4	195	10
100/134/045	4.5	6.2	215	10
100/134/050	5.0	6.8	220	10
100/134/055	5.5	7.5	230	10
100/134/060	6.0	8.2	240	10
100/134/065	6.5	8.8	250	10
100/134/070	7.0	9.6	275	10
100/134/075	7.5	10.2	290	10



<b>North Polar Polyurethane Endotracheal Tube (Nasal)</b> Uncuffed preformed tracheal tube, magill tip				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/130/060	6.0	8.2	290	10
100/130/065	6.5	8.9	300	10
100/130/070	7.0	9.6	310	10
100/130/080	8.0	10.9	330	10

## Speciality Tubes

### ORDERING INFORMATION



<b>Extra-long Croup Endotracheal Tube (Oral)</b>				
Tracheal tube for intubation where there is inflammation of the larynx and trachea, murphy eye				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/112/025	2.5	3.5	220	10
100/112/030	3.0	4.2	250	10
100/112/045	4.5	6.2	330	10



<b>Cole-style Endotracheal Tube (Oral)</b>					
Tracheal tube with stepped design for placement and seal, magill tip					
Item Number	I.D. (mm)	Proximal O.D. (mm)	Distal O.D. (mm)	Length (mm)	Units per Case
100/430/150	1.5	5.5	3.0	130	10
100/430/200	2.0	5.8	3.5	130	10
100/430/250	2.5	6.5	4.0	135	10
100/430/300	3.0	6.5	4.5	135	10
100/430/350	3.5	7.5	5.0	148	10
100/430/400	4.0	7.5	5.5	148	10



<b>Microlaryngeal Ivory Endotracheal Tube (Oral/Nasal)</b>					
Murphy eye					
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø	Units per Case
100/190/050	5.0	7.2	315	25	10

# Thoracic Airways

## ORDERING INFORMATION

Endobronchial Tube						
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Tracheal Cuff Ø	Bronchial Cuff Ø	Units per Case
Left-sided endobronchial tube						
198-28L	4.6	10.3	390	23	12	1
198-32L	5.4	11.6	390	26	13	1
198-35L	6.0	12.8	400	27	18	1
198-37L	6.3	13.3	400	28	18	1
198-39L	6.6	14.2	400	29	23	1
198-41L	7.0	14.8	400	31	23	1
Right-sided endobronchial tube						
197-28R	4.6	10.3	390	23	13	1
197-32R	5.4	11.6	390	26	13	1
197-35R	6.0	12.8	400	27	21	1
197-37R	6.3	13.3	400	28	21	1
197-39R	6.6	14.2	400	29	21	1
197-41R	7.0	14.8	400	31	21	1



## Laryngeal Masks

### ORDERING INFORMATION



<b>Silicone Laryngeal Mask</b> Single piece laryngeal mask with aperture bars					
Item Number	Size (mm)	I.D. Tube (mm)	O.D. Tube (mm)	Max Cuff Volume (ml)	Units per Case
100/222/200	2.0	7.0	11.5	10	10
100/222/250	2.5	8.3	13.3	14	10
100/222/300	3.0	10.0	15.0	20	10
100/222/400	4.0	10.0	15.0	30	10
100/222/500	5.0	11.5	16.5	40	10



<b>Sodasorb CO<sub>2</sub> Absorbents</b>		
Product Code	Description	Units Per Case
008860	Canister Pack (Refill Bag, 1.64 L)	12
008850	Bottle (5 L)	5
008865	LF CO <sub>2</sub> Canister Pack (Refill Bag, 1.64 L)	12
008875	LF CO <sub>2</sub> Pre Pack, Canister Disks (1.3 L)	12
008855	LF CO <sub>2</sub> Pail (1 gallon, 4.3 L or 4.2 kg)	3

## Accessories

### ORDERING INFORMATION

#### Laryngotracheal Anaesthesia Kit

For patients with difficult airways

Item Number	Syringe	Catheter	Units per Case
F031-61	5cc	with Luer Lock	20

#### Guedel Airway

Disposable without colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/320/100	0	60	NA	50
100/320/010	1	70	NA	50
100/320/020	2	90	NA	50
100/320/030	3	100	NA	50
100/320/040	4	110	NA	50

#### Guedel Airway

Disposable with colour coding



Item Number	Size	Length (mm)	Colour Code	Units per Case
100/322/040	000	40	Pink	10
100/322/050	00	50	Dark Blue	10
100/322/060	0	60	Black	10
100/322/070	1	70	White	10
100/322/080	2	80	Green	10
100/322/090	3	90	Yellow	10
100/322/100	4	100	Red	10
100/322/110	5	110	Light Blue	10
100/322/120	6	120	Orange	10

#### Berman Airway

Sterile and disposable with two channels for suctioning




Item Number	Size	Length (mm)	Units per Case
100/324/040	000	40	10
100/324/050	00	50	10
100/324/060	0	60	10
100/324/070	1	70	10
100/324/080	2	80	10
100/324/090	3	90	10
100/324/100	4	100	10
100/324/110	5	110	10

ORDERING INFORMATION



Nasopharyngeal Airway				
Item Number	I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
100/210/060	6.0	8.8	135	10
100/210/070	7.0	10.2	155	10
100/210/080	8.0	11.6	175	10
100/210/090	9.0	13.1	195	10



Intubation Stylet Standard and extra length				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Standard length stylets				
100/120/100	2.5 – 4.5	2.0	254	10
100/120/200	5.0 – 8.0	4.0	335	10
100/120/300	8.5 – 11.0	5.0	365	10
Extra length stylets				
100/121/200	5.0 – 8.0	4.0	673	10
100/121/300	8.5 – 11.0	5.0	693	10



Bougie Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable bougies				
14-504-17	6.0 – 11.0	5.0	600	1
Single-use bougies				
100/123/515	6.0 – 11.0	5.0	700	10



Exchange Guide Single-use and reusable versions				
Item Number	Recommended tube I.D. (mm)	O.D. (mm)	Length (mm)	Units per Case
Reusable exchange guides				
14-504-76	4.0 – 5.5	3.3	700	1
14-504-68	6.0 – 11.0	5.0	700	1
Single-use exchange guides				
100/125/005	2.0 – 3.5	1.7	500	10
100/125/010	4.0 – 5.5	3.3	700	10



## ORDERING INFORMATION



Cuff Inflator Pressure Gauge		
Item Number	Item Description	Units per Case
100/568/000	Cuff Inflator Pressure Gauge with Connecting Tube	1



Paediatric Endotracheal Tube Holder		
Item Number	Recommended tube I.D. (mm)	Units per Case
Paediatric Tracheal Tube Holder		
100/603/025	2.5	10
100/603/030	3.0	10
100/603/035	3.5	10
100/603/040	4.0	10
100/603/045	4.5	10
69-1200	-	100
Siliconized Tracheal Tube with Holder		
100/611/030	3.0	10
100/611/045	4.5	10



100/257/000

Connectors		
Item Number	Item Description	Units per Case
100/250/001	Reusable swivel 15 mm connector	10
100/251/001	Catheter mount adaptor female 15 mm	10
100/255/150	Double swivel connector	10
100/257/000	Fibreoptic bronchoscope swivel connector	10
100/590/000	Catheter Mount without double swivel connector	20
100/594/000	Catheter Mount with double swivel connector	10



100/594/000

14-504-17 ..... 16	100/112/025 ..... 12	100/136/055 ..... 10	100/179/070 ..... 6	100/322/080 ..... 15
14-504-68 ..... 16	100/112/030 ..... 12	100/136/060 ..... 10	100/179/075 ..... 6	100/322/100 ..... 15
14-504-76 ..... 16	100/112/045 ..... 12	100/136/065 ..... 10	100/179/080 ..... 6	100/322/110 ..... 15
35W040 ..... 9	100/120/100 ..... 16	100/136/070 ..... 10	100/179/085 ..... 6	100/322/120 ..... 15
35W045 ..... 9	100/120/200 ..... 16	100/136/075 ..... 10	100/179/090 ..... 6	100/324/040 ..... 15
35W050 ..... 9	100/120/300 ..... 16	100/136/080 ..... 10	100/189/060 ..... 5	100/324/050 ..... 15
35W055 ..... 9	100/121/200 ..... 16	100/136/085 ..... 10	100/189/065 ..... 5	100/324/060 ..... 15
35W060 ..... 9	100/121/300 ..... 16	100/136/090 ..... 10	100/189/070 ..... 5	100/324/070 ..... 15
35W065 ..... 9	100/123/515 ..... 16	100/141/025 ..... 7	100/189/075 ..... 5	100/324/080 ..... 15
35W070 ..... 9	100/125/005 ..... 16	100/141/030 ..... 7	100/189/080 ..... 5	100/324/090 ..... 15
69-1200 ..... 17	100/125/010 ..... 16	100/141/035 ..... 7	100/189/085 ..... 5	100/324/100 ..... 15
100/105/025 ..... 8	100/130/060 ..... 11	100/141/040 ..... 7	100/189/090 ..... 5	100/324/110 ..... 15
100/105/030 ..... 8	100/130/065 ..... 11	100/141/045 ..... 7	100/190/050 ..... 12	100/430/150 ..... 12
100/105/035 ..... 8	100/130/070 ..... 11	100/141/050 ..... 7	100/199/050 ..... 5	100/430/200 ..... 12
100/105/040 ..... 8	100/130/080 ..... 11	100/141/055 ..... 7	100/199/055 ..... 5	100/430/250 ..... 12
100/105/045 ..... 8	100/131/060 ..... 10	100/141/060 ..... 7	100/199/060 ..... 5	100/430/300 ..... 12
100/105/050 ..... 8	100/131/065 ..... 10	100/141/065 ..... 7	100/199/065 ..... 5	100/430/350 ..... 12
100/105/055 ..... 8	100/131/070 ..... 10	100/141/070 ..... 7	100/199/070 ..... 5	100/430/400 ..... 12
100/105/060 ..... 8	100/131/075 ..... 10	100/141/075 ..... 7	100/199/075 ..... 5	100/568/000 ..... 17
100/105/065 ..... 8	100/131/080 ..... 10	100/141/080 ..... 7	100/199/080 ..... 5	100/590/000 ..... 17
100/105/070 ..... 8	100/131/085 ..... 10	100/141/085 ..... 7	100/199/085 ..... 5	100/594/000 ..... 17
100/105/080 ..... 8	100/131/090 ..... 10	100/141/090 ..... 7	100/199/090 ..... 5	100/603/025 ..... 17
100/110/050 ..... 9	100/132/030 ..... 11	100/150/050 ..... 6	100/199/095 ..... 5	100/603/030 ..... 17
100/110/055 ..... 9	100/132/035 ..... 11	100/150/055 ..... 6	100/199/100 ..... 5	100/603/035 ..... 17
100/110/060 ..... 9	100/132/040 ..... 11	100/150/060 ..... 6	100/210/060 ..... 16	100/603/040 ..... 17
100/110/065 ..... 9	100/132/045 ..... 11	100/150/065 ..... 6	100/210/070 ..... 16	100/603/045 ..... 17
100/110/070 ..... 9	100/132/050 ..... 11	100/150/070 ..... 6	100/210/080 ..... 16	100/611/030 ..... 17
100/110/075 ..... 9	100/132/055 ..... 11	100/150/075 ..... 6	100/210/090 ..... 16	100/611/045 ..... 17
100/110/080 ..... 9	100/132/060 ..... 11	100/150/080 ..... 6	100/222/200 ..... 14	197-28R ..... 13
100/110/085 ..... 9	100/132/065 ..... 11	100/150/085 ..... 6	100/222/250 ..... 14	197-32R ..... 13
100/110/090 ..... 9	100/132/070 ..... 11	100/150/090 ..... 6	100/222/300 ..... 14	197-35R ..... 13
100/110/095 ..... 9	100/133/060 ..... 10	100/150/095 ..... 6	100/222/400 ..... 14	197-37R ..... 13
100/111/020 ..... 7	100/133/065 ..... 10	100/166/050 ..... 5	100/222/500 ..... 14	197-39R ..... 13
100/111/025 ..... 7	100/133/070 ..... 10	100/166/055 ..... 5	100/250/001 ..... 17	197-41R ..... 13
100/111/030 ..... 7	100/133/075 ..... 10	100/166/060 ..... 5	100/251/001 ..... 17	198-28L ..... 13
100/111/035 ..... 7	100/133/080 ..... 10	100/166/065 ..... 5	100/255/150 ..... 17	198-32L ..... 13
100/111/040 ..... 7	100/134/030 ..... 11	100/166/070 ..... 5	100/257/000 ..... 17	198-35L ..... 13
100/111/045 ..... 7	100/134/035 ..... 11	100/166/075 ..... 5	100/320/010 ..... 15	198-37L ..... 13
100/111/050 ..... 7	100/134/040 ..... 11	100/166/080 ..... 5	100/320/020 ..... 15	198-39L ..... 13
100/111/055 ..... 7	100/134/045 ..... 11	100/166/085 ..... 5	100/320/030 ..... 15	198-41L ..... 13
100/111/060 ..... 7	100/134/050 ..... 11	100/166/090 ..... 5	100/320/040 ..... 15	
100/111/065 ..... 7	100/134/055 ..... 11	100/166/095 ..... 5	100/320/090 ..... 15	<b>F</b>
100/111/070 ..... 7	100/134/060 ..... 11	100/166/100 ..... 5	100/320/100 ..... 15	F031-61 ..... 15
100/111/075 ..... 7	100/134/065 ..... 11	100/179/050 ..... 6	100/322/040 ..... 15	
100/111/080 ..... 7	100/134/070 ..... 11	100/179/055 ..... 6	100/322/050 ..... 15	
100/111/085 ..... 7	100/134/075 ..... 11	100/179/060 ..... 6	100/322/060 ..... 15	
100/111/090 ..... 7	100/136/050 ..... 10	100/179/065 ..... 6	100/322/070 ..... 15	



For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical Australasia Pty. Ltd.**  
Suite 2.03, 97 Waterloo Road  
Macquarie Park NSW 2113 Australia  
Tel AU: 1800 654 949  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Distributed by:**  
USL Medical  
494 Rosebank Road  
Avondale, 1026 New Zealand  
Tel NZ: 0800 658 814

MHYTCA-1295

**smiths medical**

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. SACETT, Aire-Cuf, Bivona, PressureEasy, SecureEasy, SODASORB and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2019 Smiths Medical. All rights reserved. SM196408GBAUS-062019