

smiths medical
jelco®



Safety Blood Collection Systems




SAFETY BLOOD COLLECTION SYSTEMS

Smiths Medical provides a full line of advanced blood collection safety products for such clinical applications as venipuncture, winged needle draws, syringe draws with needleless transfer, and capillary sampling.

Multi-Sampling Blood Collection Needles

Multi-Sampling Blood Collection Needles			
Code	Description	Hub Color	Units per case
110201020006	20G x 1" Blood Collection Needle (10 Boxes of 100)	Yellow	1000
110201020007	20G x 1½" Blood Collection Needle (10 Boxes of 100)	Yellow	1000
110201020004	21G x 1" Blood Collection Needle (10 Boxes of 100)	Teal	1000
110201020005	21G x 1½" Blood Collection Needle (10 Boxes of 100)	Teal	1000
110201020002	22G x 1" Blood Collection Needle (10 Boxes of 100)	Black	1000
110201020003	22G x 1½" Blood Collection Needle (10 Boxes of 100)	Black	1000



Saf-T Wing® Blood Collection Sets

SAFETY FOR DELICATE VENOUS BLOOD DRAWS

- End-to-end protection
- Minimizes unwanted exposures by preventing needle removal and holder re-use



Saf-T Wing® Blood Collection Sets			
Code	Description	Hub Color	Units per case
972106	21G x ¾" with 6" Tubing (4 Boxes of 50)	Green	200
972112	21G x ¾" with 12" Tubing (4 Boxes of 50)	Green	200
972306	23G x ¾" with 6" Tubing (4 Boxes of 50)	Blue	200
972312	23G x ¾" with 12" Tubing (4 Boxes of 50)	Blue	200
972506	25G x ¾" with 6" Tubing (4 Boxes of 50)	Orange	200
972512	25G x ¾" with 12" Tubing (4 Boxes of 50)	Orange	200
982106	21G x ¾" with 6" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Green	200
982112	21G x ¾" with 12" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Green	200
982306	23G x ¾" with 6" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Blue	200
982312	23G x ¾" with 12" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Blue	200
982506	25G x ¾" with 6" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Orange	200
982512	25G x ¾" with 12" Tubing and Saf-T Holder® Needle Protection Device (4 Boxes of 50)	Orange	200

Venipuncture Needle-Pro® Safety Devices

SAFETY FOR VENOUS BLOOD DRAWS

- Minimizes unwanted exposures by preventing needle removal and holder re-use
- Activates with a simple one-handed technique against a firm surface
- Needle cannot be removed once safety is activated
- Compatible with Multi-Sampling Blood Collection Needles



Venipuncture Needle-Pro® Safety Devices		
Code	Description	Units per case
4140	Individual Bags (40 Bags of 25)	1000
4141	Bulk Pack (2 Bags of 500)	1000
4142	Dispenser Box (5 Boxes of 100)	500

Blood Draw Hypodermic Needle-Pro® Safety Devices

SAFETY FOR VENOUS BLOOD DRAWS

- Safety device allows bevel rotation
- Safety feature activation consistent with all Needle-Pro® Safety Devices
- Compatible with Luer slip and Luer lock syringes



Blood Draw Hypodermic Needle-Pro® Safety Devices			
Code	Description	Hub Color	Units per case
4270	20G x 1 1/2" Needle	Yellow	800
4271	20G x 1" Needle	Yellow	800
4272	21G x 1 1/2" Needle	Green	800
4273	21G x 1" Needle	Green	800
4274	22G x 1 1/2" Needle	Grey	800
4275	22G x 1" Needle	Grey	800
4276	23G x 1" Needle	Blue	800
4277	25G x 5/8" Needle	Orange	800
4278	25G x 1" Needle	Orange	800

Saf-T Holder® and Blood Transfer Devices

SAFETY FOR VENOUS BLOOD DRAWS

- Back-end needle and Luer permanently attached to the holder
- Provide needleless transfer after a syringe draw
- Available separately or with Saf-T Wing® Blood Collection Sets




Multi-Sampling Blood Collection Needles

Code	Description	Units per case
96000	Saf-T Holder® Needle Protection Device (4 Boxes of 50)	200
96000S	Short Saf-T Holder® Needle Protection Device (4 Boxes of 50)	200
96002	Saf-T Holder® Needle Protection Device with Female Luer Adapter (4 Boxes of 50)	200
96002S	Short Saf-T Holder® Needle Protection Device with Female Luer Adapter (4 Boxes of 50)	200
96004	Saf-T Holder® Blood Culture Needle Protection Device with Male Luer	50
96005	Saf-T Holder® Blood Culture Needle Protection Device with Female Luer	50
96007	Pediatric Tube Insert (10 Boxes of 10)	100

All Saf-T Holder® Needle Protection Devices incorporate a multi-Luer adapter

Safe-T-Lance® Plus Safety Lancets

Code	Description	Blade Width	Needle Gauge	Hub Color	Units per case
1005-1	Maximum low	1.5 mm	n/a	Purple	3600
1008-1	High Flow	n/a	18G	Yellow	3600
1014-1	Medium Flow	n/a	21G	Green	3600
1020-1	Low Flow	n/a	25G	Blue	3600
1025-1	Micro Flow	n/a	28G	Light Blue	3600
1028-1	Pediatric	1.5 mm	n/a	Pink	3600



neoheel™ Safety Lancets

Description	Item Number	Incision Depth	Incision Length	Blood Volume	Colour	Units per Case
Maximum low	1050M	0.65 mm	1.4 mm	Medium	Green	200
High Flow	1051P	0.85 mm	1.75 mm	Medium	Orange	200
Medium Flow	1052N	1.00 mm	2.5 mm	Medium	Yellow	200
Low Flow	1053T	2.00 mm	3.00 mm	High	Light Blue	200



For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN ALL GEOGRAPHIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Phone: 1-219-989-3297
Toll-Free USA 1-866-216-8808
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Jelco and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. There is no connection between Smiths Medical and Intuitive Surgical. The product referenced is CE-marked. ©2019 Smiths Medical. All rights reserved. SS195213EN-092019

MKEECA-0226

**Rx
ONLY**

smiths medical