

SOL-M™ Blunt Fill Needle and Blunt Fill Needle with 5 Micron Filter

A Class IIa Medical Device that meets the requirements of the 93/42/EEC
Sterile hypodermic needle for single use (ISO 7864:1993)

SOL-M™ Blunt Fill Needle – BN1815

Sol-M Blunt Fill Needle is used in conjunction with a syringe to penetrate rubber-stoppered medicine vials and pre-slit septum covering injection sites.

Features and Benefits

18G cannula facilitates the draw-up of medications

45° bevel to reduce the risk of needlestick injuries and to reduce vial stopper coring

Cannula is pre-lubricated to reduce septum insertion forces

Not made with natural rubber latex. DEHP/PVC free

Conventional standard hub to fit standard Luer lock syringes

Packaging

Shipper 1000 units per case (545 mm x 188 mm x 116 mm)

Shelf Pack 100 each (105 mm x 84 mm x 86 mm)

Expiration 5 years



SOL-M™ Blunt Fill Needle with 5 micron filter – BN1815F

Sol-M Blunt Fill Needle with filter features a 5 micron filter membrane for use when preparing medications from ampoule.

Features and Benefits

18G cannula facilitates draw-up of medications

45° bevel to reduce the risk of needlestick injuries and to reduce vial stopper coring

Cannula is pre-lubricated to reduce insertion forces

Conventional standard hub to fit standard Luer lock syringes

Contains 5 micron filter membrane to capture 98% of particulate matter


Not made with natural rubber latex. DEHP/PVC free

Packaging

Shipper 1000 units per case (545 mm x 188 mm x 116 mm)

Shelf Pack 100 each (105 mm x 84 mm x 86 mm)

Expiration 5 years



5 Micron Filter

The Blunt Fill Needle is intended for use in combination with a standard Luer lock or Luer slip syringe for draw-up of medication. MDD Classification IIa.

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN, 55442, USA
Tel: 1-614-210-7300
Toll-Free USA 1-800-258-5361
www.smiths-medical.com

Manufactured by:
SOL-Millennium®
Pudong, Shanghai (China) 200120

MKEECA-0217

CE Rx
0473 ONLY

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Jelco and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. SS195236EN-042019